Biopsy pads

Intended use

For in vitro diagnostic use only. Biopsy pads are support material intended for holding tissue biopsies in place.

Introduction

Trajan biopsy pads are drawn from fully reticulated cellular polyester urethane foam which allows optimum reagent transfer around tissue specimen whilst securing it. The foam is soft and flexible with fine pores, designed to prevent loss of specimen during processing (e.g., fixation, dehydration, infiltration). The size is compatible with all routine size tissue cassettes used in histopathology laboratories.

Important

⚠ Precaution

- Single use – do not reuse.
- Store in a dry place.
- Do not store in direct sunlight.
- For laboratory use only by trained personnel and professional use.
- Some leaching of dye may occur.

Instructions

1. Trajan biopsy pads should be gently separated prior to use. They are compatible with all routine sized tissue cassettes used in histopathology laboratories. The pads prevent tissue from being lost during processing.
2. To avoid artefacts, biopsy pads should be pre-soaked with 10% neutral buffered formalin or a fixative of choice before use.
3. One or two biopsy pads should be used to line the tissue cassette and secure the specimen. The cassette can be secured with either a plastic or metal lid.
4. When the tissue processing is complete, carefully open the cassette and extract tissue whilst separating the pads gently.
5. After tissue is extracted, dispose of biopsy pads as per laboratory and local regulations.

Specifications

Dimensions (length x width)	31 mm x 26 mm
Thickness	2.5 mm
Material	100% reticulated polyester urethane

Storage and handling

- Store the product in original packaging, in a dry environment, preferably protected from direct sources of heat and light.
- Store at room temperature – Recommended storage temperature between 2°C and 30°C.
- Avoid large variations in temperature conditions during storage.
- Do not use after use-before date.
- Shelf-life is indicated by the use-before date, printed on outer packaging.
- Trajan biopsy pads resist temperatures from -40°C to 121°C and are autoclavable.

Care and maintenance

If tearing of biopsy pads occurs, or if suboptimal performance is noticed and cannot be attributed to variations in laboratory procedures, please contact your regional Trajan representative. Some leaching of dye may occur during processing.

Information and support

Visit www.trajanscimed.com or contact techsupport@trajanscimed.com

Specifications are subject to change without notice.

	Trajan Scientific Australia Pty Ltd (Sponsor) 7 Argent Place, Ringwood, VIC 3134, AU. Tel: +61 (0) 3 9874 8577
EC REP	Axel Semrau GmbH Stefansbecke 42, 45549 Sprockhövel, Germany
US AGENT	Trajan Scientific Americas Inc 1421 W Wells Branch Parkway, Suite 108, Austin TX 78660, USA Tel: +1 512 837 7190
Made in Australia	

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