
Product Manual

StemTAG™ Alkaline Phosphatase Activity Assay Kit (Colorimetric)

Catalog Number

CBA-301

100 assays

FOR RESEARCH USE ONLY
Not for use in diagnostic procedures



CELL BIOLABS, INC.
Creating Solutions for Life Science Research

Introduction

Embryonic stem (ES) cells are continuous proliferating stem cell lines of embryonic origin first isolated from the inner cell mass (ICM). Two distinguishing features of ES cells are their ability to be maintained indefinitely in an undifferentiated state and their potential to develop into any cell within the body. Based on previous methods developed for mouse ES cells, human ES cell lines were first established by Dr. James Thomson and colleagues. Like mouse ES cells, human ES cells express high levels of membrane alkaline phosphatase (AP) and Oct-4, a transcriptional factor critical to ICM and germline formation. However, unlike mouse ES cells, hES cells do not express stage-specific embryonic antigen (SSEA-1). In addition, prolonged propagation of hES cells is typically achieved by coculture with primary mouse embryonic fibroblasts (MEFs) serving as feeder cells. Human ES cell lines are not able to maintain their undifferentiated state in the absence of supporting feeder layer cells, even when exogenous cytokines such as leukemia inhibitory factor (LIF) and gelatin-coated plates are used.

Marker Name	Mouse ES Cells	Mouse EG Cells	Human ES Cells	Human EG Cells	Human EC Cells
AP	√	√	√	√	√
SSEA-1	√	√	–	√	–
SSEA-4	–	–	√	√	√
TRA-1-60	–	–	√	√	√
TRA-1-81	–	–	√	√	√
Oct-4	√	√	√	unknown	√
ES Cell = Embryonic stem cell EG Cell = Embryonic germ cell EC Cell = Embryonic carcinoma cell					

Table 1. Comparison of Mouse and Human Pluripotent Stem Cells.

Although stem cells from different origins require different growth conditions for self-renewal and display different cell surface markers (see Table 1), AP is the most widely used stem cell marker. The StemTAG™ Alkaline Phosphatase Activity Assay Kit provides an efficient system for monitoring ES cell undifferentiation/ differentiation through AP activity by quantitative assay.

Related Products

1. CBA-300: StemTAG™ Alkaline Phosphatase Staining Kit
2. CBA-312: MEF Feeder Cells (Puromycin-resistant)
3. CBA-316: SNL Feeder Cells
4. CBA-320: CytoSelect™ 96-Well Hematopoietic Colony Forming Cell Assay

Kit Components (shipped at room temperature)

1. StemTAG™ AP Activity Assay Substrate (Part No. 30004): One bottle – 5 mL
2. Cell Lysis Buffer (Part No. 30005): One bottle – 20 mL
3. 10X Stop Solution (Part No. 30006): One bottle – 10 mL
4. AP Activity Assay Standard (Part No. 30007): One tube – 1 mL of 5 mM p-Nitrophenol

Materials Not Supplied

1. Human or Mouse Embryonic Stem Cells and Culture Medium
2. 1X PBS
3. 1X PBST (1X PBS containing 0.05% Tween-20)

Storage

Store all components at 4°C.

Preparation of Reagents

- 1X Stop Solution: Prepare a 1X Stop Solution by diluting the provided 10X stock 1:10 in deionized water. Store the diluted solution at room temperature.

Preparation of Standard Curve

1. Prepare a 10-fold dilution of the AP Activity Assay Standard (5 mM pNP) with 1X Stop Solution. For example, in a microtube, add 100 μ L of the AP Activity Assay Standard to 900 μ L of 1X Stop Solution, mixing well.
2. Prepare 2-fold serial dilutions of the AP Activity Assay Standard solution with 1X Stop Solution. For example, label ten microtubes #1 to #10, add 0.5 mL of 1X Stop Solution to each tube. Transfer 0.5 mL of the 10-fold diluted AP Assay Standard Solution (0.5 mM final) to tube #1, mix well and transfer 0.5 mL of the mixture to tube #2. Repeat until tube #9, and use tube #10 as blank.
3. Transfer 150 μ L of each dilution, in duplicate, to a 96-well plate, read the absorbance of each well at 405 nm.

Assay Protocol

1. Culture mouse ES cells in medium containing LIF; alternatively, culture human ES cells on a MEF feeder layer.
2. Gently aspirate the medium from the ES cells and wash the cells twice with cold PBS. Aspirate the wash solutions.
3. Lyse the cells in Cell Lysis Buffer (0.5 mL for a 35 mm dish).
4. Incubate for 10 minutes at 4°C, remove the solution and spin down the cell debris at 12,000 X g for 10 minutes. Save the supernatant as cell lysate. Perform a BCA assay or other protein assay to determine the protein concentration of the cell lysate.
5. Add 50 μ L of cell lysate to a 96-well plate. In addition, prepare blank wells that contain 50 μ L Cell Lysis Buffer. We recommend testing samples in triplicate.
6. Initiate the reaction by adding 50 μ L of StemTAG™ AP Activity Assay Substrate. Incubate for 10-30 minutes at 37°C.
7. Stop the reaction by adding 50 μ L of 1X Stop Solution and mix by placing the plate on an orbital plate shaker for 30 seconds.
8. Read the absorbance of each well at 405 nm.

Example of Results

The following figures demonstrate typical results with the StemTAG™ Alkaline Phosphatase Activity Assay Kit. One should use the data below for reference only. This data should not be used to interpret actual results.




Figure 1: AP staining of ES Cells. Murine embryonic stem cells (ES-D3) are maintained in an undifferentiated stage on gelatin-coated dishes in the presence of LIF, as indicated by the high AP activity. To induce differentiation, LIF was withdrawn over a period of several days; various differentiation events were observed (cells became flattened and enlarged with reduced proliferation). At the end of day 5, AP staining of undifferentiated cells was performed with the StemTAG™ Alkaline Phosphatase Staining Kit (Cat # CBA-300).




Figure 2: pNP Standard Curve and AP Activity Assay. **A:** A serial 2-fold dilution of pNP standard was prepared in 1X Stop Solution, and the absorbance of each dilution was measured at 405 nm. **B:** Mouse embryonic D3 cells were grown in the presence or absence of LIF for 5 days. 10 μ g of cell lysate was assayed for AP activity according to the Activity Assay Instructions.

References

1. Wobus AM, Holzhausen H, Jäkel P et al. (1984) *Exp Cell Res* 152:212–219.
2. Thomson JA, Itskovitz-Eldor J, Shapiro SS et al. (1998) *Science* 282:1145–1147.
3. Smith AG, Nichols J, Robertson M et al. (1992) *Dev Biol* 151:339–351.
4. Reubinoff BE, Pera MF, Fong CY et al. (2000) *Nat Biotechnol* 18:399–404.

Recent Product Citations

1. Sanz-Horta, R. et al. (2023). Polycaprolactone with multiscale porosity and patterned surface topography prepared using sacrificial 3D printed moulds: Towards tailor-made scaffolds. *Biomater Adv.* **151**:213465. doi: 10.1016/j.bioadv.2023.213465.
2. Choi, S. et al. (2020). Biochemical activity of magnesium ions on human osteoblast migration. *Biochem Biophys Res Commun.* doi: 10.1016/j.bbrc.2020.07.057.
3. Chang, P.H. et al. (2020). Chitosan 3D cell culture system promotes naïve-like features of human induced pluripotent stem cells: A novel tool to sustain pluripotency and facilitate differentiation. *Biomaterials.* doi: 10.1016/j.biomaterials.2020.120575.
4. Ito, K. et al. (2020). MicroRNA-204 regulates osteogenic induction in dental follicle cells. *J Dent Sci.* doi: 10.1016/j.jds.2019.11.004.
5. Nam, Y.J. et al. (2020). CRH receptor antagonists from *Pulsatilla chinensis* prevent CRH-induced premature catagen transition in human hair follicles. *J Cosmet Dermatol.* doi: 10.1111/jocd.13328.
6. Escobar, A. et al. (2019). Mesoporous Titania Coatings with carboxylated Pores for Complexation and slow Delivery of Strontium for osteogenic Induction. *Appl Surf Sci.* doi: 10.1016/j.apsusc.2019.145172.
7. Escobar, A. et al. (2019). Strontium Titanate (SrTiO₃) Mesoporous Coatings for Enhanced Strontium Delivery and Osseointegration on Bone Implants. *Adv. Eng. Mater.* doi:10.1002/adem.201801210.
8. Li, J. et al. (2019). Osteogenic capacity and cytotherapeutic potential of periodontal ligament cells for periodontal regeneration in vitro and in vivo. *PeerJ.* 7:e6589. doi: 10.7717/peerj.6589.
9. Escobar, A. et al. (2019). Antibacterial Mesoporous Titania Films with Embedded Gentamicin and Surface Modified with Bone Morphogenetic Protein 2 to Promote Osseointegration in Bone Implants. *Advanced Materials Interfaces.* 1801648. doi:10.1002/admi.201801648.
10. Cheng, J. et al. (2019). Stilbene glycoside protects osteoblasts against oxidative damage via Nrf2/HO-1 and NF-κB signaling pathways. *Arch Med Sci.* **15**(1):196-203. doi: 10.5114/aoms.2018.79937.
11. Guo, Y.C. et al. (2018). Ubiquitin-specific protease USP34 controls osteogenic differentiation and bone formation by regulating BMP2 signaling. *EMBO J.* **37**(20). pii: e99398. doi: 10.15252/embj.201899398.
12. Xiong, S. et al. (2018). Immunization with Na⁺/K⁺ ATPase DR peptide prevents bone loss in an ovariectomized rat osteoporosis model. *Biochem Pharmacol.* **156**:281-290. doi: 10.1016/j.bcp.2018.08.024.
13. Camp, E. et al. (2018). miRNA-376c-3p Mediates TWIST-1 Inhibition of Bone Marrow-Derived Stromal Cell Osteogenesis and Can Reduce Aberrant Bone Formation of TWIST-1 Haploinsufficient Calvarial Cells. *Stem Cells Dev.* **27**(23):1621-1633. doi: 10.1089/scd.2018.0083.
14. Yang, F. et al. (2018). Fatty acids modulate the expression levels of key proteins for cholesterol absorption in Caco-2 monolayer. *Lipids Health Dis.* **17**(1):32. doi: 10.1186/s12944-018-0675-y.

15. Abueva, C. D. G. et al. (2018). Multi-channel biphasic calcium phosphate granules as cell carrier capable of supporting osteogenic priming of mesenchymal stem cells. *Materials & Design*. **141**:142–149. doi: 10.1016/j.matdes.2017.12.040.
16. IMAI, K. et al. (2017). Influence of Fluoride Contamination on Titanium Surface on Cell Viability and Cell Differentiation of ES-D3 Cells. *J Oral Tissue Engin*. **15**(1):35-40. doi: 10.11223/jarde.15.35.
17. IMAI, K. et al. (2017). Study of ES Cell Differentiation using Three-dimensional Culture with Silica Fiber. *Nano Biomedicine*. **9**(2):55-60. doi: 10.11344/nano.9.55.
18. Kamiya, N., et al. (2017). Targeted disruption of NF1 in osteocyte increases FGF23 and osteoid with osteomalacia-like bone phenotype. *J Bone Miner Res*. doi: 10.1002/jbmr.3155.
19. Jin, H. et al. (2016). Increased activity of TNAP compensates for reduced adenosine production and promotes ectopic calcification in the genetic disease ACDC. *Sci. Signal*. **9**:ra121.
20. Lee, H. Y. et al. (2016). Porcine placenta hydrolysates enhance osteoblast differentiation through their antioxidant activity and effects on ER stress. *BMC Complement Altern Med*. doi:10.1186/s12906-016-1274-y.
21. Choi, H. Y. et al. (2016). Efficient mRNA delivery with graphene oxide-polyethylenimine for generation of footprint-free human induced pluripotent stem cells. *J Control Release*. **235**:222-235.
22. Pengjam, Y. et al. (2016). Anthraquinone glycoside aloin induces osteogenic initiation of MC3T3-E1 cells: Involvement of MAPK mediated wnt and bmp signaling. *Biomol Ther*. **24**:123-131.
23. Yue, Y. et al. (2015). Safe and bodywide muscle transduction in young adult Duchenne muscular dystrophy dogs with adeno-associated virus. *Hum Mol Genet*. doi:10.1093/hmg/ddv310.
24. Pino-Barrio, M. J. et al. (2015). V-myc immortalizes human neural stem cells in the absence of pluripotency-associated traits. *PLoS One*. **10**:e0118499.
25. Pan, X. et al. (2015). AAV-8 is more efficient than AAV-9 in transducing neonatal dog heart. *Hum Gene Ther Methods*. doi:10.1089/hgtb.2014.128.
26. Salem, O. et al. (2014). Naproxen affects osteogenesis of human mesenchymal stem cells via regulation of Indian hedgehog signaling molecules. *Arthritis Res Ther*. **16**:R152.
27. Guo, L. et al. (2014). Effects of erythropoietin on osteoblast proliferation and function. *Clin Exp Med*. **14**:69-76.
28. Dong, Y. et al. (2014). NOTCH-mediated maintenance and expansion of human bone marrow stromal/stem cells: a technology designed for orthopedic regenerative medicine. *Stem Cells Transl Med*. **3**:1456-1466.
29. Dixon, J.E. et al. (2014). Combined hydrogels that switch human pluripotent stem cells from self-renewal to differentiation. *PNAS* **111**:5580-5585.
30. Moussalli, M.J. et al. (2011). Mechanistic contribution of ubiquitous 15-lipoxygenase-1 expression loss in cancer cells to terminal cell differentiation evasion. *Cancer Prevention Research* **4**:1961-1972.

Warranty

These products are warranted to perform as described in their labeling and in Cell Biolabs literature when used in accordance with their instructions. THERE ARE NO WARRANTIES THAT EXTEND BEYOND THIS EXPRESSED WARRANTY AND CELL BIOLABS DISCLAIMS ANY IMPLIED WARRANTY OF MERCHANTABILITY OR WARRANTY OF FITNESS FOR PARTICULAR PURPOSE. CELL BIOLABS' sole obligation and purchaser's exclusive remedy for breach of this warranty shall be, at the option of CELL BIOLABS, to repair or replace the products. In no event shall CELL BIOLABS be liable for any proximate, incidental or consequential damages in connection with the products.

Contact Information

Cell Biolabs, Inc.
5628 Copley Drive
San Diego, CA 92111
Worldwide: +1 858-271-6500
USA Toll-Free: 1-888-CBL-0505
E-mail: tech@cellbiolabs.com
www.cellbiolabs.com

©2010-2024: Cell Biolabs, Inc. - All rights reserved. No part of these works may be reproduced in any form without permissions in writing.