

## OH-Pen <Lipid Radical Inhibitor>

Catalog NO. FDV-0043

Research use only, not for human or animal therapeutic or diagnostic use.

日本語版はこちらからダウンロードできます。  
①弊社ウェブサイトよりWebページ番号検索にて【70873】で検索  
②QRコードより



### Product Background

**Lipid peroxidation (LPO)** is one of the several degradation processes of lipids under oxidative stress (Figure 1). Primary products in LPO are lipid radicals and there are two major initiators to induce LPO process, pro-oxidants and lipid oxidative enzymes including lipoxygenase (LOX) and cytochrome P450 (CYP). 1) For pro-oxidant-induced LPO, lipids containing unsaturated fatty acid, especially polyunsaturated fatty acids (PUFAs), are attacked by pro-oxidants including reactive oxygen species (ROS) and form lipid-derived radicals. Lipid radical ( $L\cdot$ ) can be easily oxidized to lipid peroxy radical ( $LOO\cdot$ ). Unstable  $LOO\cdot$  immediately extracts a hydrogen from another lipid molecule generating a lipid hydroperoxide (LOOH) and a new lipid radical ( $L\cdot$ ). 2) Another pathway enzyme-induced LPO, lipids containing PUFAs are oxidized to lipid hydroperoxides (LOOH), which decomposes to lipid peroxy radicals  $LOO\cdot$  or alkoxy radicals  $LO\cdot$  by metal ions ( $Fe^{2+}$  etc.). Once lipid radical is produced by the above two processes, lipid radicals expand the radical chain reaction (radical propagation step). In the termination reaction, antioxidants donate a hydrogen atom to the lipid peroxy radical ( $LOO\cdot$ ) species resulting in the formation of many different aldehydes including malondialdehyde (MDA), acrolein, propanal, hexanal, and 4-hydroxynonenal (4-HNE). These aldehydes are cytotoxic because reactive aldehydes attack biomolecules (proteins, DNA/RNA, etc.) to form secondary products. These reactive aldehydes are considered causative factors of organ injury, ferroptosis and ER-stress. To understand the molecular mechanism and physiological relevance of LPO, lipid radical-specific inhibitors are very powerful tools. OH-Pen is a unique lipid radical-specific inhibitor originally developed by Dr. Ken-ichi Yamada, Kyushu University, and will not react with reactive oxygen species. Funakoshi also has lipid radical-specific detector, **LipiRADICAL Green** (#FDV-0042).




Figure 1. Overview of lipid radicals in LPO pathway

## Description

Catalog Number: FDV-0043  
Size: 0.1 mg  
Formulation: C<sub>13</sub>H<sub>26</sub>NO<sub>2</sub> ·  
Chemical structure: See Fig. 2  
Molecular weight: 228.19g/mol  
Solubility: Soluble in DMSO




Figure 2. Chemical structure

## Reconstitution and Storage

Reconstitution: stock solution recommended concentration 1-10 mM in 100% DMSO.

Storage :

Store powder at -20°C.

After reconstitution in DMSO, aliquot and store at -20 °C, avoid repeated freeze-thaw cycles.

## Application data


### Inhibition of nitrosamine-induced carcinogenesis by OH-Pen

Rats received diethylnitrosamine (DEN, 100 mg/kg body weight), which is a well-known hepatic procarcinogen. Subsequently, rats received OH-Pen (2.5 μmol/kg body weight) by intraperitoneal injection after 1 hour DEN administration. For the acute model and chronic model, livers were dissected after 24 hours and 12 weeks DEN administration, respectively.

(Upper panel) Livers from chronic hepatocellular carcinoma model and total foci number

(Middle panel) Quantification of LPO-derived aldehydes in acute model livers.

(Lower panel) Quantification of tissue damage markers. In all panels, OH-Pen clearly suppressed DEN-induced hepatocellular carcinoma.



## Reference

1. Yamada *et al.*, *Nat. Chem. Biol.*, **12**, 608-613 (2016) Fluorescence probes to detect lipid-derived radicals.
2. Matsuoka *et al.*, *Anal. Chem.*, **92**, 6993-7002, (2020) Method for structural determination of lipid-derived radicals

## Related products

### LipiRADICAL Green <Lipid Radical Detection Reagent>

LipiRADICAL Green is a specific fluorescent dye for lipid-derived radicals which are the most upstream factor of lipid peroxidation (LPO). LipiRADICAL Green can be applied into both *in vitro* assay and cell-based assay to monitor lipid radical productions.

Catalog No. FDV-0042


Size 0.1 mg

Features

- Recommended Ex/Em:~480 nm / 540 nm
- Enable to detect very unstable lipid-derived radicals
- Compatible with *in vitro* assay and in cell-based assay
- An innovative reagent for comprehensive identification of lipid-derived radicals by lipidomics


#### Application note 1; Cell-based detection of lipid radicals induced by diethylnitrosamine (DEN)

Hepal-6 cells were treated with 1  $\mu$ M of “LipiRADICAL Green” for 20 min and washed twice with PBS. For inducing an LPO signal, the cells were co-treated with diethylnitrosamine (DEN) and “LipiRADICAL Green”, an LPO initiator. Immediately after DEN addition, the cells were observed by confocal microscopy (Ex.458 nm/ Em. 490-674 nm) for 20 min with 2 min interval. The fluorescent signal of “LipiRADICAL Green” from the DEN-treated cells clearly increased.



#### Application note 2; *in vitro* detection of lipid radicals derived from LDL

Purified low-density lipoprotein (LDL, 20  $\mu$ g protein/mL) was mixed with pro-oxidants hemin or AAPH and “LipiRADICAL Green” and the green fluorescence (Ex. 470 nm/ Em 530 nm) was measured for 60 min at 37°C. Both hemin and AAPH increased green fluorescence indicating the production of lipid radicals from LDL particles in a time-dependent manner.



### AcroleinRED <Cell-based Acrolein Detection Reagent>


Acrolein is a LPO downstream aldehyde and one of the most toxic oxidative stress marker. AcroleinRED is the world first cell-based acrolein detection reagent.

Catalog No. FDV-0022

Size 0.5 mg

Features

- Easy and quick protocol
- Enable to monitor acrolein production under live cells with various stimulations



### CellFluor™ GST <Cell-based GST Activity Assay Reagent >


Glutathione *S*-Transferases (GSTs) are major detoxication enzyme family which neutralize LPO-derived toxic aldehydes. CellFluor™ GST is a novel fluorescent probe for monitoring wide GST members' activity both in cell and *in vitro*. CellFluor™ GST releases green fluorophore rhodamine 110 upon GST activities. This probe has cell-permeability and can detect intracellular GST activity.

Catalog No. FDV-0031

Size 0.1 μmol

Features

- Easy and quick protocol
- Broad specificity for various GST family members
- Ex/Em: 496 nm/520 nm (commercial FITC filters are available)



#### Disclaimer/免責事項

This product has been commercialized by Funakoshi Co., Ltd. based on the results of academic research, and the advertisement text, figures and manuals (hereinafter "Product information") have been prepared based on published research reports on November, 2020. The academic interpretation at the time of creation of the Product Information may change in accordance with future developments in the relevant research field and expansion of various scientific findings, and the latest version and certainty of the Product Information are not guaranteed. The specifications of this product and the Product Information are subject to change without notice. Please contact us for the latest information.

本製品は学術研究成果を基にフナコシ株式会社が製品化したもので、2020年11月時点における公開研究報告を基に広告文章およびマニュアル(以下、製品資料)を作成しています。今後の当該研究分野の発展および各種学術知見の拡大にともない、製品資料作成時の学術的解釈が変更になる可能性があり、最新性・確実性を保証するものではありません。また、本製品の仕様および製品資料を予告なく変更する場合がございます。最新の情報に関しましては、弊社までご確認いただきますようお願い申し上げます。

E-mail Newsletter Sign Up

Japanese English

funakoshi  
FRONTIERS IN LIFE SCIENCE  
URL: <http://funakoshi.co.jp>  
9-7 Hongo 2-Chome, Bunkyo-ku, Tokyo 113-0033