

Operation Manual

Version 1.0

Neo  **Biotech**

NB-12-0035-27

Microplate Reader

PC Software



Thanks for purchasing our Microplate reader. This Operation manual describes PC software functions and how to operate. Please read it carefully before operation and keep this operation manual for later use.

Initial Inspection

Please check the instrument, as well as all included accessories according to the packing list when you first open the packaging, if you find any damaged or missing, please contact distributor or manufacturer.

Contents

Chapter 1	Introduction	1
Chapter 2	Features	2
Chapter 3	Instrument Overview	3
Chapter 4	Software Installation	5
Section 1	Installation Environment.....	5
Section 2	Installation Procedures	5
2.1	Software Installation.....	5
2.2	Environment Installation	7
2.3	DriverSetup Installation	9
Chapter 5	Operation Instruction	11
Section 1	Startup Interface	11
Section 2	Main Interface.....	12
2.1	Title	13
2.2	Menu	13
2.3	Toolbar	13
2.4	Protocol	13
2.5	Status	13
2.6	Parameter info.....	13
2.7	Protocol tree	13
2.8	Control bar.....	14
Section 3	Information	15
3.1	Document Management.....	16
3.2	System settings	22
Section 4	Layout	24
4.1	Layout Parameters.....	24
4.2	Plate layout parameters information	25
Section 5	Step	29
5.1	End-point method.....	29
Section 6	Results.....	31
6.1	Standard Curve	31
6.2	Logit/Log	32
6.3	Spline	34
6.4	Qualitative Analysis.....	35
Section 7	Report	37
7.1	Report exporting.....	37

Chapter 1 Introduction

This microplate reader is an instrument for EIA test, measuring concentration, absorbance, measuring positive or negative reaction between antibody and antigen in samples by reading enzymatic color change – Enzyme Linked Immunosorbent Assay (ELISA).

Key Features and Benefits:

- 1) Easy-to-use controls: 7 inch, color touch screen and 3 external keys.
- 2) Operating system allows acquisition, editing and saving of data.
- 3) Can be used independently from a computer, data can be saved and transferred by USB drive.
- 4) 8 position optical filter wheel; 4 standard optical filters are included, additional optional filters are available.
- 5) 8 channel vertical optical path, zero dispersion single mode fiber measurement system, automatic plate well center position function.
- 6) 96-well visual layout allows easy setting of blank, sample, positive/negative control, quality control and multi-value control.
- 7) Multi-choice tests on single plate.
- 8) Single or double wavelength measurement.
- 9) Self-checking optical path, top reading, mechanical motion.
- 10) Adjustable Plate shaking function, time and speed.
- 11) Quartz Halogen, energy saving light source.
- 12) Multiple ports for data export.
- 13) Measurement results can be exported as .CSV file, compatible with MS Excel.

Chapter 2 Features

Operating conditions:

Ambient operating temperature: 4°C ~ 45°C

Relative humidity: ≤70%



Input power: AC100~240V 50/60Hz 2A

Model Parameter	Microplate reader
Light source	6V, 10W, Quartz-Halogen lamp
Light source est. life	300 hours
Wavelength	400~750nm
Optical filter	Four standard filters included: 405,450,492,630nm. Filter wheel holds up to 8 filters.
Detector Type	Silicon photodetector
Read-out range	0.000-4.000 Abs
Resolution	0.001Abs
Linearity	±1% (0~2.000A) , ±2% (2.000~4.000A)
Precision	CV≤0.2% (0-3Abs); CV≤1% (3-4Abs)
Accuracy (450nm)	±0.005Abs (0 ~2.0 Abs), ±1% (2.0~3.0 Abs), ±1.5% (3.0~4.0 Abs)
Channel deviation	<0.01Abs at 450nm
Measurement Speed	Single wavelength <15s/96well plate, double wavelength <28s/96well plate (Normal speed setting)
Dimension (WXDXH)	295x440x225mm / 11.5x17x9 inches
Weight (kg)	10kg

Chapter 3 Instrument Overview

Before first use of this instrument, please read this chapter carefully.

Structure



There are three input buttons: "Start" "Stop" and "Plate in/out"

USB Ports: For connection of keyboard, mouse and USB drive. (the upper port is for printer connection)

Data port: for PC connection. (contact Neo Biotech or distributor for updated software and features)

Internet port: for PC connection. (contact Neo Biotech or distributor for updated software and features) Rear door panel: can opened to replace light source and access optical filter wheel.

Knob: To lock and unlock the rear door.

Chapter 4 Software Installation

Section 1 Installation Environment

Operation system: Windows7, 64byte, as Fig 1.





Fig 1

Section 2 Installation Procedures

2.1 Software Installation

1. Click icon  AMR-100.exe .
2. Choose installation path way, then click “Next”, as Fig 2.




Fig 2

3. Please choose "Create a desktop shortcut" if needed, as Fig 3.




Fig 3

4. Then "Next" to installation interface as Fig 4 and 5.




Fig 4




Fig 5

2.2 Environment Installation

1. Interface of environment installation will pop out after the above installation finished, choose "I agree to the license terms and conditions" then click "Install" to next step, as Fig 6.




Fig 6

2. Then click "Install" to enter interface as Fig 7 and Fig 8 below.




Fig 7




Fig 8

- Please confirm if the installation is successful in “Control Panel” which path is “control board”→ “Programs”→ “Programs and Features ” as highlighted in Fig 9.




Fig 9

2.3 DriverSetup Installation

- DriverSetup(X64) will pop out once the above installation finished, as Fig 10.




Fig 10

- Click “INSTALL” and then the interface will pop out as Fig 11.



Fig 11

*Attention: This software maybe can not work properly if there is any Chinese character in the path or the path is too long.

Chapter 5 Operation Instruction

Section 1 Startup Interface

Startup interface as below Fig 12.



Fig 12

1. Load and verify required resource files.
2. Animations presentation.
3. Check the connection and acquire basic information.
4. Automatically match the proper port during the connection.
5. Jump to the next interface.

Section 2 Main Interface

Detecting interface as Fig 13 which is composed by 8 parts, please see below table.




Fig 13

No.	Name	Function
①	Title	Including three functions: minimization, maximization and shut off.
②	Menu	Switching among "Shortcut", "Start", "Layout", "Step", "Results" and "Report"
③	Toolbar	Sub-bar of the Menu
④	Protocol	Switching among protocols
⑤	Status	Control connection with Microplate Reader
⑥	Parameter info.	Displays setting information and results
⑦	Protocol tree	Including "Info", "Layout", "Step", "Results" and "Report".
⑧	Control bar	Controlling plate in/out and protocol running

2.1 Title

This part including three functions : minimize, maximum and shut off.

2.2 Menu

Please see below table:

Content	Function
Shortcut	Including "New", "Open", "Save" and "Export".
Start	Including "New", "Open", "Save", "Save as", "Settings", "Language"and "Help".
Layout	Click to Layout interface
Step	Click to Step interface
Results	Click to Result interface
Report	Can export report in excel format.

2.3 Toolbar

For protocols operation, like "Demo1" "Demo2" "Demo3" etc.

2.4 Protocol

Displays opened protocols, Demo1, Demo2 and Demo3 are three default protocols.

2.5 Status

Two buttons included: "Connect/Disconnect" and "Logout", if click "Connect" in the main interface, it will connect to the Microplate Reader.

*** Remark: If connected, both "Disconnect" button in the upper-right corner and  in the lower-left corner will appear, or "Connect" and  will appear when instrument disconnected.**

2.6 Parameter info.

Including setting information and results.

2.7 Protocol tree

Displays the current protocol and newly added steps.

Click protocol tree can switch to selected steps quickly.

2.8 Control bar

Including two buttons “Plate in/out” and “Run”.

Section 3 Information

It will enter into detection interface which including “Info”, “Layout”, “Step”, “Results” and “Report” after the software started, “Info” is the default interface as Fig 14.




Fig 14

Info interface is composed by two parts: shortcut and start, mainly includes “New”, “Open”, “Save”, “Export”, “Path” “System” etc.

3.1 Document Management

3.1.1 New

Click “New” in “Shortcut” or “Start” to create a new protocol as Fig 15 below. Only protocol that without test data is with icon T, protocol that with test data is without icon T.




Fig 15

3.1.2 Open

Click “Open” in “Shortcut” or “Start”, default path is “Pro_Save” file under installation directory as Fig 16.




Fig 16

3.1.3 Save

- Save files by “Save” or “Save as” button, default path is “Pro_Save” file under installation directory

which is also can be changed as Fig 17.




Fig 17

- Click “Save” in “Start” menu, protocol saved in “Pro_Save” file in installation director, then click “OK” to finish. Please note: the saving path can not be changed.




Fig 18

3.1.4 Export

Click “Export” in “Shortcut” menu, default export path is “EXCEL_save” file in installation directory, the path can be changed as Fig 19.



Fig 19

3.1.5 Path

- Click “Path” button to open the path of the last time exporting as Fig 20.




Fig 20

3.2 System settings

3.2.1 Settings

The absorbance and concentration can be set, only for their decimal place(Default value of absorbance is three decimal while two for concentration), as below Fig 21.




Fig 21

3.2.2 Language

Language can be switched between Chinese and English. Restart needed once language changed.




Fig 22

Section 4 Layout

4.1 Layout Parameters

Two options in Layout interface: “New” and “Delete”. For now only one 96-well plate allowed to be newly added.




Fig 23

4.2 Plate layout parameters information

1) Plate layout interface is composed by two parts: "Plate Layout" and "Type"(as Fig 24). "Plate" button is for positions layout while "List" for protocol list and result list, see below Fig 25 and Fig 26.




Fig 24




Fig 25



Fig 26

2) Plat layout setting: in the “Layout” interface, choose the type needed, the button will turn gray(as Fig 27), then click wells in “Play Layout” area. If the type of a well need to be changed, choose clear first, then change the corresponding well. There are 7 types of well: Blank, Standard, QC, Unknown, Negative, Positive and Clear.




Fig 27

3) About plate type:

- Blank: used for blank control, you can set several wells as blank.
- Standard: Standard sample well in dark green are for setting standard curve, you can set several wells with the same number of standard sample, operator can be "+", "-", "×", "÷" and "step by" can be any integer as Fig 28.

Blank	Group:
Standard	Group
QC	Name:
Unknown	Standard
Negative	Replicates:
Positive	1 column * 1 rows
Clear	Concentration:
	Conc.: 1 Step by: +1
	series value
	First: 1
	Operator: +
	Step by: 1
	Unit: Unit

Fig 28

- QC: QC control positions are in dark yellow, the setting procedure is the same as that for standard sample position.
- Unknown: In light blue, several positions can be set as unknown. Click the right place of the "Unknown" button, it can be named freely, the prefix part can accept numbers, letters, even Chinese characters, but for the suffix part, only numbers accepted(except 0). Dilution ratio also can be set, default value is 1:200 as below Fig 29.

Remark: If sample dilution ratio changed, all sample positions dilution ratio in plate layout also will be changed.

- Negative: Negative control position in purple, you can set several positions as negative control.
- Positive: Positive control position in light yellow, also several positions can be set as positive control.
- Clear: positions in gray without samples.

<input type="checkbox"/> Blank	Group:
<input type="checkbox"/> Standard	Group prefix suffix
<input type="checkbox"/> QC	Name: <input type="text"/> 1 <input type="text"/>
<input type="checkbox"/> Unknown	Replicates:
<input type="checkbox"/> Negative	1 column * 1 rows
<input type="checkbox"/> Positive	Dilution:
<input type="checkbox"/> Clear	1:200
	<input type="text" value="1: 200"/>

Fig 29

Section 5 Step

“Step” interface as below:




Fig 30

5.1 End-point method

End-point method including two sections: Wavelength and Shake as Fig 31, Fig 32 below.



Fig 31

Wavelength	Shake Settings:
Shake	<input type="checkbox"/> Shake
	Shake Time: <input type="text" value="00"/> : <input type="text" value="00"/> : <input type="text" value="01"/> HH:MM:SS
	Shake Rate:
	<input checked="" type="radio"/> High <input type="radio"/> Medium <input type="radio"/> Low
	Shake Type:
	<input checked="" type="radio"/> Continue

Fig 32

Note: "WL1" and "WL2" can not be with the same wavelength simultaneously.

Section 6 Results

Result interface as below Fig 33.




Fig 33

6.1 Standard Curve

Click “standard curve” in protocol tree to enter into standard curve interface.




Fig 34

You can switch the fitting type by clicking the radio button as Fig 35 below.




Fig 35

Fitting type and functions as below chart:

Fitting type	Function
Liner	For liner fitting
4PL	4PL fitting
Quadratic	For quadratic polynomial fitting
Cubic	Cubic polynomial fitting
Quartic	Quartic polynomial fitting
Logit/log	Logit/log fitting
Spline	Cubic polynomial fitting between two points

6.2 Logit/Log

Curve fitting conditions:

1. There are must be one sample concentration set to "0", and position of the maximum absorbance point must be corresponding to the "0" concentration point.
***Note: "0" concentration must be the maximum absorbance point, or there will be something wrong with the calculation.**
2. There are must be two or more standard samples with different X values (except 0 concentration point).
3. The curve is linear fitting of logit/log, x-axis stands for linear fitting while y-axis absorbance.

Formula: $\ln \frac{p}{1-p} = k \log_{10}(x) + b$ ($p=y/y_0$, y_0 stands for absorbance value of standard sample at 0 concentration)

4. IC50 is the corresponding concentration value when the absorbance is half of the maximum absorbance of standard sample.

5. For this curve fitting, radio button “Log X” and “Log Y” can not be selected as Fig 36.




Fig 36


6.3 Spline

Standard sample is node which combines with the two adjacent nodes consist the Spline. (Except "0" concentration point)

***Note: "0" concentration must be the maximum absorbance point, or there will be something wrong with the calculation.**

Curve fitting conditions:

1. There are must be one sample concentration set to "0", and position of the maximum absorbance point must be corresponding to the "0" concentration point.
2. IC50 is the corresponding concentration value when the absorbance is half of the maximum absorbance of standard sample.
3. In the curving fitting mode, radio button "Log X" and "Log Y" are in gray as below.



IC50: 10.83954

Name: Cover Save Log X Log Y

CH: Orun

Fitting Type:

Liner 4PL logit/log

Ext. METH Quadratic Spline

Factor Cubic Quartic

Limit: 200

Positive: 100 Low Positive: 0

Fig 37

6.4 Qualitative Analysis

Qualitative analysis is mainly used to verify positive, negative and low positive, only logit/log and Spline fitting type can be used.

Fig 38

Qualitative analysis functions as below. Low positive function is open when the checkbox was chosen, otherwise, it will be closed.

Name	Function
Limit	Limit is determined by input value
Positive	Positive is also determined by input value
Low Positive	Also determined by input value

Click "List" button, please see result interface Fig 39. When sample concentration higher than the positive input value, qualitative result will be positive; when sample concentration lower than the positive input value and higher than low positive input value; qualitative result will be low positive(the checkbox of low positive should be chosen); or it will be negative if in other states.

Plate	Well	Group	Type	Name	Concentration	Dilution	N/P
Platel	B9		Clear		0	1.00	
Platel	B10		Clear		0	1.00	
Platel	B11		Clear		0	1.00	
Platel	B12		Clear		0	1.00	
Platel	C1	Group	Standard	Standard3	6.00000	1.00	
Platel	C2	Group	Unknown	5	> 6.000000	1.00	Positive → Positive
Platel	C3	Group	Unknown	13	4.86371	1.00	→ Negative
Platel	C4		Clear		0	1.00	
Platel	C5		Clear		0	1.00	
Platel	C6		Clear		0	1.00	
Platel	C7		Clear		0	1.00	
Platel	C8		Clear		0	1.00	
Platel	C9		Clear		0	1.00	
Platel	C10		Clear		0	1.00	
Platel	C11		Clear		0	1.00	
Platel	C12		Clear		0	1.00	
Platel	D1	Group	Standard	Standard4	18.00000	1.00	
Platel	D2	Group	Unknown	6	5.22532	1.00	Low Positive → Low Positive
Platel	D3	Group	Unknown	14	5.78200	1.00	Positive → Positive
Platel	D4		Clear		0	1.00	
Platel	D5		Clear		0	1.00	
Platel	D6		Clear		0	1.00	

Fig 39

Section 7 Report

7.1 Report exporting

Click “Excel” in Report interface, default path is “EXCEL_save” under installation directory, also you can change the path as you like as Fig 40.



Fig 40