

Cell

Cell Biology Products

Bio

logy

Distributed by
cliniSciences

Mycoplasma 01

TransDetect[®] Luciferase Mycoplasma Detection Kit
TransDetect[®] PCR Mycoplasma Detection Kit
TransDetect[®] qPCR Mycoplasma Detection Kit
TransSafe[™] Mycoplasma Elimination Reagent (TransMyco Plus) Kit
TransSafe[™] Mycoplasma Prevention Reagent

Cell Isolation and Culture 05

Human Peripheral Blood Lymphocyte Separation Solution
RBC Lysis Buffer (1×)
Recombinant Trypsin-EDTA Solution (1×)
ArtMedia[®] Human T Cell SerumFree Medium (GMP Grade, with/without phenolic red)
ArtMedia[®] Human NK Cell SerumFree Expansion Kit
Human MSC Characterization Kit
TransStem[®] Chemically Defined Xeno-free Cell Cryopreservation Medium-Protein Free
TransStem[®] Chemically Defined Xeno-free Cell Cryopreservation Medium III-DMSO Free, Protein Free
TransStem[®] Natural Killer Cell Cryopreservation Medium—Protein Free

Cell Transfection 16

TransIntro[®] EL Transfection Reagent
TransIntro[®] PL Transfection Reagent
TransIntro[®] mRNA Transfection Reagent

Cell Detection 20

Cell Viability/Proliferation Detection
TransDetect[®] Cell Counting Kit (CCK)
TransDetect[®] Cell LIVE/DEAD Viability/ Cytotoxicity Detection Kit
TransDetect[®] Luminescent Cell Viability Detection Kit
Cell Apoptosis Detection
TransDetect[®] Annexin V-FITC/PI Cell Apoptosis Detection Kit
TransDetect[®] In Situ Fluorescein TUNEL Cell Apoptosis Detection Kit
Luciferase Reporter Assay
TransDetect[®] Single-Luciferase (Firefly) Reporter Assay Kit
TransDetect[®] Single-Luciferase (Renilla) Reporter Assay Kit
TransDetect[®] Double-Luciferase Reporter Assay Kit
TransDetect[®] Dual-Luc Pro Luciferase Reporter Assay Kit
TransDetect[®] Bright-Luc Firefly Luciferase Reporter Assay Kit

GMP-Grade Cytokines 31

Distributed by
clinisciences

Mycoplasma



TransDetect[®] Luciferase Mycoplasma Detection Kit (FM301)

01

- Simple & time-saving operation, performs two luminescence readings and delivers results in 20 minutes.
- High sensitivity, detecting as low as 50 CFU/ml, effectively avoiding false-negative results.
- Excellent inter-well reproducibility for more accurate result determination.

High sensitivity

Sample	Product Result	Company L				TransGen			
		A	B	B/A	Result	A	B	B/A	Result
Negative control		462	371	0.8	Negative	172	113	0.66	Negative
Stock solution		400	50419	126.05	Positive	188	29177	155.2	Positive
8		360	33386	92.74	Positive	191	18283	95.72	Positive
64		403	9978	24.76	Positive	184	4714	25.62	Positive
512		398	5827	14.64	Positive	202	3081	15.25	Positive
4096		485	1543	3.18	Positive	186	776	4.17	Positive
32768		359	366	1.02	Suspected positive	179	153	0.86	Suspected positive
262144		263	184	0.7	Negative	181	107	0.59	Negative

Serially diluted positive samples were detected using reagents from TransGen and Company L. The results demonstrated comparable performance between the two products, with both detecting highly diluted positive samples.

High accuracy

Sample	Product Result	Company L				TransGen			
		A	B	B/A	Result	A	B	B/A	Result
Positive control		198	19613	99.06	Positive	120	12460	103.83	Positive
DMEM medium		296	259	0.88	Negative	102	76	0.75	Negative
1640 medium		423	324	0.77	Negative	141	82	0.58	Negative
Cell culture supernatant 1		1524	379	0.25	Negative	1086	133	0.12	Negative
Cell culture supernatant 2		742	286	0.39	Negative	311	82	0.26	Negative
Cell culture supernatant 3		934	2670	24.72	Positive	89	2643	29.7	Positive
Cell culture supernatant 4		1079	583	0.54	Negative	356	151	0.42	Negative
Cell culture supernatant 5		745	388	0.52	Negative	321	102	0.32	Negative


Across a range of different sample types, TransGen's product demonstrated high consistency with Company L, which is considered the industry gold standard.

TransDetect[®] PCR Mycoplasma Detection Kit (FM311)

02

This product detects the presence of mycoplasma contamination in biological materials such as cultured cells by PCR method. Highly specific primers have been designed to amplify a fragment of 16S rRNA coding DNA that is conserved across all commonly-known mycoplasma species. The kit includes an optimized supermix and primer, ultrapure water and positive control template. Using this kit, cell culture supernatants can be tested directly without DNA extraction. The kit provides a very easy-to-use, simple, fast (within 2 hours), specific and sensitive PCR-based mycoplasma detection method.

- High sensitivity: Able to detect as low as 20 copies of mycoplasma genome.
- High specificity: Only detect mycoplasma DNA, not eukaryotic or bacterial DNA.
- Easy to use: No need to extract genomic DNA, suitable for the detection of a large number of cell samples.
- Positive control and negative control to ensure the accuracy of PCR test results.



Sensitivity of mycoplasma detection by PCR

Distributed by
CliniSciences



TransDetect[®] qPCR Mycoplasma Detection Kit (FM321)

03

- Compliant with the requirements of nucleic acid amplification technology (NAT) for mycoplasma detection as stipulated in the European Pharmacopoeia (EP 2.6.7) and the Japanese Pharmacopoeia (JP G3).
- Wide coverage: Validated for the detection of 13 strains of mycoplasma and plasmid DNA, with database comparisons covering at least 129 species of mycoplasma DNA sequences.
- High sensitivity: Validated mycoplasma strains meet 95% detection requirements, with a detection limit of up to 10 CFU/mL; plasmid standards can reach 1 copy/ μ L.
- High specificity: No detection interference from commonly occurring bacteria or production cell genomes, among others.
- dUTP/UDG contamination prevention system: Effectively prevents contamination of PCR products, ensuring accurate data.
- High stability: Product can undergo at least 15 freeze-thaw cycles without significant impact on performance.
- Broad instrument compatibility: Compatible with mainstream qPCR instruments, ensuring consistent detection results.

High sensitivity, detection limit is 10 CFU/mL


Mycoplasma Species	10 CFU/mL Detection Rate
Mycoplasma gallisepticum	24/24
Mycoplasma synoviae	24/24
Mycoplasma arginini	24/24
Spiroplasma citri	24/24
Mycoplasma fermentans	24/24
Acholeplasma laidlawii	24/24
Mycoplasma hyorhinis	24/24
Mycoplasma orale	24/24
Mycoplasma pneumoniae	24/24
Mycoplasma salivarium	24/24



- 10 CFU/mL Fermentation-derived mycoplasma standards
- NTC, NSC (Negative sample extraction product)
- Positive control
- Internal reference

High specificity, no cross-reaction with non-mycoplasma species

Prokaryotic Species	Eukaryotic Species	Prokaryotic Species
Escherichia coli	Pichia pastoris	Lactobacillus acidophilus
Mycobacterium smegmatis	Chinese Hamster Ovary	Streptococcus pneumoniae
Staphylococcus aureus	Sp2/O cell (mouse)	Streptococcus salivarius
Rhodococcus erythropolis	293T cell	Enterococcus faecalis
Micrococcus luteus	Sf9 cell (Spodoptera frugiperda)	Bacillus subtilis



- Internal reference signal from all samples
- Mycoplasma signal for each non-mycoplasma species
- Mycoplasma signal for positive controls



Distributed by
CliniSciences

TransSafe™ Mycoplasma Elimination Reagent (TransMyco Plus) Kit (FM421)

04

This product contains two types of antibiotic components to effectively eliminate mycoplasma contamination. The first type interferes with protein expression of mycoplasma while the second type interferes with DNA replication of mycoplasma. This product has minimum impact on eukaryotic cells with lower cytotoxicity to maintain high cell viability.

- High mycoplasma removal efficiency
- Active at low working concentration. Low cytotoxicity. Suitable for a wide range of cell types
- Easy to use. Add to the culture medium directly.




TransSafe™ Mycoplasma Prevention Reagent (FM501)

05

This product is a mixture of three mycoplasma antibiotics including macrolides, tetracyclines and fluoroquinolone derivatives, which can effectively interfere with protein expression and DNA replication of mycoplasma. This product can effectively prevent mycoplasma contamination during cell culture process, and can also prevent common gram-positive and gram-negative bacteria. It can be used as a better substitute of common penicillin and streptomycin antibiotics in daily cell culture, with low concentration and almost no cytotoxicity.

- Maximum preventive effect on mycoplasma contamination
- Active at low working concentration. Low cytotoxicity. Suitable for a wide range of cell types.
- Easy to use. Add to the culture medium directly.
- An alternative to Penicillin-Streptomycin with wide range of prevention from common gram-positive and gram-negative bacteria.



Preventive effect on mycoplasma contamination determined by PCR

Distributed by
CliniSciences

Cell Isolation and Culture


B

01


Human Peripheral Blood Lymphocyte Separation Solution (FB102)

- The osmotic molarity of this solution is similar to that of human peripheral blood cells.
- The quantity of separated lymphocytes is greater than 1×10^6 cells/mL whole blood.
- The post-separation viability of lymphocytes > 95%.
- This product is ready-to-use, and the lymphocytes separated under sterile conditions can be used for in vitro culture and immunological assays.

Separation diagram



The cell morphology is well preserved



The yield of cells is greater than 1×10^6 cells/mL whole blood, with a viability of over 95%.

Testing Parameters	Donor 1	Donor 2	Donor3	Donor 4	Donor 5
Number of cells harvested per mL of whole blood	1.51×10^6	1.78×10^6	1.12×10^6	1.33×10^6	1.5×10^6
Lymphocyte viability	98.22%	98.59%	97.08%	98.41%	99.28%

RBC Lysis Buffer (1×) (FB101)

Red Blood Cell (RBC) Lysis Buffer (1×) is designed for the simple and rapid lysis of red blood cells. Due to osmotic stress, it can rapidly rupture non-nucleated red blood cells from blood or tissue samples with minimal effects on lymphocytes or other nucleated cells. This product has been treated with sterile filtration. After processing with RBC Lysis Buffer (1×), the blood or tissue samples can be applied in a variety of routine analysis, tests, primary cell cultures, etc.

- Simple operation with consistent lysis efficiency.
- Broad applicability, suitable for routine analysis, detection, and culture of various cell types and tissues.

02

Distributed by
CliniSciences

Recombinant Trypsin-EDTA Solution (1×) (FG302)

03


- Safe: Xeno-free, no antibiotics, no phenol red.
- Mild reaction: The treatment is mild, resulting in minimal cellular damage, and digestion for 20 minutes has no effect on cell morphology and proliferation rate.
- Convenient: Activity of trypsin can be reduced by dilution without serum termination.

Compared to conventional porcine-derived trypsin – safer and gentler

	Porcine-derived trypsin	Recombinant Trypsin-EDTA Solution (1×)
Animal-derived	Yes	No
Risk of viral contamination	Yes	No
Drug applications	Difficult	Easy
Storage conditions	-20 C	2-8 C , ready-to-use
Digestion termination	Serum, serum-containing culture media, or trypsin inhibitors are required to terminate the digestion process.	The activity can be reduced through dilution with basic culture media or serum-free media, avoiding the need for serum to terminate the reaction.


Gentle digestion with high cell viability and high cell recovery rate

hUC-MSCs were digested for 3 minutes



	TransGen-Porcine-derived trypsin	TransGen-Recombinant trypsin	Company T-Recombinant trypsin
Digestion time	2 min	2-3 min	2-3 min
Viability(%)	93.09	96.1	93.68
Recovered cell count	2.52×10^5	2.83×10^5	2.56×10^5

hUC-MSCs were digested for 20 minutes




	TransGen-Porcine-derived trypsin	TransGen-Recombinant trypsin	Company T-Recombinant trypsin
Digestion time	20 min	20 min	20 min
Viability(%)	83.01	89.01	85.71
Recovered cell count	2.46×10^5	2.45×10^5	2.36×10^5

Distributed by
CliniSciences


Cells exhibited normal morphology and maintained proliferative capacity after either 3- or 20-minute digestion

hUC-MSCs were seeded and cultured for 3 days after 3-minute digestion



	TransGen-Porcine-derived trypsin	TransGen-Recombinant trypsin	Company T-Recombinant trypsin
Digestion time	2 min	2-3 min	2-3 min
Cell proliferation count	8.43×10^5	8.84×10^5	8.85×10^5

hUC-MSCs were seeded and cultured for 3 days after 20-minute digestion



	TransGen-Porcine-derived trypsin	TransGen-Recombinant trypsin	Company T-Recombinant trypsin
Digestion time	20 min	20 min	20 min
Cell proliferation count	5.25×10^5	8.34×10^5	8.65×10^5


ArtMedia[®] Human T Cell Serum-Free Medium (GMP Grade, with/without phenolic red)

04

- Independent research and development
- Serum-free and xeno-free medium
- Suitable for rapid expansion and high-density culture of T-cells
- Can be paired with a variety of T cells activation methods (antibody soluble method, antibody-coated method, magnetic bead conjugate antibody method)
- Production and management with GMP standards


Distributed by
CliniSciences

Good T-cell Morphology




Scale: 100 μm


Supported robust and rapid expansion of T cells, resulting in a 200 to 300-fold increase over a continuous 11- to 13-day culture period.



CD3⁺ T cells comprised >98% of the population, and the proportion of CD4⁺CD8⁺ T cells was sustained.




Strong proliferation capacity, high cell viability



Distributed by
ClinSciences

Enhanced IFN- γ Secretion and Tumor Killing



ArtMedia[®] Human NK Cell Serum-Free Expansion Kit (MK101)

05

1. Safety and Stability

Xeno-Free & Chemically Defined: Contains no animal-derived components and has a fully defined chemical composition.
 Serum-Free: Requires only a small amount of inactivated autologous plasma.
 Cytokine Kit Format: A pure cytokine reagent kit; does not require the use of feeder cells.
 Contains phenol red. Formulated without antibiotics.

2. Powerful Expansion Capability

After 14 days: Total cell count expands 100-250 fold; NK cell count expands 1,000-3,000 fold.
 After 21 days: Total cell count expands 200-500 fold; NK cell count expands 3,000-8,000 fold.

3. High NK Cell Purity

After 14 days: Purity of CD3⁺CD56⁺ cells \geq 90%; Purity of CD3⁺CD56⁺CD16⁺ cells \geq 85%.
 After 21 days: Purity of CD3⁺CD56⁺ cells \geq 95%; Purity of CD3⁺CD56⁺CD16⁺ cells \geq 90%.

4. Potent NK Cell Cytotoxicity


Demonstrates high tumor-killing activity against cancer cell lines.

5. Operational Excellence and Cost-Effectiveness

Simple & Scalable: Easy to use and highly suitable for industrial-scale production.
 High Cost-Effectiveness: Offers superior performance at an optimized cost.
 GMP Compliance: Manufactured and managed in accordance with GMP standards.


20 mL fresh heparinized peripheral blood cultured in a 1 L culture system

Cell Morphology



Scale Bar: 100 μ m


High Cell Viability



High fold expansion

Product Used	Time	Medium Volume (mL)	Cell Number($\times 10^7$ cells)	Total Cell Fold Expansion	NK Cell Fold Expansion
TransGen	D0	10	1.94	1	1
	D5	70	3.49	1.8	/
	D7	200	21	10.82	136.94
	D9	550	53.5	27.58	373.42
	D11	1000	95.8	49.38	684.89
	D14	1000	215.2	110.93	1569.18
Company A	D0	20	1.94	1	1
	D5	80	5.29	2.73	/
	D7	275	23.46	12.09	102.26
	D9	525	59.4	30.62	311.5
	D11	1000	124	63.92	745.86
	D14	1000	203.6	104.95	1303.71

High Purity




30 mL of fresh heparin-anticoagulated peripheral blood cultured in a 2 L system

High cell viability, high fold expansion, high purity

Time	Medium Volume (mL)	Cell Viability(%)	Cell Number($\times 10^7$ cells)	Total Cell Fold Expansion	NK Cell Fold Expansion	CD3 CD56 ⁺	CD3 CD16 ⁺	CD3 CD56 ⁺ CD16 ⁺
D0	20	91.89	3.6	1	1	10.02	10.64	8.97
D5	160	95	10.75	2.99	/	/	/	/
D7	500	95.74	62.61	17.39	146.98	84.68	54.74	51.95
D9	1200	92.83	195.6	54.33	496.48	91.56	87.41	86.25
D11	2000	91.84	408.8	113.56	1085.92	95.82	94.36	91.84
D14	2000	94.13	720	200	1973.85	98.89	96.76	96.15


Distributed by
CliniSciences

Potent and Highly Efficient Killing Activity



30 mL fresh sodium citrate-anticoagulated umbilical cord blood

The 20-day culture yielded ~10 billion highly viable (>90%) and potently cytotoxic cells




14-Day Cultured NK Cell Purity

Frozen PBMCs: 2.12x10⁷ cells in 1 L culture system

High cell viability, high fold expansion, high purity

Time	Medium Volume (mL)	Cell Viability (%)	Cell Number (x10 ⁷ cells)	Total Cell Fold Expansion	NK Cell Fold Expansion	CD3 ⁺ CD56 ⁺	CD3 ⁺ CD16 ⁺	CD3 ⁺ CD56 ⁺ CD16 ⁺
D0	10	88.09	2.12	1	1	4.09	4.12	4.06
D5	70	95.56	4.17	1.97	/	/	/	/
D7	200	98.06	14.26	6.73	72.35	59.17	53.62	50.69
D9	550	98.27	60.64	28.6	389.95	71.91	70.19	70.97
D11	1000	92.16	127	59.91	1277.85	83.79	82.71	81.36
D14	1000	92.86	255	120.28	3018.77	95.31	95	95.02

Potent and Highly Efficient Killing Activity




Distributed by
CliniSciences

Human MSC Characterization Kit (HF001)

- Staining was performed using a multicolor flow cytometry panel containing four antibodies: CD73, CD105, CD90, and CD45.
- The protocol is designed to minimize the number of cells required for the assay.
- The assay is highly specific, user-friendly, and suitable for the identification of MSCs from various tissue sources.

Identification of Mesenchymal Stem Cell Surface Markers



TransStem[®] Chemically Defined Xeno-free Cell Cryopreservation Medium-Protein Free (MC102)

- Safety: Serum-free, xeno-free, protein-free, chemically defined.
- High efficiency: Significantly enhances cell revival rates, with revival rates exceeding 90% for mesenchymal stem cells.
- Stability: Does not affect the phenotype and growth status of stem cells, while maintaining their multipotency.
- Convenience: Ready-to-use product, no need for programmed cooling, can be directly transferred to liquid nitrogen for long-term storage after overnight storage at -80°C.
- Reliability: Suitable for the cryopreservation of various cells, including stem cells, immune cells, tumor cells, etc.

Differences from traditional serum-containing cell cryo-preservation solutions

Differences	Traditional serum-containing cell cryopreservation solution	Serum-free, protein-free cell cryopreservation solution (containing DMSO)
Serum	Yes	No
Protein	Yes	No
Clear chemical composition	Uncertain	Certain
Cooling program	Require	Not required
Cryopreservation solution preparation	Use immediately after preparation	Ready to use, store at 4°C
Revival rate	Normal	Outstanding (support high concentration cryopreservation)
Difference between batches	High	Low
Risk of carrying virus	Yes	No
Applicable to clinical research	No, containing animal serum and protein	Yes, clear chemical composition, no animal-derived substances, protein-free

Distributed by
CliniSciences



TransStem[®] Chemically Defined Xeno-free Cell Cryopreservation Medium III-DMSO Free, Protein Free (MC131) (Completely Using Pharmaceutical Grade Materials)

- Safety: Serum-free, xeno-free, protein-free, DMSO-free, clear chemical component, completely using pharmaceutical grade materials, with complete toxicity test.
- High efficiency: Effectively improves cell recovery viability, with over 90% recovery rate for most cells (including mesenchymal stem cells and immune cells).
- Stability: Does not affect stem cell phenotype and growth status, while maintaining stem cell pluripotency.
- Convenience: Ready-to-use product, no programmed cooling required, can be directly transferred to liquid nitrogen for long-term storage after overnight storage at -80°C.
- Reliability: Suitable for cryopreservation of various cells including stem cells, immune cells, tumor cells, and transformed cell lines.

Differences from traditional serum-containing cell cryopreservation solutions

Differences	Traditional serum-containing cell cryopreservation solution	Serum-free, protein-free, DMSO-free cell cryopreservation solution
Serum	Yes	No
Protein	Yes	No
DMSO	Yes	No
Clear chemical composition	Uncertain	Certain
Cooling program	Require	Not required
Cryo-preservation solution preparation	Use immediately after preparation	Ready to use, store at 4°C
Revival rate	Normal	Outstanding (support high concentration cryo-preservation)
Difference between batches	High	Extremely low
Risk of carrying virus	Yes	No
Apply for pharmaceutical excipients	Apply for pharmaceutical excipients	Apply for pharmaceutical excipients
Applicable to clinical research	No, containing animal serum and protein	Yes, clear chemical composition, xeno-free, protein-free

CD3+T cells cryopreservation




Distributed by
CliniSciences


TransStem[®] Natural Killer Cell Cryopreservation Medium—Protein Free (MC105)

- Safety: Serum-Free, Animal-Component Free, Protein-Free, Chemically Defined, Pharmaceutical-Grade Ingredients, Ultra-Low Endotoxin, < 0.1 EU/mL.
- High efficiency: Effectively increases post-thaw cell viability and supports high-density NK cell cryopreservation. Validated for >85% viability and >80% recovery rate after 3 years of storage.
- Good Stability: minimize post-thaw apoptosis, ensuring sustained cell health and function after resuscitation.
- Ready-to-Use: Compatible with controlled rate freezing protocols and supports direct storage at -80°C.


>80% Viability at 12h Post-Thaw




Normal Surface Marker Phenotype After 1-Year Cryopreservation



>85% Viability After 3 Years



High-Density Cryopreservation Support



Cell Transfection





TransIntro™ EL Transfection Reagent (FT201)

- High transfection efficiency: Improve the transfection efficiency of a variety of difficult-to-transfect cells.
- Low cytotoxicity: Minimum effect on cellular physiological processes.
- Simple operation: No need to change the medium after transfection, tolerance with serum and antibiotics.
- Broad range of applications: Applicable for both DNA and RNA.

01

Successfully transfected cell types with TransIntro™ EL Transfection Reagent

A549	COS-1	Hep G2	MEF	PANC-1
B16-F10	DLD-1	HL-60	MIA PaCa-2	PT67
BHK-21	HCT-116	K562	Neuro-2a	SGC-7901
BTC	HEK-293	L929	NIH/3T3	SH-SY5Y
C2C12CEF	HEK-293T	NRK	P815	Porcine epithelial granulosa cells
CHO	HeLa	Vero	U2OS	Bovine fibroblasts
CEF	COS-7	STO	MARC-145	Pilose antler stem cells
MCF-7				



siRNA transfection effect with HEK-293 cells detected by Western Blot


TransIntro[®] PL Transfection Reagent (FT301)

02

- Highly efficient and stable transfection for tumor cell lines, such as lung cancer, colorectal cancer, liver cancer, pancreatic cancer, facilitating the establishment of stable and reliable cell models.
- Low cytotoxicity and minimum effect on cellular physiological processes.
- Tolerance with antibiotics.

Successfully transfected cell types with TransIntro[®] PL Transfection Reagent

Tumor cell lines			Neural cell lines	Other cell lines		Primary cell lines
A549	HCT 116	MCF-7	Neuro-2a	BHK21	NIH3T3	HUVEC
B16-F10	Huh7	MIA PaCa-2	SH-SY5Y	CHO	Vero	MEF
DLD-1	HT-29	PANC-1		HEK-293	L929	MCM
HeLa	SGC-7901	SMMC-7721		HEK-293T	Marc 145	
Hep G2	SW480			BTC		




Distributed by
CliniSciences

TransIntro[®] mRNA Transfection Reagent (FT501)

03


- High transfection efficiency: Efficiently transfect difficult-to-transfect cells.
- Low cytotoxicity: It has little effect on cell physiological processes.
- Easy operation: There is no need to change culture medium after transfection, and it can be added directly to cells in culture medium with serum and antibiotics.

High transfection efficiency




Hela cells were transfected with 1 kb GFP-mRNA using transfection reagents from TransGen, Company Y, and Company T. At 48 hours post-transfection, cell morphology was observed under a fluorescence microscope, and transfection efficiency was measured using flow cytometry. The results showed that cells transfected with the TransGen reagent exhibited better morphology and higher transfection efficiency.

Transfection with mRNA of different lengths



Hela cells were transfected with 1 kb and 2 kb GFP-mRNA using the TransGen reagent. At 24 hours post-transfection, the cell status remained excellent under fluorescence microscopy, and the transfection efficiency detected by flow cytometry reached over 86%.

High cell viability



Different cell types were transfected with 1 kb GFP-mRNA using reagents from TransGen, Company Y, and Company T. Cell viability was assessed 24 hours post-transfection using a cell counter. The results demonstrated that the TransGen reagent achieved higher cell viability compared to the others.

Cells successfully transfected with TransIntro[®] mRNA Transfection Reagent

A549	CHO	HeLa	Vero	HCT116	MEF	SH-SY5Y	Neuro-2a	Jurkat
C2C12CEF	CEF	MCF-7	HEK-293	HEK-293T	Hep G2	HL-60	NIH/3T3	

Distributed by
CliniSciences

Cell Detection



Cell Viability/Proliferation Detection


01

TransDetect[®] Cell Counting Kit (CCK) (FC101)

Designed for cell proliferation assays as well as cytotoxicity assays by utilizing a water-soluble tetrazolium salt. The faster the cell proliferation, the lower cytotoxicity and the more cells, the darker the color. The color development has a good linear relationship with the number of cells.

- Fast detection and high sensitivity.
- Low cytotoxicity.
- Broad linear range.
- Good stability and high repeatability.

Schematic diagram of CCK




High Sensitivity

HeLa cells were seeded in a 96-well plate at densities of 0, 5,000, 10,000, 15,000, 20,000, and 25,000 cells per well and incubated at 37°C for 3 hours. After removing the cell culture supernatant, 110 µl of culture medium pre-mixed with CCK reagent (100 µl of DMEM with 10% FBS + 10 µl CCK reagent) was added to each well. The cells were then returned to the 37°C incubator for further culture. At 1 hour and 3 hours after adding CCK, the plate was placed in a microplate reader to measure the absorbance at a wavelength of 450 nm. A standard curve was plotted. The results showed a strong linear relationship between the absorbance and the cell number, with $R^2 > 0.99$.


Distributed by
CliniSciences

CCK sensitivity test



Low Cytotoxicity

Based on the above experiment, after the 3-hour measurement, the culture medium in the wells seeded with 10,000 cells was replaced with fresh DMEM medium containing 10% FBS. Within 24 hours after adding CCK, the cell morphology in the culture wells was observed. As shown in the figure, the cells were in good condition, indicating low cytotoxicity of FC101 and that the cells remained viable for subsequent experiments after the assay.




Good Stability

Experimental data demonstrate that the key performance parameters of this CCK product exhibited no significant alterations after three freeze-thaw cycles at -20°C, thereby meeting all predefined quality specifications.

Freeze-Thaw Cycles	R ² Value (Indicator of Linearity)	CV% @ 5000 cells (Indicator of Precision)
0 (Fresh)	0.995	4.5%
1 Cycle	0.993	5.1%
2 Cycles	0.990	6.2%
3 Cycles	0.985	7.8%

Application Example: bFGF Cytokine Activity Assay

NIH-3T3 cells were treated with varying concentrations of the growth factor bFGF (0.001, 0.005, 0.01, 0.1, 0.2, 0.5, 1, 2, 5, 10, 50 ng/ml) for 48 hours, followed by detection with the CCK reagent. The EC₅₀ of bFGF was determined to be 1.350 ng/ml.




Distributed by
CliniSciences

02

TransDetect[®] Cell LIVE/DEAD Viability/ Cytotoxicity Detection Kit (FC301)


- Simple operation, high specificity of fluorescence
- Low cytotoxicity
- Quantification or sorting can be performed by flow cytometry

Detection by Flow Cytometry



The HeLa living cells and the fixed permeabilized dead cells were mixed at a ratio of 4:1, and the results were analyzed by flow cytometry after double staining.

Detection by Fluorescence Microscopy



Fluorescence microscopy results of double stained Hep G2 cells


Distributed by
Clinisciences

03

TransDetect[®] Luminescent Cell Viability Detection Kit (FC401)

- Simple operation: Single-component formulation enables rapid lysis in just 10 minutes, with compatibility for high-throughput detection.
- High sensitivity: Capable of detecting as few as 5 cells.
- Excellent stability: Half-life of up to 3 hours; maintains strong linearity even after 7 days at room temperature.
- Wide linear range: Demonstrates reliable linearity across a broad range (5–100,000 cells).

Schematic diagram of the detection principle






Schematic diagram of luminescent cell viability assay principle Procedure

Procedure




Wide linear range, high sensitivity



Using TransGen's product for cell viability detection across varying cell concentrations, the results demonstrated excellent linearity in fluorescence values within a range of 5–100,000 cells. The assay exhibited a broad linear range and high sensitivity, with a detection limit as low as 5 cells.

Distributed by
Clinisciences


Excellent stability and prolonged half-life




Using TransGen's product to analyze 10,000 cells, fluorescence intensity was measured at different time points. The results demonstrated high stability of the product, maintaining 98% of fluorescence signal intensity within 1 hour and 75% within 2 hours, with a half-life exceeding 3 hours.

Comparison with competitors

Wide linear range, high sensitivity




Excellent stability and prolonged half-life



Cell viability assays were performed using TransGen, Company P, and Company V products across varying cell concentrations. The results demonstrated that TransGen's product exhibited superior linearity ($R^2 > 0.99$) within a broad range of 5-50,000 cells, along with significantly higher fluorescence intensity compared to competitors.

Cell viability was assessed using TransGen, Company P, and Company V products with fluorescence measurements taken at multiple time points. The results demonstrated significantly high stability for TransGen's product.

Application example - cytokine activity analysis



	TransGen	Company P
log(agonist) vs. response -- Variable slope (four parameters)		
Best-fit values		
Bottom	751.7	984.3
Top	5915	6465
LogEC50	0.1055	0.1197
HillSlope	1.904	1.542
EC ₅₀	1.275	1.317
Span	5164	5480


The product's performance was benchmarked against international competitor Company P in detecting recombinant human IL-2 bioactivity. Calculated EC₅₀ values were 1.275 ng/mL (TransGen) and 1.317 ng/mL (Company P).

Cell Apoptosis Detection

TransDetect[®] Annexin V-FITC/PI Cell Apoptosis Detection Kit (FA101)

01

- Detect apoptosis in both adherent and suspension cells using flow cytometry or fluorescence microscopy with Annexin V-FITC as the fluorescent probe.
- Combine Annexin V-FITC with propidium iodide (PI) to distinguish cells at different apoptotic stages.




Flow cytometry analysis of camptothecin-induced apoptosis in Jurkat T-lymphoma cells

TransDetect[®] In Situ Fluorescein TUNEL Cell Apoptosis Detection Kit (FA201)

02

- Low toxicity: Free of highly toxic organic arsenic-based buffers required in traditional TUNEL assays.
- High sensitivity & specificity: Optimized ratio of labeled to unlabeled substrates ensures precise detection.
- Simple one-step protocol: Just mix TdT with Labeling Solution—no complex pipetting or extended incubations.
- Broad compatibility & flexible detection.

High sensitivity, low background staining




Detection result of HeLa cells after 10-minute DNase I treatment


Detection results of HeLa cells after 24-hour camptothecin treatment

Distributed by
Clinisciences

Broad applicability



Good stability



Luciferase Reporter Assay

TransDetect[®] Single-Luciferase (Firefly) Reporter Assay Kit (FR101)

01

TransDetect[®] Single-Luciferase (Renilla) Reporter Assay Kit (FR111)

02


Distributed by
clinisciences

TransDetect[®] Double-Luciferase Reporter Assay Kit (FR201)

03

The activity of firefly luciferase reporter gene was detected with luciferin as the substrate, and the activity of *Renilla* luciferase reporter gene was detected with coelenterin as the substrate.

- Rapid detection
- High sensitivity
- Wide detection range
- No interference with cellular endogenous activity




TransDetect[®] Dual-Luc Pro Luciferase Reporter Assay Kit (FR203)

04

- Stable luminescent signal: a half-life of approximately 2 hours (20-25°C), making it suitable for high-throughput assays.
- Easy to use: No need to remove the culture medium—simply add the detection reagent directly.
- More reliable data: Higher termination efficiency, no background signal interference, and more accurate results.
- Precise detection: Normalized using Renilla luciferase as an internal control gene to correct for well-to-well variations in cell number and transfection efficiency.


High fluorescence intensity with long half-life

For Firefly Luciferase




Distributed by
clinisciences

For Renilla Luciferase



Continuously monitor the luminescence intensity and signal stability of both firefly luciferase and Renilla luciferase for 3 hours.

Firefly luciferase exhibits effective luminescence quenching, delivering more reliable data




Detection of residual signal after quenching firefly luciferase luminescence

Application example

Table 1 Plasmid Dosage for Transfection

96-well plates/well	Fluc plasmid	Rluc plasmid
Control	60 ng	300 ng
Fluc high concentration	120 ng	300 ng
Fluc low concentration	30 ng	300 ng



Comparison between TransGen products and competitor Company P in transfecting 293T cells with varying plasmid amounts.


05

TransDetect[®] Bright-Luc Firefly Luciferase Reporter Assay Kit

(FR104)

- Simple Operation: One-step protocol—no washing or cell harvesting required.
- Fast Detection: Complete all steps in just 5–10 minutes.
- High Stability: Extended half-life of 40 minutes.
- Bright Signal: High fluorescence intensity, ideal for luminescence-sensitive assays.
- Batch processing: Designed for rapid small-batch runs (up to 10 plates).

High fluorescence intensity with excellent signal stability (40-minute half-life)



Comparison of luminescence intensity in 293T cells stably expressing firefly luciferase: TransGen vs. competitor products

Distributed by
CliniSciences


GMP-Grade Cytokines




Recombinant Human IL-2 Protein (GMP Grade) (PM101)

- GMP-grade quality management system.
- Endotoxin level below 10 EU/mg.
- High purity, high activity, and high stability.
- High batch-to-batch consistency.

Purity > 98% (SDS-PAGE & SEC-HPLC)



Biological activity >1×10⁷ IU/mg



Product Name	Cat. No.	Expression Host	Specification	Application		
				Immune Cell Culture	Stem Cell Culture	Organoid Culture
Recombinant Human IL-2 Protein (GMP Grade)	PM101	CHO Cells	1×10 ⁶ IU/100 µg 1×10 ⁷ IU/1 mg	★		
Recombinant Human IL-15 Protein (GMP Grade)	PM102	CHO Cells	10 µg 50 µg 1 mg	★		
Recombinant Human IL-12 Protein (GMP Grade)	PM103	CHO Cells		★		
Recombinant Human FGF basic Protein (GMP Grade)	PM104	<i>E.coli</i>			★	★
Recombinant Human VEGF165 Protein (GMP Grade)	PM105	CHO Cells			★	★
Recombinant Human IL-1 alpha Protein (GMP Grade)	PM106	<i>E.coli</i>		★		
Recombinant Human IL-4 Protein (GMP Grade)	PM107	CHO Cells		★		★
Recombinant Human IL-6 Protein (GMP Grade)	PM108	CHO Cells		★		
Recombinant Human IL-7 Protein (GMP Grade)	PM109	CHO Cells		★		
Recombinant Human IL-18 Protein (GMP Grade)	PM110	CHO Cells		★		
Recombinant Human IL-21 Protein (GMP Grade)	PM111	CHO Cells		★		
Recombinant Human Activin A Protein (GMP Grade)	PM112	CHO Cells			★	★
Recombinant Human TGF-beta 1 Protein (GMP Grade)	PM113	CHO Cells			★	★
Recombinant Human IFN-gamma Protein (GMP Grade)	PM114	CHO Cells		★		
Recombinant Human GM-CSF Protein (GMP Grade)	PM115	CHO Cells		★		
Recombinant Human TNF-alpha Protein (GMP Grade)	PM116	CHO Cells		★		
Recombinant Human EGF Protein (GMP Grade)	PM117	CHO Cells			★	★
Recombinant Human HGF Protein (GMP Grade)	PM118	CHO Cells			★	★
Recombinant Human PDGF-BB Protein (GMP Grade)	PM119	CHO Cells			★	★
Recombinant Human SCF Protein (GMP Grade)	PM120	CHO Cells			★	★
Recombinant Human Vitronectin Protein (GMP Grade)	PM121	CHO Cells			★	★
Recombinant Human Fibronectin Protein (GMP Grade)	PM122	CHO Cells		★	★	★

Distributed by
CliniSciences

Mycoplasma

Product Name	Cat. No.	Specifications
<i>TransDetect</i> [®] Luciferase Mycoplasma Detection Kit	FM301-01	25 rxns
	FM301-02	50 rxns
<i>TransDetect</i> [®] PCR Mycoplasma Detection Kit	FM311-01	100 rxns
<i>TransDetect</i> [®] qPCR Mycoplasma Detection Kit	FM321-00	50 rxns
	FM321-01	100 rxns
<i>TransSafe</i> [™] Mycoplasma Elimination Reagent (TransMyco Plus) Kit	FM421-01	10 ml
	FM421-02	20 ml
<i>TransSafe</i> [™] Mycoplasma Prevention Reagent	FM501-01	1 ml
	FM501-02	5 ml

Cell Isolation and Culture

Product Name	Cat. No.	Specifications
Human Peripheral Blood Lymphocyte Separation Solution	FB102-02	200 ml
RBC Lysis Buffer (1×)	FB101-01	100 ml
Recombinant Trypsin-EDTA Solution (1×)	FG302-01	100 ml
	FG302-02	500 ml
<i>TransSerum</i> [®] Chemically Defined Xeno-free Cell Cryopreservation Medium-Protein Free	MC102-01	100 ml
<i>TransSerum</i> [®] Stem Cell Cryopreservation Medium—Protein Free	MC103-01	100 ml
<i>TransSerum</i> [®] Peripheral Blood Mononuclear Cell Cryopreservation Medium—Protein Free	MC104-01	100 ml
<i>TransSerum</i> [®] Natural Killer Cell Cryopreservation Medium—Protein Free	MC131-01	5 ml
	MC131-02	50 ml
	MC131-03	100 ml
<i>ArtMedia</i> [®] Human T Cell Serum-Free Medium (GMP Grade)	MT101-01	1000 ml
<i>ArtMedia</i> [®] Human T Cell Serum-Free Medium, no phenol red (GMP Grade)	MT102-01	1000 ml
<i>ArtMedia</i> [®] Human NK Cell Serum-Free Expansion Kit	MK101-01	1 L
	MK101-02	2 L
Human MSC Characterization Kit	HF001-01	25 tests
	HF001-02	50 tests

Cell Transfection

Product Name	Cat. No.	Specifications
<i>TransIntro</i> [®] EL Transfection Reagent	FT201-01	0.75 ml
	FT201-02	2×0.75 ml
<i>TransIntro</i> [®] PL Transfection Reagent	FT301-01	0.75 ml
	FT301-02	2×0.75 ml
<i>TransIntro</i> [®] PEI Transfection Reagent (GMP Grade)	FT401-01	1 ml
		10 ml
<i>TransIntro</i> [®] mRNA Transfection Reagent	FT501-01	0.2 ml
	FT501-02	1 ml

Cell Detection


Product Name	Cat. No.	Specifications
<i>TransDetect</i> [®] Cell Counting Kit (CCK)	FC101-01	1 ml
	FC101-02	5 ml
	FC101-03	2×5 ml
	FC101-04	6×5 ml
<i>TransDetect</i> [®] Cell LIVE/DEAD Viability/Cytotoxicity Detection Kit	FC301-01	100 rxns
<i>TransDetect</i> [®] Luminescent Cell Viability Detection Kit	FC401-01	10 ml
	FC401-02	100 ml
<i>TransDetect</i> [®] Annexin V-FITC/PI Cell Apoptosis Detection Kit	FA101-01	25 rxns
	FA101-02	50 rxns
<i>TransDetect</i> [®] <i>In Situ</i> Fluorescein TUNEL Cell Apoptosis Detection Kit	FA201-01	25 rxns
	FA201-02	50 rxns
<i>TransDetect</i> [®] Single-Luciferase (Firefly) Reporter Assay Kit	FR101-01	50 rxns
	FR101-02	200 rxns
<i>TransDetect</i> [®] Single-Luciferase (<i>Renilla</i>) Reporter Assay Kit	FR111-01	50 rxns
	FR111-02	200 rxns
<i>TransDetect</i> [®] Double-Luciferase Reporter Assay Kit	FR201-01	50 rxns
	FR201-02	200 rxns
<i>TransDetect</i> [®] Dual-Luc Pro Luciferase Reporter Assay Kit	FR203-01	100 rxns
	FR203-02	1000 rxns
<i>TransDetect</i> [®] Bio-Luc Firefly Luciferase Reporter Assay Kit	FR103-01	10 ml
	FR103-02	100 ml
	FR103-03	100 ml
<i>TransDetect</i> [®] Bright-Luc Firefly Luciferase Reporter Assay Kit	FR104-01	10 ml
	FR104-02	100 ml

How to navigate through our website?

Homepage: When you first visit our website, you'll land on the homepage. This is where you can find an overview of our services, products, or any important announcements. From here, you can navigate to other sections of the website.



Search Bar: Look for specific Products and equipment



“Remember, if you ever get lost or have questions, don't hesitate to reach out to our customer support team. They're here to help you address any concerns you may have.”

clinisciences

Austria

Company: CliniSciences GmbH
Address: Sternwartestrasse 76, A-1180
Wien - Austria
Telephone: +43 720 115 580
Fax: +43 720 115 577
Email: osterreich@clinisciences.com
Web: <https://www.clinisciences.com>



Belgium

Company: CliniSciences S.R.L
Address: Avenue Stalingrad 52, 1000
Brussels - Belgium
Telephone: +32 2 31 50 800
Fax: +32 2 31 50 801
Email: belgium@clinisciences.com
Web: <https://www.clinisciences.com>



Denmark

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: denmark@clinisciences.com
Web: <https://www.clinisciences.com>



Finland

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: suomi@clinisciences.com
Web: <https://www.clinisciences.com>



France

Company: CliniSciences S.A.S
Address: 74 Rue des Suisses, 92000
Nanterre- France
Telephone: +33 9 77 40 09 09
Fax: +33 9 77 40 10 11
Email: info@clinisciences.com
Web: <https://www.clinisciences.com>



Germany

Company: Clinisciences GmbH
Address: Wilhelm-Mauser-Str. 41-43,
50827 Köln - Germany
Telephone: +49 221 9498 320
Fax: +49 221 9498 325
Email: deutschland@clinisciences.com
Web: <https://www.clinisciences.com>



Iceland

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: island@clinisciences.com
Web: <https://www.clinisciences.com>



Ireland

Company: CliniSciences Limited
Address: Ground Floor, 71 lower Baggot street
Dublin D02 P593 - Ireland
Telephone: +353 1 6971 146
Fax: +353 1 6971 147
Email: ireland@clinisciences.com
Web: <https://www.clinisciences.com>



Italy

Company: CliniSciences S.r.l
Address: Via Maremmana inferiore 378
Roma 00012 Guidonia Montecelio - Italy
Telephone: +39 06 94 80 56 71
Fax: +39 06 94 80 00 21
Email: italia@clinisciences.com
Web: <https://www.clinisciences.com>



Netherlands

Company: CliniSciences B.V.
Address: Kraaijenhoffstraat 137A,
1018RG Amsterdam, - Netherlands
Telephone: +31 85 2082 351
Fax: +31 85 2082 353
Email: nederland@clinisciences.com
Web: <https://www.clinisciences.com>



Norway

Company: CliniSciences AS
Address: c/o MerVerdi Munkerudtunet 10
1164 Oslo - Norway
Telephone: +47 21 988 882
Email: norge@clinisciences.com
Web: <https://www.clinisciences.com>



Poland

Company: CliniSciences sp.Z.o.o.
Address: ul. Rotmistrza Witolda Pileckiego 67
lok. 200 - 02-781 Warszawa -Poland
Telephone: +48 22 307 0535
Fax: +48 22 307 0532
Email: polska@clinisciences.com
Web: <https://www.clinisciences.com>



Portugal

Company: Quimigen Unipessoal LDA
Address: Rua Almada Negreiros, Lote 5, Loja 14,
2615-275 Alverca Do Ribatejo - Portugal
Telephone: +351 30 8808 050
Fax: +351 30 8808 052
Email: info@quimigen.com
Web: <https://www.quimigen.pt>



Spain

Company: CliniSciences Lab Solutions
Address: C/ Hermanos del Moral 13
(Bajo E), 28019, Madrid - Spain
Telephone: +34 916 750 700
Fax: +34 91 269 40 74
Email: espana@clinisciences.com
Web: <https://www.clinisciences.com>



Sweden

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: sverige@clinisciences.com
Web: <https://www.clinisciences.com>



Switzerland

Company: CliniSciences AG
Address: Fracht Ost Flughafen Kloten
CH-8058 Zürich - Switzerland
Telephone: +41 (044) 805 76 81
Fax: +41 (044) 805 76 75
Email: switzerland@clinisciences.com
Web: <https://www.clinisciences.com>



UK

Company: CliniSciences Limited
Address: 11 Progress Business center, Whittle
Parkway, SL1 6DQ Slough- United Kingdom
Telephone: +44 (0)1753 866 511
or +44 (0) 330 684 0982
Fax: +44 (0)1753 208 899
Email: uk@clinisciences.com
Web: <https://www.clinisciences.com>



USA

Company: CliniSciences LLC
Address: c/o Carr Riggs Ingram,
500 Grand Boulevard, Suite 210 Miramar
Beach, FL 32550- USA
Telephone: +1 850 650 7790
Fax: +1 850 650 4383
Email: usa@clinisciences.com
Web: <https://www.clinisciences.com>

