

# MYCO

Mycoplasma  
Detection Products

# PLASMA


# DETECTION

Distributed by  
**cliniSciences**

# TransDetect<sup>®</sup> Luciferase Mycoplasma Detection Kit

01

This product is designed using an enzyme unique to mycoplasma. This enzyme degrades MycoDetect Substrate and converts ADP to ATP. Luciferase catalyzes the oxidation of luciferin to produce bioluminescence in the presence of ATP. Cell cultures can be characterized for the presence of mycoplasma contamination by detecting bioluminescence. As this assay can only detect bioactive mycoplasma, the results will be more accurate than that of the PCR assay. The method features high sensitivity, simple operation and time saving.



	TransGen	Company L
Original sample	+	+
1:5 dilution	+	+
1:25 dilution	+	+
1:125 dilution	+	+/-
1:625 dilution	+/-	-
1:3125 dilution	-	-
Negative control	-	-


+ Sample contaminated with mycoplasma  
 +/- Need to check and verify  
 - Samples not contaminated with mycoplasma

# TransDetect<sup>®</sup> PCR Mycoplasma Detection Kit

02

This product detects the presence of mycoplasma contamination in biological materials such as cultured cells by PCR method. Highly specific primers have been designed to amplify a fragment of 16S rRNA coding DNA that is conserved across all commonly-known mycoplasma species. The kit includes an optimized supermix and primer, ultrapure water and positive control template. Using this kit, cell culture supernatants can be tested directly without DNA extraction. The kit provides a very easy to use, simple, fast (within 2 hours), specific and sensitive PCR-based mycoplasma detection method.

- High sensitivity: Able to detect as low as 20 copies of mycoplasma genome.
- High specificity: Only detect mycoplasma DNA, not eukaryotic and bacterial DNA.
- Easy to use, no need to extract genomic DNA, suitable for the detection of a large number of cell samples.
- Positive control and negative control to ensure the accuracy of PCR test results.



Mycoplasma PCR detection sensitivity

Distributed by  
**CliniSciences**



# TransDetect<sup>®</sup> qPCR Mycoplasma Detection Kit

# 03

- Compliant with the requirements of nucleic acid amplification technology (NAT) for mycoplasma detection as stipulated in the European Pharmacopoeia (EP 2.6.7) and the Japanese Pharmacopoeia (JP G3).
- Wide Coverage: Validated for the detection of 13 strains of mycoplasma and plasmid DNA, with database comparisons covering at least 129 species of mycoplasma DNA sequences.
- High Sensitivity: Validated mycoplasma strains meet 95% detection requirements, with a detection limit of up to 10 CFU/mL; plasmid standards can reach 1 copy/μL.
- High Specificity: No detection interference from commonly occurring bacteria or production cell genomes, among others.
- dUTP/UDG Contamination Prevention System: Effectively prevents contamination of PCR products, ensuring accurate data.
- High Stability: Product can undergo at least 15 freeze-thaw cycles without significant impact on performance.
- Broad Instrument Compatibility: Compatible with mainstream qPCR instruments, ensuring consistent detection results.

## High sensitivity, detection limit is 10 CFU/mL

Mycoplasma species	10 CFU/mL detection rate
Mycoplasma gallisepticum	24/24
Mycoplasma synoviae	24/24
Mycoplasma arginini	24/24
Spiroplasma citri	24/24
Mycoplasma fermentans	24/24
Acholeplasma laidlawii	24/24
Mycoplasma hyorhinis	24/24
Mycoplasma orale	24/24
Mycoplasma pneumoniae	24/24
Mycoplasma salivarium	24/24




- 10 CFU/mL Fermentation-derived mycoplasma standards
- NTC, NSC (Negative sample extraction product)
- Positive control
- Internal reference

Distributed by  
**CliniSciences**

**High specificity, no cross-reaction with non-mycoplasma species**

Prokaryotic Species	Prokaryotic Species	eukaryotic Species
Lactobacillus acidophilus	Escherichia coli	Pichia pastoris
Streptococcus pneumoniae	mycobacterium smegmatis	Chinese Hamster Ovary
Streptococcus salivarius	Staphylococcus aureus	Sp2/O call (rat)
Enterococcus faecalis	Rhodococcus erythropolis	293T cell
Bacillus subtilis	Micrococcus luteus	Sf9 cell (Spodoptera frugiperda)





# TransSafe™ Mycoplasma Elimination Reagent (TransMyco Plus) Kit

This product contains two types of antibiotic components, which can inhibit DNA replication and interfere with protein expression of mycoplasma, respectively, so that it can be used to eliminate mycoplasma contamination effectively. It exhibits minimal cytotoxicity on eukaryotic cells and can maintain relatively high cell viability.

- High mycoplasma removal efficiency.
- Active at low working concentration. Low cytotoxicity. Suitable for a wide range of cell types.
- Easy to use. Add to the culture medium directly.

04




Distributed by  
**CliniSciences**

# TransSafe™ Mycoplasma Prevention Reagent

This product is a mixture of three mycoplasma antibiotics including macrolides, tetracyclines and fluoroquinolone derivatives, which can effectively interfere with protein expression and DNA replication of mycoplasma. This product can effectively prevent mycoplasma contamination during cell culture process, and can also prevent common gram-positive and gram-negative bacteria. It can be used as a better substitute for common penicillin and streptomycin antibiotics in daily cell culture, with low concentration and almost no cytotoxicity.

- Maximum preventive effect on mycoplasma contamination .
- Low cytotoxicity. Suitable for a wide range of cell types.
- Easy to use. Add to the culture medium directly.
- An alternative to Penicillin-Streptomycin with wide range of prevention from common gram-positive and gram-negative bacteria.

# 05



PCR detected mycoplasma prevention

Product Name	Cat.No.	Specification
TransDetect® Luciferase Mycoplasma Detection Kit	FM301-01	25 rxns
	FM301-02	50 rxns
TransDetect® qPCR Mycoplasma Detection Kit	FM321-00	50 rxns
	FM321-01	100 rxns
TransDetect® PCR Mycoplasma Detection Kit	FM311-01	100 rxns
TransSafe™ Mycoplasma Elimination Reagent (TransMyco Plus) Kit	FM421-01	10 ml
	FM421-02	20 ml
TransSafe™ Mycoplasma Prevention Reagent	FM501-01	1 ml
	FM501-02	5 ml

Distributed by  
**CliniSciences**

# How to navigate through our website?

**Homepage:** When you first visit our website, you'll land on the homepage. This is where you can find an overview of our services, products, or any important announcements. From here, you can navigate to other sections of the website.

Callouts on the homepage:

- Consult our catalog and services
- Choose your Country
- Log in to your account
- Choose your language
- View Cart
- Search for a product by description, reference, CAS, Number, target...

The homepage features the Clinisciences logo, a search bar, a navigation menu (Catalog, Research areas, Diagnostic, Custom services, CRO services, Technical resources, Special offer, Brand names, Contact), and a featured product section for 'upcyte® cells Expanded primary cells'. Below this is a 'Categories' section with various product categories like Cell-based assay kits, Culture flasks, etc.

## Search Bar: Look for specific Products and equipment

Search Bar: Search for a product (description, reference, CAS number, target ...)

Search result : 39126007 product found

Search by selecting either RUO or IVD products

Refine your search by choosing from the following options

Apply Filter

Cat#	Description	Cond.
NB-45-00042-100	Super Ni-NTA Agarose Resin	100ml
NB-45-00042-25	Super Ni-NTA Agarose Resin	25ml
NB-45-00058-4	Proteus 1 -step Batch Mini Spin Column Pack	40pc
NB-12-6001-3	NeoLine pipette 2-20 µl	1unit
NB-12-0023C	Mini Centrifuge N500C @10,000rpm (Including 6x1.5/2.0ml angle rotor)	1pcs
NB-03-0160	Proteinase K (Powder)	100mg
NB-60-0001	NeoPrep mini	50columns
NB-12-8001-19	Combs for NeoPRO4 mini (1.5mm, 15 wells)	5pieces
NB-12-8001-20	Spacer glasses flat for NeoPRO4 mini (0.75mm, 100*83mm)	5pieces
NB-12-8001-04	Short glasses flat for NeoPRO4 mini (1.0mm, 100*73mm)	10pieces

***“Remember, if you ever get lost or have questions, don't hesitate to reach out to our customer support team. They're here to help you address any concerns you may have.”***

Distributed by  
**Clinisciences**

Distributed by  
**CliniSciences**

# clinisciences

## Austria

Company: CliniSciences GmbH  
Address: Sternwartestrasse 76, A-1180  
Wien - Austria  
Telephone: +43 720 115 580  
Fax: +43 720 115 577  
Email: [osterreich@clinisciences.com](mailto:osterreich@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Belgium

Company: CliniSciences S.R.L  
Address: Avenue Stalingrad 52, 1000  
Brussels - Belgium  
Telephone: +32 2 31 50 800  
Fax: +32 2 31 50 801  
Email: [belgium@clinisciences.com](mailto:belgium@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Denmark

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [denmark@clinisciences.com](mailto:denmark@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Finland

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [suomi@clinisciences.com](mailto:suomi@clinisciences.com)  
Web: <https://www.clinisciences.com>



## France

Company: CliniSciences S.A.S  
Address: 74 Rue des Suisses, 92000  
Nanterre- France  
Telephone: +33 9 77 40 09 09  
Fax: +33 9 77 40 10 11  
Email: [info@clinisciences.com](mailto:info@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Germany

Company: Clinisciences GmbH  
Address: Wilhelm-Mauser-Str. 41-43,  
50827 Köln - Germany  
Telephone: +49 221 9498 320  
Fax: +49 221 9498 325  
Email: [deutschland@clinisciences.com](mailto:deutschland@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Iceland

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [island@clinisciences.com](mailto:island@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Ireland

Company: CliniSciences Limited  
Address: Ground Floor, 71 lower Baggot street  
Dublin D02 P593 - Ireland  
Telephone: +353 1 6971 146  
Fax: +353 1 6971 147  
Email: [ireland@clinisciences.com](mailto:ireland@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Italy

Company: CliniSciences S.r.l  
Address: Via Maremmana inferiore 378  
Roma 00012 Guidonia Montecelio - Italy  
Telephone: +39 06 94 80 56 71  
Fax: +39 06 94 80 00 21  
Email: [italia@clinisciences.com](mailto:italia@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Netherlands

Company: CliniSciences B.V.  
Address: Kraaijenhoffstraat 137A,  
1018RG Amsterdam, - Netherlands  
Telephone: +31 85 2082 351  
Fax: +31 85 2082 353  
Email: [nederland@clinisciences.com](mailto:nederland@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Norway

Company: CliniSciences AS  
Address: c/o MerVerdi Munkerudtunet 10  
1164 Oslo - Norway  
Telephone: +47 21 988 882  
Email: [norge@clinisciences.com](mailto:norge@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Poland

Company: CliniSciences sp.Z.o.o.  
Address: ul. Rotmistrza Witolda Pileckiego 67  
lok. 200 - 02-781 Warszawa -Poland  
Telephone: +48 22 307 0535  
Fax: +48 22 307 0532  
Email: [polska@clinisciences.com](mailto:polska@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Portugal

Company: Quimigen Unipessoal LDA  
Address: Rua Almada Negreiros, Lote 5, Loja 14,  
2615-275 Alverca Do Ribatejo - Portugal  
Telephone: +351 30 8808 050  
Fax: +351 30 8808 052  
Email: [info@quimigen.com](mailto:info@quimigen.com)  
Web: <https://www.quimigen.pt>



## Spain

Company: CliniSciences Lab Solutions  
Address: C/ Hermanos del Moral 13  
(Bajo E), 28019, Madrid - Spain  
Telephone: +34 916 750 700  
Fax: +34 91 269 40 74  
Email: [espana@clinisciences.com](mailto:espana@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Sweden

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [sverige@clinisciences.com](mailto:sverige@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Switzerland

Company: CliniSciences AG  
Address: Fracht Ost Flughafen Kloten  
CH-8058 Zürich - Switzerland  
Telephone: +41 (044) 805 76 81  
Fax: +41 (044) 805 76 75  
Email: [switzerland@clinisciences.com](mailto:switzerland@clinisciences.com)  
Web: <https://www.clinisciences.com>



## UK

Company: CliniSciences Limited  
Address: 11 Progress Business center, Whittle  
Parkway, SL1 6DQ Slough- United Kingdom  
Telephone: +44 (0)1753 866 511  
or +44 (0) 330 684 0982  
Fax: +44 (0)1753 208 899  
Email: [uk@clinisciences.com](mailto:uk@clinisciences.com)  
Web: <https://www.clinisciences.com>



## USA

Company: CliniSciences LLC  
Address: c/o Carr Riggs Ingram,  
500 Grand Boulevard, Suite 210 Miramar  
Beach, FL 32550- USA  
Telephone: +1 850 650 7790  
Fax: +1 850 650 4383  
Email: [usa@clinisciences.com](mailto:usa@clinisciences.com)  
Web: <https://www.clinisciences.com>

