



Jotbody EV Platform: Targeted Nucleic Acid Delivery via Extracellular Vesicles in Metastatic Breast Cancer

Distributed by
cliniSciences



Introduction

Extracellular vesicles (EVs) are naturally released nanoparticles which act as messenger in cell to cell communication. EVs deliver functional biomolecules from donor to recipient cells. With their natural biocompatibility, tissue-homing capabilities, and ability to transport functional molecules, EVs are an efficient and adaptable platform for drug delivery. However, naïve EVs suffer from non-specific uptake, which limits their therapeutic potential.

To transform EVs into a smart drug delivery system, We developed a platform for EVs' nucleic loading and surface modification, empowering EVs with the ability to deliver therapeutic nucleic acid and other molecules into targeted cells. As a proof of concept, antisense oligonucleotide (ASO) anti oncogene miR-125 was loaded into red blood cells derived EVs (RBCEVs) and a single domain antibody (sdAb) against EGFR was conjugated onto RBCEVs surface. The treatment successfully suppressed tumor growth in a mouse model of metastatic breast cancer, demonstrating the potential of Jotbody's platform for enhancing targeted nucleic acid-based therapies using EVs.

- **Targeted Delivery:** Jotbody's platform allows for conjugation of proteins, antibodies, antibody fragments, and single-domain antibodies.
- **Nucleic Acid Loading:** Jotbody's platform allows for effective loading of nucleic acids, including small RNAs and plasmids.
- **Enhanced Therapeutic Efficacy:** Enables precise and efficient delivery of genetic therapies directly to target cells, improving therapeutic outcomes while minimizing off target effects.
- **Safety and Biocompatibility:** Extensive testing confirms that our engineered EVs can be administered repeatedly without adverse effects.
- **Range of application:** The platform is dynamic and translatable to EVs of any source.

Jotbody's Smart Drug Delivery Platform



Methods

Every section of method follows the latest guidelines suggested by Minimal information for studies of extracellular vesicles (MISEV).

EV Isolation/Characterization

- Isolated RBCEVs from healthy human donors.
- Ultracentrifugation, Size exclusion chromatography (SEC), Particle Metrix ZetaView, Transmission electron microscopy (TEM), Western blotting for protein detection.

EV Functionalization

- ASO or siRNA were transfected into RBCEV using Jotbody platform.
- Alpacas derived sdAbs (VHH) were covalently conjugated onto RBCEVs surface using Jotbody platform.

Results

Nanobody conjugated EVs show increased accumulation in target cells




Figure A&B

VHH conjugation did not affect the naïve EVs structure, shape and size as confirmed by TEM and NTA.

Figure C

EGFR VHH was determinant to significantly increased the uptake of EVs into CA1a cells (EGFR+) but not into control cells MOLM13 (EGFR -) as evaluated by flow cytometry analysis.

Targeted EVs specifically deliver the therapeutic cargo into target cells

Nucleic Acid Targeted Delivery




Figure A

Surface modification with EGFR nanobody specifically drives FAM-ASO loaded EVs into EGFR-positive CA1a cells.

Figure B

Surface modification with EGFR nanobody specifically drives siRNA loaded EVs into target cells resulting in reduction of GFP expression.

C




Figure C

EGFR-targeted EVs demonstrated the ability to efficiently deliver siRNAs, resulting in a 70% reduction in GFP signal.

In Vitro Therapeutic Effects




Figure D

Comparison of Ct values revealed that each EV contained on average 450 miR-125b ASOs.




Figure E




Figure E

EGFR-targeted EVs loaded with ASO (125b) significantly decrease the expression of the oncogenic miRNA 125 in both CA1a and 4T1 breast cancer cells.

Figure F

The decrease in miRNA 125 expression led to a reduction in cell viability in both CA1a and 4T1 cells.

Functionalized EVs suppressed th tumor growth in mouse model of breast metastatic cancer

In vivo Targeting Activity




Figure A

We created an hEGFR-positive lung metastatic breast cancer model by systemically injecting 4T1-tdTomato-hEGFR cells into NSG-SGM3 mice for intratracheal EVs treatment.

Figure B

Upon intratracheal administration, EGFR-targeted EVs showed significantly higher levels of uptake by tumor cells as compared to untargeted EVs.

In vivo Therapeutic Effect




Figure C

EGFR-targeted miR-125b ASO-loaded EVs induce a significantly greater knockdown in tumor cells compared to healthy lung cells.

Figure D

EGFR-VHH-EVs loaded with miR-b ASO exhibited superior tumor-suppressive effects, leading to a significantly reduced percentage of tumor cells in the lungs compared to control groups.

Figure E

EGFR-VHH-EVs loaded with miR-125b ASO showed superior tumor suppression, significantly reducing lung tumor area compared to the control groups.




Figure F

Tumor nodules in the EGFR-targeted miR-125b-EVs treatment were significantly smaller and did not show high levels of colocalization with CD31, a well-established marker for blood vessel density in tumor tissue.

Conclusion

This modular and biocompatible approach to EV engineering presents an innovative drug delivery platform for the targeted delivery of diverse therapeutics. The method allows for stable conjugation of targeting moieties and therapeutic payloads, improving specificity, efficiency, and stability. The engineered EVs show significant potential in reducing side effects, overcoming cellular barriers, and delivering therapeutic molecules directly to diseased cells. This breakthrough paves the way for safe, scalable, and versatile EV-based drug delivery solutions, setting a new benchmark in nanomedicine.

Reference

Surface-engineered extracellular vesicles for targeted delivery of therapeutic RNAs and peptides for cancer therapy", *Theranostics* (2022); 12(7):3288-3315

Other related publication to learn further about Jotbody

Extracellular Vesicle Surface Display Enhances the Therapeutic Efficacy and Safety Profile of Cancer Immunotherapy, *Molecular Therapy* (2024) 32: 1-22

Red Blood Cell-Derived Extracellular Vesicles Display Endogenous Antiviral Effects and Enhance the Efficacy of Antiviral Oligonucleotide Therapy, *ACS Nano* 2023, 17, 21, 21639–21661

Covalent Conjugation of Extracellular Vesicles with Peptides and Nanobodies for Targeted Therapeutic Delivery", *Journal of Extracellular Vesicles*, 10.4 (2021): e12057

clinisciences

Austria

Company: CliniSciences GmbH
Address: Sternwartestrasse 76, A-1180
Wien - Austria
Telephone: +43 720 115 580
Fax: +43 720 115 577
Email: oesterreich@clinisciences.com
Web: <https://www.clinisciences.com>



Belgium

Company: CliniSciences S.R.L
Address: Avenue Stalingrad 52, 1000
Brussels - Belgium
Telephone: +32 2 31 50 800
Fax: +32 2 31 50 801
Email: belgium@clinisciences.com
Web: <https://www.clinisciences.com>



Denmark

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: denmark@clinisciences.com
Web: <https://www.clinisciences.com>



Finland

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: suomi@clinisciences.com
Web: <https://www.clinisciences.com>



France

Company: CliniSciences S.A.S
Address: 74 Rue des Suisses, 92000
Nanterre- France
Telephone: +33 9 77 40 09 09
Fax: +33 9 77 40 10 11
Email: info@clinisciences.com
Web: <https://www.clinisciences.com>



Germany

Company: Clinisciences GmbH
Address: Wilhelm-Mausier-Str. 41-43,
50827 Köln - Germany
Telephone: +49 221 9498 320
Fax: +49 221 9498 325
Email: deutschland@clinisciences.com
Web: <https://www.clinisciences.com>



Iceland

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: island@clinisciences.com
Web: <https://www.clinisciences.com>



Ireland

Company: CliniSciences Limited
Address: Ground Floor, 71 lower Baggot street
Dublin D02 P593 - Ireland
Telephone: +353 1 6971 146
Fax: +353 1 6971 147
Email: ireland@clinisciences.com
Web: <https://www.clinisciences.com>



Italy

Company: CliniSciences S.r.l
Address: Via Maremmana inferiore 378
Roma 00012 Guidonia Montecelio - Italy
Telephone: +39 06 94 80 56 71
Fax: +39 06 94 80 00 21
Email: italia@clinisciences.com
Web: <https://www.clinisciences.com>



Netherlands

Company: CliniSciences B.V.
Address: Kraaijenhoffstraat 137A,
1018RG Amsterdam, - Netherlands
Telephone: +31 85 2082 351
Fax: +31 85 2082 353
Email: nederland@clinisciences.com
Web: <https://www.clinisciences.com>



Norway

Company: CliniSciences AS
Address: c/o MerVerdi Munkerudtunet 10
1164 Oslo - Norway
Telephone: +47 21 988 882
Email: norge@clinisciences.com
Web: <https://www.clinisciences.com>



Poland

Company: CliniSciences sp.Z.o.o.
Address: ul. Rotmistrza Witolda Pileckiego 67
lok. 200 - 02-781 Warszawa -Poland
Telephone: +48 22 307 0535
Fax: +48 22 307 0532
Email: polska@clinisciences.com
Web: <https://www.clinisciences.com>



Portugal

Company: Quimigen Unipessoal LDA
Address: Rua Almada Negreiros, Lote 5, Loja 14,
2615-275 Alverca Do Ribatejo - Portugal
Telephone: +351 30 8808 050
Fax: +351 30 8808 052
Email: info@quimigen.com
Web: <https://www.quimigen.pt>



Spain

Company: CliniSciences Lab Solutions
Address: C/ Hermanos del Moral 13
(Bajo E), 28019, Madrid - Spain
Telephone: +34 916 750 700
Fax: +34 91 269 40 74
Email: espana@clinisciences.com
Web: <https://www.clinisciences.com>



Sweden

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: sverige@clinisciences.com
Web: <https://www.clinisciences.com>



Switzerland

Company: CliniSciences AG
Address: Fracht Ost Flughafen Kloten
CH-8058 Zürich - Switzerland
Telephone: +41 (044) 805 76 81
Fax: +41 (044) 805 76 75
Email: switzerland@clinisciences.com
Web: <https://www.clinisciences.com>



UK

Company: CliniSciences Limited
Address: 11 Progress Business center, Whittle
Parkway, SL1 6DQ Slough- United Kingdom
Telephone: +44 (0)1753 866 511
or +44 (0) 330 684 0982
Fax: +44 (0)1753 208 899
Email: uk@clinisciences.com
Web: <https://www.clinisciences.com>



USA

Company: CliniSciences LLC
Address: c/o Carr Riggs Ingram,
500 Grand Boulevard, Suite 210 Miramar
Beach, FL 32550- USA
Telephone: +1 850 650 7790
Fax: +1 850 650 4383
Email: usa@clinisciences.com
Web: <https://www.clinisciences.com>

