



Empowering Life Science with Precision

Jotbody provides professional service for Extracellular Vesicles

1. About us

Established in 2020, Jotbody (HK) Limited is a state-of-the-art biotechnology company specializing in extracellular vesicles (EVs) and nanobody research based in Hong Kong Science and Technology Park.

Our team, comprising experienced scientists and dedicated board members, is committed to advancing global research and development. Since 2014, we have been deeply involved in EV research, collaborating with industry and academic partners. Furthermore, we were honored with a place on Forbes ASIA 100 to watch in 2023.

Our expertise lies in EV functionalization for therapeutic applications, which has been proven across various disease models, including a number of pre-clinical models of human diseases like cancer, neurodegenerative disorders, cardiovascular diseases, and inflammatory conditions.

Within the realm of EVs, we offer comprehensive services tailored to meet the evolving needs of researchers worldwide. Our services encompass:

- Isolation.
- Purification.
- Characterization.
- Omics analysis.
- Functionalization.
- CDMO solutions.

Our vision is clear: to elevate life science with precision through the incredible potential of extracellular vesicles. We are dedicated to fostering collaboration, igniting innovation, and upholding the highest standards of quality in all our EV-focused initiatives. By spearheading EV research and technology, we aim to be the reliable partner for researchers striving to unravel the complexities of intercellular communication and harness the therapeutic benefits of EVs.

Distributed by
CliniSciences

2. DIA Proteomics/ Phosphoproteomics Service

Introduction

DIA (Data-independent acquisition) is a new mass spectrometry technique developed in recent years and belongs to the non-labeled proteomics methods. With data-independent scanning mode, it divides the entire scanning range of the mass spectrum into several windows, and then detects all the ions in each window, so it can obtain the information of all the ions in the samples without any omission or discrepancy. This technique can reduce the number of missing values in the sample assay while improving quantitative accuracy and reproducibility, enabling highly stable and accurate proteomic quantification in large sample queues.

We can provide 3 kinds of DIA services including Direct-DIA, GPF-DIA, and Fractionated DDA-base DIA. Samples provided for analysis include and are not limited to cells, tissues, paraffin sections, exosomes, etc. We also provide Proteomics Preprocessing Kit (#EV01-03-01) and Phosphopeptide Enrichment Kit (#EV01-05-01) for your proteomics experiments.

Procedure




FIGURE 1. THE FLOWCHART OF THE 3 DIA DETECTION TECHNIQUES

Advantages

- Full Scan Mode: a non-dependent data acquisition mode that greatly reduces high abundance interference and comprehensively acquires high and low abundance signals, improving data coverage by nearly 40 %.
- High reproducibility and strong stability: nearly 40% increase in repetition for identifying large number of samples, nearly 1X increase in quantitative precision.
- High quantitative accuracy: quantitative power is closed to the gold standard (SRM/MRM targeting technology).

Application

It can be applied to build up the information base of biological samples, discover biomarkers and etc. in various biological fields, such as biotechnology, agriculture, forestry, husbandry and medicine.

3. Label free Proteomics

Introduction

Label-Free shotgun is an isotope labeling-independent technique for non-labeled proteome quantification. It can analyze peptide fragments of proteins by LC-MS/MS, without the need to use expensive stable isotope tags as internal standards, but only need to analyze mass spectrometry data generated during large-scale identification of proteins, and compare the signal intensities of peptide fragments in different samples, so as to quantify the corresponding proteins of peptides. Label Free quantitative analysis technology can directly obtain the relative quantitative information of the protein represented by a peptide by comparing the peak area of the same peptide in the mass spectrometry.

Procedure




FIGURE 2. THE FLOWCHART OF LABEL FREE TECHNOLOGY (CREATED WITH MED PEER (WWW.MEDPEER.CN))

Advantages

- Cost-effective: no need for expensive isotope labeling reagents and short experimental cycles.
- Widely applicable: The "presence or absence of differential proteins" can be compared between species.
- Flexible and convenient: Unlimited number of samples, flexible and convenient.

Application

It can be applied for protein identification and quantification, screening and identification of differential proteins in various fields, such as agriculture, forestry, husbandry, basic medicine and clinical diagnosis and biomedicine.

4. Identification of interacting proteins

Introduction

LC-MS/MS protein identification technology can be used to identify proteins in purified solution samples, such as IP, Co-IP, pull-down, etc., and obtain a series of proteins interacting with the target protein.

Proteins are first digested into peptides by trypsin, and the peptides are ionized in the mass spectrometer to carry a certain amount of charge, which can be analyzed by the detector to obtain the mass-to-charge ratio (m/z) of each peptide (Primary mass spectrometry). To obtain more detailed sequence information of the peptides, some of the peptides were fragmented and analyzed again to produce secondary mass spectrometry. Finally, the mass spectrometry search software was used to select the appropriate protein database, and all the mass spectrometry data obtained were analyzed, while the identified proteins were sorted in the form of scores.

We also provide Co-IP kit (#EV 03-01-02) and one-stop service from sample preparation, IP enrichment, pretreatment and mass spectrometry.

Procedure



FIGURE 3. THE FLOWCHART FOR LC-MS/MS PROTEIN IDENTIFICATION

Advantages

- Widely applicable: It's suitable for protein identification of complex samples, can identify multiple proteins simultaneously.
- High sensitivity and high accuracy: It can capture comprehensive and detailed peptide information.

Application

It is applied for the identification of protein solutions prepared by IP, Co-IP, and pull-down methods in a variety of basic research in biotechnology, agriculture, forestry, husbandry, and medicine.

5. SDS-PAGE enzymatic assay

Introduction

It can identify proteins by LC-MS/MS for glue spot samples (2D glue dug spots), glue strip samples (SDS-PAGE samples), and WB membrane samples.

First, the protein bands were digested into peptide mixtures with enzyme, which were separated by HPLC to obtain a number of simpler fractions that were analyzed in a high-resolution mass spectrometer. Based on the secondary mass spectrometry information and the corresponding database, the software can obtain the exact sequence of peptides by combining match scoring and mismatch filtering, and then stitch the complete sequence to realize the identification of proteins.

Procedure



FIGURE 4. THE FLOWCHART OF SDS- PAGE ENZYMATICAL ASSAY

Advantages

- Easy and fast: 2D glue spot detection is fast and sample pre-treatment is easy.
- High requirement for samples: protein need to be pure or relatively homogenous.

Application

It is applied to the identification of SDS-PAGE adhesive strips prepared by IP, Co-IP, and pull-down methods in a variety of basic research in biotechnology, agriculture, forestry, husbandry, and medicine.

How to navigate through our website?

Homepage: When you first visit our website, you'll land on the homepage. This is where you can find an overview of our services, products, or any important announcements. From here, you can navigate to other sections of the website.



Search for a product by description, reference, CAS, Number, target...

Search Bar: Look for specific Products and equipment

Search for a product (description, reference, CAS number, target ...)

Search by selecting either RUO or IVD products

Refine your search by choosing from the following options

Apply Filter

Refine your search :

RUO CE / IVD

Conjugate (100) +

Reactivity species (100) +

Class (66) +

Host species (100) +

Cell name (2) +

Application (100) +

Clone (100) +

APPLY FILTERS

REINITIALIZE

SHARE THIS SEARCH

Search result : 39126007 product found

Cat#	Description	Cond.
NB-45-00042-100	Super NI-NTA Agarose Resin	100ml
NB-45-00042-25	Super NI-NTA Agarose Resin	25ml
NB-45-00058-4	Proteus 1 -step Batch Mini Spin Column Pack	40pc
NB-12-6001-3	NeoLine pipette 2-20 µl	1unit
NB-12-0023C	Mini Centrifuge N500C @10,000rpm (Including 6x1.5/2.0ml angle rotor)	1pcs
NB-03-0160	Proteinase K (Powder)	100mg
NB-60-0001	NeoPrep mini	50columns
NB-12-8001-19	Combs for NeoPRO4 mini (1.5mm, 15 wells)	5pieces
NB-12-8001-20	Spacer glasses flat for NeoPRO4 mini (0.75mm, 100*83mm)	5pieces
NB-12-8001-04	Short glasses flat for NeoPRO4 mini (1.0mm, 100*73mm)	10pieces

“Remember, if you ever get lost or have questions, don’t hesitate to reach out to our customer support team. They’re here to help you address any concerns you may have.”

clinisciences

Austria

Company: CliniSciences GmbH
Address: Sternwartestrasse 76, A-1180
Wien - Austria
Telephone: +43 720 115 580
Fax: +43 720 115 577
Email: oesterreich@clinisciences.com
Web: <https://www.clinisciences.com>



Belgium

Company: CliniSciences S.R.L
Address: Avenue Stalingrad 52, 1000
Brussels - Belgium
Telephone: +32 2 31 50 800
Fax: +32 2 31 50 801
Email: belgium@clinisciences.com
Web: <https://www.clinisciences.com>



Denmark

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: denmark@clinisciences.com
Web: <https://www.clinisciences.com>



Finland

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: suomi@clinisciences.com
Web: <https://www.clinisciences.com>



France

Company: CliniSciences S.A.S
Address: 74 Rue des Suisses, 92000
Nanterre- France
Telephone: +33 9 77 40 09 09
Fax: +33 9 77 40 10 11
Email: info@clinisciences.com
Web: <https://www.clinisciences.com>



Germany

Company: Clinisciences GmbH
Address: Wilhelm-Mausier-Str. 41-43,
50827 Köln - Germany
Telephone: +49 221 9498 320
Fax: +49 221 9498 325
Email: deutschland@clinisciences.com
Web: <https://www.clinisciences.com>



Iceland

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: island@clinisciences.com
Web: <https://www.clinisciences.com>



Ireland

Company: CliniSciences Limited
Address: Ground Floor, 71 lower Baggot street
Dublin D02 P593 - Ireland
Telephone: +353 1 6971 146
Fax: +353 1 6971 147
Email: ireland@clinisciences.com
Web: <https://www.clinisciences.com>



Italy

Company: CliniSciences S.r.l
Address: Via Maremmana inferiore 378
Roma 00012 Guidonia Montecelio - Italy
Telephone: +39 06 94 80 56 71
Fax: +39 06 94 80 00 21
Email: italia@clinisciences.com
Web: <https://www.clinisciences.com>



Netherlands

Company: CliniSciences B.V.
Address: Kraaijenhoffstraat 137A,
1018RG Amsterdam, - Netherlands
Telephone: +31 85 2082 351
Fax: +31 85 2082 353
Email: nederland@clinisciences.com
Web: <https://www.clinisciences.com>



Norway

Company: CliniSciences AS
Address: c/o MerVerdi Munkerudtunet 10
1164 Oslo - Norway
Telephone: +47 21 988 882
Email: norge@clinisciences.com
Web: <https://www.clinisciences.com>



Poland

Company: CliniSciences sp.Z.o.o.
Address: ul. Rotmistrza Witolda Pileckiego 67
lok. 200 - 02-781 Warszawa -Poland
Telephone: +48 22 307 0535
Fax: +48 22 307 0532
Email: polska@clinisciences.com
Web: <https://www.clinisciences.com>



Portugal

Company: Quimigen Unipessoal LDA
Address: Rua Almada Negreiros, Lote 5, Loja 14,
2615-275 Alverca Do Ribatejo - Portugal
Telephone: +351 30 8808 050
Fax: +351 30 8808 052
Email: info@quimigen.com
Web: <https://www.quimigen.pt>



Spain

Company: CliniSciences Lab Solutions
Address: C/ Hermanos del Moral 13
(Bajo E), 28019, Madrid - Spain
Telephone: +34 916 750 700
Fax: +34 91 269 40 74
Email: espana@clinisciences.com
Web: <https://www.clinisciences.com>



Sweden

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: sverige@clinisciences.com
Web: <https://www.clinisciences.com>



Switzerland

Company: CliniSciences AG
Address: Fracht Ost Flughafen Kloten
CH-8058 Zürich - Switzerland
Telephone: +41 (044) 805 76 81
Fax: +41 (044) 805 76 75
Email: switzerland@clinisciences.com
Web: <https://www.clinisciences.com>



UK

Company: CliniSciences Limited
Address: 11 Progress Business center, Whittle
Parkway, SL1 6DQ Slough- United Kingdom
Telephone: +44 (0)1753 866 511
or +44 (0) 330 684 0982
Fax: +44 (0)1753 208 899
Email: uk@clinisciences.com
Web: <https://www.clinisciences.com>



USA

Company: CliniSciences LLC
Address: c/o Carr Riggs Ingram,
500 Grand Boulevard, Suite 210 Miramar
Beach, FL 32550- USA
Telephone: +1 850 650 7790
Fax: +1 850 650 4383
Email: usa@clinisciences.com
Web: <https://www.clinisciences.com>



