

# Next Generation Sequencing

Next Generation Sequencing

Providing Innovative Reagents for Life Sciences since 2006

**Distributed by**  
**CliniSciences Group**

# CONTENTS

02

13

24

35

Distributed by  
**CliniSciences Group**

# DNA Library Preparation

*TransNGS*<sup>®</sup> DNA Library Prep Kit for Illumina<sup>®</sup>  
*TransNGS*<sup>®</sup> DNA Library Prep Kit for MGI<sup>®</sup>  
*TransNGS*<sup>®</sup> Fragmentase DNA Library Prep Kit for Illumina<sup>®</sup>  
*TransNGS*<sup>®</sup> Fragmentase DNA Library Prep Kit for MGI<sup>®</sup>  
*TransNGS*<sup>®</sup> Circularization Kit For MGI<sup>®</sup>  
*MagicPure*<sup>®</sup> Size Selection DNA Beads

# RNA Library Preparation

*TransNGS*<sup>®</sup> rRNA Depletion Kit (Human/Mouse/Rat)  
*TransNGS*<sup>®</sup> Magic rRNA Depletion Kit (Bacteria)  
*MagicPure*<sup>®</sup> mRNA Kit  
*TransNGS*<sup>®</sup> RNA-Seq Library Prep Kit for Illumina<sup>®</sup>  
*TransNGS*<sup>®</sup> Fast RNA-Seq Library Prep Kit for Illumina<sup>®</sup>  
*TransNGS*<sup>®</sup> Fast Stranded RNA-Seq Library Prep Kit for MGI<sup>®</sup>

# Pathogenic Microorganism Detection

*TransGuard*<sup>®</sup> Disposable Virus Sampling Tube  
*TransGuard*<sup>®</sup> Fecal DNA Sampling Tube  
*TransGuard*<sup>®</sup> Buccal Swab DNA Preservation Buffer  
*TransNGS*<sup>®</sup> Host DNA Depletion Kit  
*EasyPure*<sup>®</sup> Microbiome DNA Isolation Kit  
*MagicPure*<sup>®</sup> Microbiome DNA Isolation Kit  
*EasyPure*<sup>®</sup> Viral DNA/RNA Kit  
*MagicPure*<sup>®</sup> Fly 96 Viral DNA/RNA Kit  
*MagicPure*<sup>®</sup> Stool and Soil Genomic DNA Kit  
*TransNGS*<sup>®</sup> Fragmentase DNA Library Prep Kit for Illumina<sup>®</sup>

# Modifying Enzymes

T4 DNA Ligase for NGS  
DNA polymerase I Klenow Fragment  
T4 DNA Polymerase  
T4 Polynucleotide Kinase  
*TransNGS*<sup>®</sup> Library Amplification SuperMix

# TRANSGEN

# DNA Library Prep



Distributed by  
**CliniSciences Group**

# 01

# Common Library Preparation

## TransNGS<sup>®</sup> DNA Library Prep Kit for Illumina<sup>®</sup> (KP201)


### Features

- High library conversion rate.
- Suitable for wide range of sequencing methods: whole genome sequencing, target gene sequencing, meta-genomic sequencing, immunoprecipitation sequencing, exome sequencing, or other targeted capture sequencing.

### Applications


- Whole genome sequencing.
- Target gene sequencing.
- Exon sequencing / other targeted capture sequencing.
- Metagenomic sequencing.
- Co-immunoprecipitation sequencing.

### Schematic Diagram of Library Preparation



The IS position is indicated by a dotted line, which means that some libraries do not have this Index


### Operation flow chart




### Comparison Between Competing Products

TransGen (KI341) and Company V products were used for DNA library construction using 1ng, 10ng, 100ng, and 1000ng HeLa cell and smear DNA. The results showed that the peak and productivity of the library, the quality of sequencing result, and GC distribution generated by TransGen is as good as competing products.

Library yield



Library peak pattern





Sequencing data quality

Sample	Raw reads	Clean reads	Effective Rate(%)	Error Rate(%)	Q20(%)	Q30(%)	GC Content(%)
T-1	162024748	153782002	94.91	0.03	96.17	91.18	42.11
T-10	162024748	154921836	95.62	0.03	97.85	92.78	43.26
T-100	162024748	157370662	97.13	0.03	97.93	93.71	40.63
T-1000	162024748	158685938	97.94	0.03	97.67	92.14	40.2
V-1	162024748	158473152	97.81	0.03	96.45	90.79	42.82
V-10	162024748	158371574	97.75	0.03	97.33	92.86	43.42
V-100	162024748	158742128	97.97	0.03	97.6	93.22	41
V-1000	162024748	159029382	98.15	0.03	97.45	92.83	40.92


Sample name	Clean reads	Optical/PCR duplicate:	Unmapped reads:	Total mapped	Uniquely mapped	Average depth(X)	Coverage_at least_1X(%)	Coverage_at least_4X(%)
T-1	153782002	39022950 (25.38%)	1856324 (1.21%)	112902728 (73.42%)	103278748 (67.16%)	5.38	92.13	61.32
T-10	154921836	31328626 (20.22%)	981921 (0.63%)	122611289 (79.14%)	113816737 (73.47%)	6.04	92.44	65.52
T-100	157370662	29733483 (18.89%)	841402 (0.53%)	126795777 (80.57%)	117990902 (74.98%)	6.15	94.61	72.17
T-1000	158685938	30462223 (19.20%)	1070225 (0.67%)	127153490 (80.13%)	116148077 (73.19%)	5.94	94.54	70.71
V-1	158473152	38619248 (24.37%)	1009004 (0.64%)	118844900 (74.99%)	110511798 (69.74%)	5.71	93.59	65.88
V-10	158371574	32042078 (20.23%)	855095 (0.54%)	125474401 (79.23%)	116888800 (73.81%)	6.11	93.57	68.04
V-100	158742128	25962879 (16.36%)	727321 (0.46%)	132051928 (83.19%)	123571752 (77.84%)	6.43	95.07	75.46
V-1000	159029382	27823753 (17.50%)	711773 (0.45%)	130493856 (82.06%)	122082443 (76.77%)	6.33	95.12	75.03

GC distribution





correlation analysis



# TransNGS<sup>®</sup> DNA Library Prep Kit for MGI<sup>®</sup> (KP221)


**Features**

- High library conversion rate.
- Suitable for wide range of sequencing methods: whole genome sequencing, target gene sequencing, exosome sequencing, metagenomic sequencing, and immunoprecipitation sequencing.
- Great data quality


**Applications**

- Whole genome sequencing.
- Target gene sequencing.
- Metagenomic sequencing.
- Co-immunoprecipitation sequencing.
- Exon sequencing / other targeted capture sequencing.

**Schematic Diagram of Library Preparation**



**Operation flow chart**




## Comparison Between Competing Products

(1) comparison of library yield using different amount of samples

TransGen and Company V products were used for dsDNA library construction using different amount of fragmented Hela cell DNA. The results showed that TransGen product generates higher efficiency in library construction than competing products.

Library yield




(2) comparison of sequencing result

TransGen and Company V products were used for dsDNA library construction using different amount of fragmented Hela cell DNA. The results showed that the quality of sequencing result, GC distribution and correlation generated by TransGen is as good as competing products.

Comparison of sequencing results

Sample	Clean reads	Optical/PCR duplicate	Unmapped reads	Total mapped	Multiple mapped	Uniquely mapped	Q20(%)	Q30(%)	GC Content(%)
T-1	42720938	684528 (1.60%)	230611 (0.54%)	41805799 (97.86%)	2456018 (5.75%)	39349781 (92.11%)	96.12	90.99	43.11
T-10	42715766	123735 (0.29%)	215356 (0.50%)	42376675 (99.21%)	2474022 (5.79%)	39902653 (93.41%)	96.19	91.14	43.19
T-100	42720042	40336 (0.09%)	395135 (0.92%)	42284571 (98.98%)	2191110 (5.13%)	40093461 (93.85%)	96.43	91.66	40.78
T-1000	42720058	48045 (0.11%)	363641 (0.85%)	42308372 (99.04%)	2225122 (5.21%)	40083250 (93.83%)	96.54	91.89	41.02
V-1	42711428	914214 (2.14%)	205217 (0.48%)	41591997 (97.38%)	2560892 (6.00%)	39031105 (91.38%)	96.27	91.34	44.00
V-10	42710566	124841 (0.29%)	189493 (0.44%)	42396232 (99.26%)	2599217 (6.09%)	39797015 (93.18%)	95.96	90.63	43.90
V-100	42710232	41899 (0.10%)	418119 (0.98%)	42250214 (98.92%)	2176671 (5.10%)	40073543 (93.83%)	96.15	91.04	40.49
V-1000	42704112	40522 (0.09%)	420547 (0.98%)	42243043 (98.92%)	2159708 (5.06%)	40083335 (93.86%)	96.06	90.82	40.48

GC distribution



# 02

## Fragmentase Library Preparation

### TransNGS<sup>®</sup> Fragmentase DNA Library Prep Kit for Illumina<sup>®</sup> (KP231)


#### Features

- High library conversion rate.
- Suitable for wide range of sequencing methods: whole genome sequencing, target gene sequencing, meta-genomic sequencing, and exome sequencing.

#### Applications


- Whole genome sequencing.
- Target gene sequencing.
- Metagenomic sequencing.
- Co-immunoprecipitation sequencing.
- Exon sequencing / other targeted capture sequencing.

#### Schematic Diagram of Library Preparation



The i5 position is indicated by a dotted line, which means that some libraries do not have this Index

#### Operation flow chart




[ ] not included in the kit

★ Locations where clip sorting can be done


### Experimental data display

(1) Comparison of library yields with different starting sample amounts  
 DNA library was constructed using TransGen UDI long adapter (KI341) with 1 ng, 10 ng, 100 ng and 500 ng restricted HeLa cell genomic DNA. The results showed that the average library yields is 1 µg.



(2) library yield of different species

TransGen products were used to construct libraries for different sample species, including corn, Staphylococcus aureus, and soil microorganisms. The results showed that the average library yields is 1 µg.




(3) Sequencing result

The results show that TransGen products construct DNA libraries efficiently with 1-500 ng samples.

### Sequencing data quality

Sample	Q20 (%)	Q30 (%)	GC Content (%)	Total mapped (%)	Coverage at-1X (%)	Coverage at-4X (%)	Duplication Rate (%)	Depth (x)
T-1	95.85	90.35	40.95	72.88	93.48	69.88	25.68	6.26
T-10	96.49	91.26	40.71	83.26	94.75	78.81	15.53	7.29
T-100	96.15	90.65	40.77	78.27	96.41	85.81	20.61	8.64
T-500	96.97	92.33	40.74	80.66	95.55	85.76	18.28	8.63

### GC distribution



# TransNGS<sup>®</sup> Fragmentase DNA Library Prep Kit for MGI<sup>®</sup> (KP241)


## Features

- High library conversion rate.
- Suitable for wide range of sequencing methods: whole genome sequencing, target gene sequencing, meta-genomic sequencing, and exome sequencing.


## Applications

- Whole genome sequencing.
- Target gene sequencing.
- Metagenomic sequencing.
- Co-immunoprecipitation sequencing.
- Exon sequencing / other targeted capture sequencing.

## Schematic Diagram of Library Preparation



## Operation flow chart



## Comparison Between Competing Products

200 ng standard NA12878 human genome was treated by TransGen and Company A products to construct a DNA library. The results showed that the TransGen product was better than competing products in terms of sequencing data quality, alignment rate, SNP and InDel detection.

### Sequencing data quality

Sample	Q20 (%)	Q30 (%)	GC Content (%)	Total mapped (%)	Coverage at-4X (%)	Coverage at-10X (%)	Duplication Rate (%)	Depth (x)
T-200-1	97.24	93.05	40.43	99.99	98.62	97.34	1.24	30.31
T-200-2	97.82	94.04	40.67	99.99	98.78	97.58	1.24	30.34
A-200-1	97.41	94.16	39.45	99.99	98.54	96.00	1.22	30.63
A-200-2	97.38	93.74	39.45	99.99	98.88	96.75	0.96	30.85

### SNP and InDel detection

Sample	Total SNPs	Ti/Tv	Precision (%)	Consistency (%)	Total InDels	Precision (%)	Consistency (%)
T-200-1	3717711	2.04	99.73	95.49	827247	96.00	74.26
T-200-2	3723268	2.04	99.74		837459	96.03	
A-200-1	3700900	2.04	99.78	95.77	813681	96.79	74.76
A-200-2	3705928	2.04	99.83		823616	96.03	

# 03

## Related Products

### TransNGS<sup>®</sup> Circularization Kit For MGI<sup>®</sup> (KC101)


#### Features

- Short operation time
- high circularization efficiency

#### Applications


- Applicable for dsDNA library preparation using MGI high-throughput sequencing platform

#### Schematic Diagram of Library Preparation




#### Comparison Between Competing Products

100pmol double strand DNA library was generated by TransGen library preparation kit under MGI platform. Single-stranded circularized libraries were prepared using TransGen and Company V products respectively, and the circularization efficiency and amplification fold were calculated using products of the same quantity. The results showed that the TransGen product was better than competing products in terms of efficiency.



### Stability of the Product

Circularization efficiency was calculated using 1 pmol MGI double-stranded DNA library and fragments with MGI adapters, circularized by TransGen products in different production batches. The consistency of circularization efficiency shows that TransGen products have good batch-to-batch stability.



# MagicPure<sup>®</sup> Size Selection DNA Beads (EC401)

### Features

- Easy operation, designed for automated workstations.


### Applications

- DNA purification
- DNA condensation
- DNA fragment sorting in high-throughput sequencing library construction

### DNA purification and size selection

DNA samples were purified and sorted by TransGen products. Purified DNA product was identified by Agilent's high-sensitivity DNA chip, and TransGen's product shows high efficiency and accuracy.

### DNA Purification




Agilent high-sensitivity DNA electropherogram

Smear - Control dissolved in Nuclease-free Water; 0.6x~1.8x - DNA sample purified by the corresponding magnetic bead ratio

## DNA size selection

Average length of sorted fragments (bp)						
190~220	220~250	250~300	300~400	400~500	500~600	600~750
1 <sup>st</sup> volume ratio (DNA Beads:DNA)						
1.0×	0.9×	0.80×	0.70×	0.60×	0.55×	0.50×
2 <sup>nd</sup> volume ratio (DNA Beads:DNA)						
0.20×	0.20×	0.20×	0.20×	0.15×	0.15×	0.15×





Agilent high-sensitivity DNA electropherogram

Smear - Control dissolved in Nuclease-free Water; 1.0+0.2-0.5+0.15 - DNA sample purified by the corresponding magnetic bead ratio

## Comparison Between Competing Products



### (1) Size selection

DNA samples with same length were sorted using products from TransGen and Company B respectively. Selected fragments were analyzed by Agilent high-sensitivity DNA chips and the concentration was measured using Qubit. The results showed that under same conditions, the length and recovery amount of the product gives no significant difference.



### (2) Stability of the Product

After stored in 37°C for 3 weeks/6 weeks, products from TransGen and Company B (2-8°C as a control), DNA samples were fragmented (0.6x+0.1x) and analyzed by Agilent high-sensitivity DNA chips. Qubit was used to calculate the product concentration. The results showed that after high-temperature destructive experiments, the fragment length distribution and recovery amount are not significantly different from competing products, which shows good stability of the product.



# TRANSGEN

# RNA Library Prep

# B

Distributed by  
**CliniSciences Group**

# 01 rRNA Removal mRNA Capture

## TransNGS<sup>®</sup> rRNA Depletion Kit (Human/Mouse/Rat) (KD101)

### Features


- Removes up to 99% of rRNA in humans/mouse/rats effectively .
- Control qPCR Primer Sets included. Able to detect changes in rRNA and non-ribosomal RNA content before and after removal

### Applications

- Human/mouse/rat total RNA (100 ng-1 µg) samples
- Completely or partially degraded (e.g. FFPE RNA) RNA samples

### Stability of the product

Different batches of products were used to remove rRNA from 1 µg of total RNA (HepG2 cells), and rRNA Primer and non-rRNA Primer were used to perform qRT-PCR before and after rRNA removal. Since  $\Delta Cq$  values of different primers among different batches maintain relatively constant, the stability of the product is good.



### Comparison Between Competing Products

Products from TransGen and Company K were used to remove rRNA in 1 µg of human (H), mouse (M) and rat (R) total RNA samples respectively, and then library construction (TransGen, KP601) was performed for sequencing and data analysis. The result shows that the performance of TransGen products is consistent with, or even better than, competing products in terms of sequencing quality, rRNA removal rate, circRNA and lncRNA detection rate.

### Sequencing Quality

Sample	rRNA Rate(%)	Clean reads	Q20 (%)	Q30 (%)	GC content (%)	Total Mapped	Uniquely Mapped
H-CK	76.4	30611582	96.36	90.37	48.1	24914766 (81.39%)	21400556 (69.91%)
H-T	3.13	33860104	97.28	91.69	44.52	30467321 (89.98%)	29898471 (88.3%)
H-K	4.1	25323776	96.97	87.72	45.3	22631858 (89.37%)	22135512 (87.41%)
M-CK	82.2	31734716	96.37	90.39	50.52	26266824 (82.77%)	23902588 (75.32%)
M-T	3.37	26782378	97.13	89.71	44.91	23541710 (87.9%)	21366981 (79.78%)
M-K	5.7	28705540	97.05	88.71	45.11	25039842 (87.23%)	22829515 (79.53%)
R-CK	76.1	29980720	96.92	91.45	47.73	23828676 (79.48%)	20344916 (67.86%)
R-T	0.6	36024365	97.09	89.21	45.01	31315980 (86.93%)	28981601 (80.45%)
R-K	1.47	38587178	97.07	88.96	45.06	33566986 (86.99%)	30862024 (79.98%)


H-CK/M-CK/R-CK : Samples without rRNA removal

rRNA removal rate

Species	Sample	rRNA Removal (%)	rRNA Rate (%)
Human	H-T	99.88	3.13
	H-K	99.54	4.1
Mouse	M-T	99.79	3.37
	M-K	98.82	5.7
Rat	R-T	99.58	0.6
	R-K	99.03	1.47

rRNA Removal is analyzed by Mirabait software, rRNA Rate is compared using nt library, where Mirabait software is more accurate.

circRNA/lncRNA Analysis



# TransNGS<sup>®</sup> Magic rRNA Depletion Kit (Bacteria) (KD401)


## Features

- Removes more than 95% of bacterial rRNA effectively.
- Suitable for single or mixed bacteria.
- Suitable for feces, soil, oral cavity and other samples.
- Little damage to mRNA.


## Applications

- Suitable for both Gram-positive and Gram-negative bacteria (removal efficiency may vary for different bacterial species).
- Suitable for 500 ng–2 µg total RNA samples.
- Suitable for RNA samples with good integrity or moderate degradation (RIN value >4)

## Schematic diagram




## Operation flow chart




## Gram-positive bacteria and Gram-negative bacteria rRNA removal efficiency

qPCR and sequencing of *S. aureus* and *E. coli* RNA samples treated by TransGen product were used to detect rRNA removal efficiency. The results showed that rRNA removal rate of TransGen product is more than 98%.

Changes in Cq values before and after rRNA removal based on qPCR




Sequencing of residual rRNA reads




### Strain compatibility

Products from TransGen were used to remove rRNA and construct RNA library and RNA sequencing of 500 ng~1 µg bacterial total RNA sample. Sequencing data was compared with the Silva database to calculate the residual ratio of rRNA. The results showed that in most of tested species, rRNA residues were within 5%, which indicate that TransGen products have good compatibility and removal efficiency.




### rRNA removal effect on complex samples

TransGen products were used to remove rRNA from soil, saliva, and fecal samples respectively, and then RNA libraries were constructed and sequenced. The result were compared with Silva database to calculate the ratio of rRNA reads. The results showed that after removal, the proportion of rRNA in complex samples was significantly reduced.



### rRNA removal performance on degraded RNA sample

Products from TransGen were used to remove rRNA and construct RNA library on 1 µg of *Staphylococcus aureus* RNA samples under different degradation (RNA samples after high temperature treatment were tested by Agilent 2100, and the RIN values were in the intervals of 4~5, 5~6, 6~7, 7~8 and >8 respectively). The result was compared with the Silva database to calculate the residual rate of rRNA reads, which shows that TransGen products are suitable for intact or moderately degraded RNA with RIN values >4.



# MagicPure<sup>®</sup> mRNA Kit (EC511)

## Features

- Easy Operation.
- High yield and purity.
- Suitable for magnetic high-throughput nucleic acid extraction equipment.

## Sample Requirement

- 0.1-10 µg purified total RNA with good integrity (RIN value ≥8).

## Applications


- Products are suitable for RT-PCR, qRT-PCR, second-generation sequencing, etc.

## Comparison Between Competing Products


Capture and Libraries Construction of Samples from Different Species

TransGen and Company V products were used to capture the total RNA of human, rat, mouse, tobacco, lily, tomato, soybean, arabidopsis, wheat and rice respectively, where 200 ng of total RNA input to human, rat and mouse; 1000 ng of total RNA input to other species. Library construction kit (TransGen, KP701) was used to construct RNA library, where library yield and peak pattern were analyzed with sequencing. Results show that TransGen products perform consistently with Company V products.

Library Yield



Peak Pattern of Library Construction



# 02

# Library Preparation

## TransNGS<sup>®</sup> RNA-Seq Library Prep Kit for Illumina<sup>®</sup> (KP601)


### Features

- High library conversion rate.
- High data quality.


### Applications

- Whole transcriptome sequencing.
- Gene expression analysis.
- Single nucleotide variation analysis.
- Variable shear detection.
- Fusion gene detection.
- Non-coding RNA and RNA precursor analysis

### Schematic Diagram of Library Construction Principle



### Operation Flow Chart




### Comparison with Competing Products

TransGen, Company V, and Company products were used for library construction from 1 µg, 100 ng of total human blood RNA samples after rRNA removal (TransGen, KD101). The results show that TransGen products perform consistently with competitors in terms of sequencing data quality, gene expression levels, and correlation.

Sequencing quality



Sample	rRNA Rate (%)	Clean reads	Q20 (%)	Q30 (%)	GC content (%)	Total Mapped	Uniquely Mapped
H (1 µg)-T	0.78	27709932	98.85	95.67	49.03	25146763 (90.75%)	23079602 (83.29%)
H (100 ng)-T	2.31	24376592	98.96	95.99	48.57	21702480 (89.03%)	20093625 (82.43%)
H (1 µg)-V	0.88	24545480	98.92	95.85	46.92	22486114 (91.61%)	20961840 (85.40%)
H (100 ng)-V	2.91	26241828	98.94	95.57	49.11	23575658 (89.84%)	21841073 (83.23%)
H (1 µg)-N	0.69	19833710	98.72	95.22	50.75	17,980,295 (90.66%)	16,545,319 (83.42%)
H (100 ng)-N	2.55	20951360	98.65	95.00	47.95	19,036,259 (90.86%)	17,730,696 (84.63%)

Compare regional distributions




Gene expression and differential analysis

FPKM interval	H (1 µg)-T	H (100 ng)-T	H (1 µg)-V	H (100 ng)-V	H (1 µg)-N	H (100 ng)-N
0~1	52,346 (71.27%)	52,782 (71.86%)	51,114 (69.59%)	50,586 (68.87%)	51,300 (70.95%)	50,587 (69.96%)
1~3	6,098 (8.30%)	5,929 (8.07%)	8,588 (11.69%)	8,842 (12.04%)	8,161 (11.29%)	8,606 (11.90%)
3~5	3,784 (5.15%)	3,470 (4.72%)	3,397 (4.62%)	3,347 (4.56%)	1,893 (2.62%)	1,975 (2.73%)
5~15	4,043 (5.50%)	4,133 (5.63%)	3,765 (5.13%)	3,640 (4.96%)	3,679 (5.09%)	3,595 (4.97%)
15~60	4,870 (6.63%)	4,830 (6.58%)	4,582 (6.24%)	4,715 (6.42%)	5,031 (6.96%)	5,185 (7.17%)
>60	2,309 (3.14%)	2,306 (3.14%)	2,004 (2.73%)	2,320 (3.16%)	2,245 (3.10%)	2,361 (3.27%)



Correlation analysis



# TransNGS<sup>®</sup> Fast RNA-Seq Library Prep Kit for Illumina<sup>®</sup> (KP701)


## Features

- Fast library preparation
- High library conversion rate
- High data quality

## Applications

- Whole transcriptome sequencing
- Gene expression analysis
- Single nucleotide variation analysis
- Variable splicing testing
- Fusion genetic testing
- Analysis of non-coding RNA and RNA precursor

## Schematic Diagram of Library Preparation



Principle chart of strand-specific library construction  
(No dUTP incorporation and UDG digestion steps in second-strand synthesis of non-strand-specific library)

## Operation Flow Chart




☐ The dotted lines are not included in this kit flow  
★ represents where the fragment can be sorted

Flow chart of strand-specific library construction  
(No dUTP incorporation and UDG digestion steps in second-strand synthesis of non-strand-specific library)


## Comparison with Competing Products

mRNA purification kit was used to capture mRNA from total RNA of human, mouse, wheat and tobacco, in which 200 ng total RNA is input for human and mouse, and 1000 ng total RNA is input for wheat and tobacco. TransGen and Company V products were used for RNA library construction. After library sorting, the library yield, library peak shape, and sequencing analysis is performed. Results show that TransGen products perform consistently with Company V products.

### Library yield



### Library peak type



### Sequencing quality

Sample name	clean ratio (%)	low quality ratio (%)	adapter ratio (%)	Q20 (%)	Q30 (%)	GC Content (%)	duplicate (%)	Total mapped (%)	Strand specific (%)
Human-T	92.59	6.88	0.53	98.24	94.09	50.05	31.63	96.42	98.62
Human-V	89.25	8.21	2.54	98.21	93.95	49.17	33.91	96.66	97.86
Mouse-T	92.44	7.09	0.47	97.77	92.89	50.53	30.12	95.65	98.30
Mouse-V	88.36	9.09	2.55	98.19	93.95	49.17	31.05	95.35	97.54
Wheat-T	92.94	6.47	0.59	98.25	94.15	54.90	26.51	93.96	98.72
Wheat-V	88.74	8.71	2.55	98.16	93.94	54.05	26.91	94.20	97.85
Tobacco-T	92.48	6.88	0.64	97.56	92.41	45.14	20.46	56.17	98.43
Tobacco-V	88.98	7.98	3.04	98.32	94.17	44.56	21.75	53.52	97.64

# TransNGS<sup>®</sup> Fast Stranded RNA-Seq Library Prep Kit for MGI<sup>®</sup> (KP801)

## Features

- Fast library preparation
- High library conversion rate
- High data quality


## Applications

- Whole transcriptome sequencing
- Gene expression analysis
- Single nucleotide variation analysis
- Variable splicing testing
- Fusion genetic testing
- Analysis of non-coding RNA and RNA precursor


## Comparison with Competing Products

mRNA purification kit was used to capture mRNA from total RNA of human, mouse, and arabidopsis, in which 200 ng total RNA is input for human and mouse, and 1000 ng total RNA is input for arabidopsis. TransGen and Company V products were used for RNA library construction. After library sorting, the library yield, library peak shape, and sequencing analysis is performed. Results show that TransGen products perform consistently with Company V products.

### Library yield



### Library peak type



### Sequencing quality

Sample	Q20 (%)	Q30 (%)	GC Content (%)	Total mapped (%)	Duplication Rate (%)
Human-T	97.82	94.19	50.53	94.22	19.45
Human-V	97.79	94.19	50.70	94.19	19.41
Mouse-T	97.81	94.18	49.44	93.68	16.46
Mouse-V	97.92	94.48	49.25	92.66	17.64
Arabidopsis-T	97.90	94.39	45.18	92.44	23.62
Arabidopsis-V	98.00	94.69	45.32	92.52	23.39

# TRANSGEN Pathogenic Microorganism Detection



Distributed by  
**CliniSciences Group**

# 01

## Sample Preservation


### TransGuard<sup>®</sup> Disposable Virus Sampling Tube (ES101)

#### Features


- Inactivate virus: To avoid secondary contamination during the process of transportation and detection.
- Stable preservation: Protect viral nucleic acids from degradation at room temperature for 1 week.
- Transport and preserve at room temperature without the need of cold chain.
- High detection rate: Compatible with a variety of virus nucleic acid extraction reagents.
- Simple to operate: No professional training is required, allowing self-sampling.

#### Stable preservation

TransGen products were used to collect and preserve the serially diluted ( $10^5$ - $10^6$  copies/ml) Newcastle disease virus. After 7 days, the nucleic acid was extracted and qRT-PCR was performed to test the preservation effect. The results show that TransGen products can effectively preserve viruses at different concentrations.



TransGen and Company D products were used to collect and preserve Newcastle disease virus. After 7 days, the nucleic acid was extracted and qRT-PCR was performed to test the preservation effect. The results show that TransGen products exhibits a better preservation effect.




### TransGuard<sup>®</sup> Fecal DNA Sampling Tube

#### Features

- Inactivate safely: To avoid secondary contamination during the process of transportation and detection.
- Stable preservation: Protect sample nucleic acids from degradation at room temperature for 3 months.
- Convenient transportation and preservation: Transport and preserve at room temperature without the need of cold chain.
- Simple to operate: No professional training is required, allowing self-sampling.

#### Stable preservation

Stool samples were preserved using products from TransGen, Company T and PBS. The preservation effect was evaluated by agarose gel electrophoresis for gDNA extracted after 0 day, 60 days and 120 days stored at 25°C. The results showed that the integrity of the sample gDNA preserved by TransGen products were better than those of Company T and PBS, with little degradation.





# TransGuard<sup>®</sup> Buccal Swab DNA Preservation Buffer (ES103)

## Features

- Inactivate safely: To avoid secondary contamination during the process of transportation and detection.
- Stable preservation: Protect sample nucleic acids from degradation at room temperature for 6 months.
- Convenient transportation and preservation: Transport and preserve at room temperature without the need of cold chain.
- Simple to operate: No professional training is required, allowing self-sampling.

## Stable preservation

Four oral swabs were collected using TransGen products. The preservation effect was evaluated by agarose gel electrophoresis for gDNA extracted after 30 days and 120 days stored at 25°C. The results show that the sample gDNA preserved by TransGen products has good integrity.



# 02

## Host Nucleic Acid Removal

### TransNGS<sup>®</sup> Host DNA Depletion Kit (EH301)


#### Features

- Simple and time-saving: simple operation, and the nucleic acid removal of the human host is completed in 50 minutes.
- High-efficiency removal: special lysate can efficiently remove human host nucleic acid and improve the detection rate of microorganisms.
- Wide range of samples: suitable for buccal/throat swabs, sputum, pleural fluid, ascites, cerebrospinal fluid, amniotic fluid and other biological fluid samples.
- Compatible with a variety of extraction platforms: compatible with column extraction and magnetic bead method and other downstream extraction methods.


#### Effective Depletion of Host Nucleic Acid

TransGen, Company V and Company Q were used to process different samples, and the corresponding the column extraction products of pathogenic microorganism nucleic acid were used to extract nucleic acid, and the host conserved genes were quantitatively detected by qPCR. The results showed that TransGen could effectively remove host nucleic acid, and the removal efficiency was significantly better than that of Company V and Company Q products.


#### Throat Swab




#### Sputum



#### Ascites



#### Pleural Fluid



# 03

## Nucleic Acid Extraction

### EasyPure<sup>®</sup> Microbiome DNA Isolation Kit (EE401)


#### Features

- The extracted nucleic acid is of high quality and meets a variety of downstream testing needs.
- Strictly control the contaminations of background bacteria in reagents to reduce the risk of false positives.
- Suitable for microbial nucleic acid extraction from the throat swab, sputum, ascites and pleural fluid samples.


#### Efficient Extraction of Multiple Samples

High and low input of different bacteria diluents were added to throat swabs, sputum, pleural effusion and ascites samples, and nucleic acid was extracted using TransGen, Company V and Company Q products, respectively. qPCR was used to detect the extraction effect. The results showed that the extraction efficiency of TransGen products for different microorganisms in different samples was better than competing products.


#### Throat Swab




#### Sputum



#### Pleural Fluid



#### Ascites



# MagicPure<sup>®</sup> Microbiome DNA Isolation Kit (EC107)


## Features

- The extracted nucleic acid is of high quality and meets a variety of downstream detection requirements.
- Strictly control the contaminations of background bacteria in reagents to reduce the risk of false positives.
- Suitable for microbial nucleic acid extraction from the throat swab, sputum, ascites and pleural fluid samples.
- Suitable for high-throughput magnetic-rod nucleic acid extractor.


## Efficient Extraction of Multiple Samples

The high and low input of different bacteria diluents were added to the throat swab, sputum, ascites and pleural fluid samples. The nucleic acid was extracted using the TransGen column-based method and magnetic bead method, and the extraction effect was detected by qPCR. The results showed that both TransGen magnetic bead method and column based method can efficiently extract samples with high and low input.


### Throat Swab




### Sputum



### Pleural Fluid



### Ascites



# EasyPure<sup>®</sup> Viral DNA/RNA Kit (ER201)


## Features

- Use a unique lysis buffer to lyse virus and silica-based spin column to specifically adsorb DNA.
- High yield and high purity.
- Suitable for extracting viral DNA/RNA from a variety of samples such as whole blood, swab and cell supernatant.


## High Extraction Efficiency

Extract COVID-19 pseudovirus RNA using products of TransGen and Company TH respectively, followed by qRT-PCR to detect the ORF1ab and N genes. The results show that TransGen's kit has high extraction efficiency.

ORF1ab



N



# MagicPure<sup>®</sup> Fly 96 Viral DNA/RNA Kit


## Features

- Fast and simple, requiring only one wash, which enhances manual operation efficiency.
- High extraction sensitivity, improving virus detection rates.
- Extract fast, 14 minutes can be completed
- Suitable for extracting viral DNA/RNA from samples such as serum, swabs, and cell supernatants.
- Utilizes silicon-based magnetic beads for specific adsorption and purification of virus DNA/RNA, compatible with various 96-channel automated nucleic acid extractors.


## High Extraction Efficiency

Using COVID-19 pseudovirus samples diluted to 10<sup>4</sup>-fold and 10<sup>6</sup>-fold, nucleic acid extraction was performed using TransGen EC331-96 and EC311-96 products, followed by qPCR (TransGen, AQ322) to assess the extraction efficiency. The results show that EC331-96 product yields a higher extraction quantity.


ORF1ab



N



RNase P



# MagicPure<sup>®</sup> Stool and Soil Genomic DNA Kit (EC801)

## Features


- DNA-specific adsorption magnetic beads, high purity, and fast extraction speed.
- Suitable for extracting various types of soil samples, efficiently removing inhibitors from the samples.

## Comparison with competitors

### Soil sample


#### (1) High yield

Extract genomic DNA from different soil samples using products of TransGen and Company O respectively. The Qubit test results showed that TransGen products yielded higher DNA concentrations in the extracted DNA from all six soil samples.




#### (2) The extracted DNA has good integrity

Extract genomic DNA from different soil samples using products of TransGen and Company O respectively. Gel electrophoresis results show that TransGen products yielded better DNA quantity and integrity.




(3) Better suited for downstream applications

Fungal ITS sequences were amplified using 50 ng of genomic DNA extracted from different soil samples with TransGen and Company O products as templates. Gel electrophoresis results show that genomic DNA extracted with TransGen products is better suited for downstream PCR detection.




16S rDNA were amplified using 4 µL genomic DNA extracted from different soil samples with TransGen and Company O products as templates. Gel electrophoresis results show that genomic DNA extracted with TransGen products is better suited for downstream PCR detection.



Stool sample

(1) The extracted DNA has good integrity

Extract genomic DNA from different soil samples (Mouse feces 200 mg, human feces 200 µL) using products of TransGen and Company T respectively. Gel electrophoresis results show that TransGen products yielded better DNA quantity and integrity.



(2) Better suited for downstream applications

Extract genomic DNA using products of TransGen and Company T respectively from human and mouse feces, and then different amounts of genomic DNA were used as templates to amplify 16S rDNA. Gel electrophoresis results show that genomic DNA extracted with TransGen products is better suited for downstream PCR detection.



# TransNGS<sup>®</sup> Fragmentase DNA Library Prep Kit for Illumina<sup>®</sup> (KP231)

## Applications

- Wide applicability to various sample types.
- High library conversion rate.


## Applications:

- Whole genome sequencing.
- Targeted gene sequencing.
- Exome sequencing/other targeted capture sequencing.
- Metagenomic sequencing.


## Comparison with competing products:

Compared to our competitor, we conducted library construction using TransGen and Company V products on soil, feces, and sputum samples. The results showed that TransGen's product performance was consistent with that of our competitor in terms of library yield and peak shape, sequencing data quality, sequencing results, and sequence assembly.

Library yield




The peak shape of the library construction



Sequencing data quality

Sample	Clean reads	rmHost reads	rmHost Bases(bp)	Percent in clean(%)	Q20(%)	Q30(%)	GC Content(%)
Feces-T	7111142	7107050	1047093115	99.94	96.76	91.76	49.2
Feces-V	7368622	7365764	1096527329	99.96	87.69	78.57	48.96
Soil-T	7157144	7156848	1054467309	100	96.55	91.44	62.46
Soil-V	7359946	7355026	1091022557	99.94	96.89	91.89	62.85
Sputum-T	7126684	5711464	833888064	80.91	96.19	90.79	43.95
Sputum-V	7369508	6109882	904691171	83.34	97.15	92.21	45.11



Sequence assembly results

Sample	Contigs	Contig bases(bp)	N50(bp)	N90(bp)	Max(bp)	Min(bp)
Feces-T	158092	116284135	842	365	42833	200
Feces-V	131004	91218323	764	360	34560	200
Soil-T	116560	57703998	465	321	18461	200
Soil-V	122432	60987964	464	324	29016	200
Sputum-T	63449	47478014	822	356	232780	200
Sputum-V	66250	53287413	916	367	271388	200

# TRANSGEN

## Modifying Enzymes



Distributed by  
**CliniSciences Group**

# T4 DNA Ligase for NGS (LL101)


## Applications

- This is primarily used for adapter ligation during NGS library construction.
- It can also be used for cloning of restriction enzyme fragments.


## Comparison with competing products

Using both TransGen and Company V products for library adapter ligation in next generation sequencing, the results indicate that the library yield and peak shape of TransGen are consistent with those of Company V products.

### Library yield



### The peak shape of the library construction



# DNA polymerase I Klenow Fragment (LE201)


## Applications

The complementary filling of the 5' overhang and excision of the 3' overhang in double-stranded DNA.

## Comparison with competing products

Using a specific quantity of single-stranded cDNA (GFP gene, approximately 700 bp) as a template, dual-strand synthesis was performed using TransGen and Company N products separately. The reaction was carried out at 25°C for 10 minutes, and the yield was measured. The results indicate that the performance of TransGen products is consistent with that of Company N products in next-generation sequencing.

### Library yield



# T4 DNA Polymerase (LP201)


## Applications

- Smoothing of DNA 5' or 3' protruding ends.
- Synthesis of labeled DNA probes through displacement reaction.
- Subcloning after single-strand deletion.
- Synthesis of the second strand during gene site-directed mutagenesis process.

## Comparison with competing products

Using a certain amount of single-stranded cDNA (GFP gene, approximately 700 bp) as a template, we performed double-stranded synthesis using TransGen products and similar competitors. The reaction was carried out at 25°C for 10 minutes, and the yield was measured. The results indicate that the performance of TransGen products is consistent with that of the competitors.

Library yield




# T4 Polynucleotide Kinase (LK101)

## Applications

- Oligonucleotide, DNA, or RNA 5' end labeling for use as probes in Southern, Northern, EMSA, DNA sequencing primers, PCR primers, etc.
- Phosphorylation of oligonucleotide, DNA, or RNA 5' ends.
- Removal of 3' end phosphate groups.

## Comparison with competing products

Using 1U and 2U TransGen products, as well as similar competitor products, a 50 ng 250 bp PCR fragment was phosphorylated. The resulting product was digested using  $\lambda$  enzyme (with phosphorylated double-stranded DNA as the optimal substrate) and the remaining product was quantified. The results indicate that the phosphorylation efficiency of TransGen products is consistent with that of the competitors.



# TransNGS<sup>®</sup> Library Amplification SuperMix (KA101)

## Features


- Ultra high fidelity.
- Low preference.
- High sensitivity and specificity.
- Hot Start.

## Applications

Next-generation sequencing library amplification.


## Stability of the Product

Different batches of products were used to amplify 25 ng of DNA libraries with different GC contents (33%~72%) for 6 cycles. After the products were purified with 1.0× magnetic beads, Qubit was used to detect the concentration. The results show that the amplification efficiency of different batches of products for templates with different GC contents is basically the same.




## Comparison between competing products

1. Use TransGen, Company N, Company K and Company V products to amplify 25 ng of DNA libraries with different GC contents (33%~72%) respectively, and amplify for 6 cycles. After the products are purified with 1.0× magnetic beads, use Qubit detects concentration. Results show that TransGen products demonstrate higher amplification efficiency and lower GC bias.



2. Use TransGen and Company N products to amplify libraries with different GC contents. The amplified libraries are DNA fragments connected to the same adapter, and sequenced using the Illumina sequencing platform. The distribution of sequencing reads with different GC contents is shown in the curve in the figure. The library amplified using TransGen products is closer to the expected distribution.



Category	Product Name	Catalog Number	Specifications
DNA Library Preparation	TransNGS® DNA Library Prep Kit for Illumina®	KP201-11/03	12 rxns / 96 rxns
	TransNGS® DNA Library Prep Kit for MGI®	KP221-01/02	12 rxns / 96 rxns
	TransNGS® Fragmentase DNA Library Prep Kit for Illumina®	KP231-01/02	12 rxns / 96 rxns
	TransNGS® Fragmentase DNA Library Prep Kit for MGI®	KP241-01/02	12 rxns / 96 rxns
	TransNGS® Circularization Kit For MGI®	KC101-01/02	12 rxns / 96 rxns
	MagicPure® Size Selection DNA Beads	EC401-01/02/03/04	1 mL / 5 mL / 60 mL / 450 mL
RNA Library Preparation	TransNGS® rRNA Depletion Kit (Human/Mouse/Rat)	KD101-11/03	12 rxns / 96 rxns
	TransNGS® Magic rRNA Depletion Kit (Bacteria)	KD401-01/02	12 rxns / 96 rxns
	MagicPure® mRNA Kit	EC511-01/02	24 rxns / 96 rxns
	TransNGS® RNA-Seq Library Prep Kit for Illumina®	KP601-01/02	12 rxns / 96 rxns
	TransNGS® Fast RNA-Seq Library Prep Kit for Illumina®	KP701-01/02	12 rxns / 96 rxns
	TransNGS® Fast Stranded RNA-Seq Library Prep Kit for MGI®	KP801-01/02	12 rxns / 96 rxns
	TransNGS® Circularization Kit For MGI®	KC101-01/02	12 rxns / 96 rxns
	MagicPure® RNA Beads	EC501-01/02/03	1 ml / 5 ml / 60 ml
Library Quantification and Adapters	TransNGS® Library Amplification SuperMix	KA101-01/02	1 mL / 5×1 mL
	TransNGS® Library Quantification Kit for Illumina®	KQ101-01/02	100 rxns / 500 rxns
	TransNGS® Library Quantification qPCR SuperMix	KQ201-01/02	1 mL / 5×1 mL
	TransNGS® Library Quantification DNA Standards (S1-S6)	KS101-21	50%, 120 µL each
	TransNGS® Library Dilution Buffer	KB101-01	5×1 mL
	TransNGS® UDI Indexed Adapter Kit for Illumina®	KI341-01/02	192 rxns / 384 rxns
	TransNGS® 384 UDI Indexed Adapter Kit for Illumina®	KI351-01/02	192 rxns / 384 rxns
	TransNGS® Indexed Adapter Kit for MGI®	KI401-S <sub>1-3</sub> -01/02	192 rxns / 384 rxns
	TransNGS® Index Primers (384) Kit for Illumina®	KI241-01/02	96 rxns / 384 rxns
Pathogenic Microorganism Detection	TransGuard® Disposable Virus Sampling Tube	ES101-01	50 rxns
	TransGuard® Fecal DNA Sampling Tube	ES102-01	50 rxns
	TransGuard® Buccal Swab DNA Preservation Buffer	ES103-01	50 mL
	TransNGS® Host DNA Depletion Kit	EH301-01	50 rxns
	EasyPure® Microbiome DNA Isolation Kit	EE401-01	50 rxns
	MagicPure® 32 Microbiome DNA Isolation Kit	EC107-32-11	32 rxns
	EasyPure® Viral DNA/RNA Kit	ER201-01/02	50 rxns / 200 rxns
	MagicPure® Fly 96 Viral DNA/RNA Kit	EC331-96	96 rxns
	MagicPure® Stool and Soil Genomic DNA Kit	EC801-11	50 rxns
Modifying Enzymes	T4 DNA Ligase (for NGS)	LL101-01/02	200 µL / 1 mL
	T4 DNA Polymerase	LP201-01/02	150 units / 5×150 units
	T4 Polynucleotide Kinase	LK101-01/02	500 units / 4×500 units
	DNA Polymerase I Klenow Fragment	LE201-01	500 units

# How to navigate through our website?

**Homepage:** When you first visit our website, you'll land on the homepage. This is where you can find an overview of our services, products, or any important announcements. From here, you can navigate to other sections of the website.

Consult our catalog and services

Choose your Country

Log in to your account

Choose your language

View Cart

Search for a product by description, reference, CAS, Number, target...

Categories

- Cell-based assay kits
- Culture flasks
- Culture supports for microscopy
- Detection kits for IHC
- FPLC Columns
- Histological special stains
- Immunoassays
- Lectins
- Next Generation Sequencing (NGS)
- Plant tissue culture
- Radiochemicals
- More categories...

## Search Bar: Look for specific Products and equipment

Search by selecting either RUO or IVD products

Refine your search by choosing from the following options

Apply Filter

Search result : 39126007 product found

Cat#	Description	Cond.
NB-45-00042-100	Super Ni-NTA Agarose Resin	100ml
NB-45-00042-25	Super Ni-NTA Agarose Resin	25ml
NB-45-00058-4	Proteus 1 -step Batch Mini Spin Column Pack	40pc
NB-12-6001-3	NeoLine pipette 2-20 µl	1unit
NB-12-0023C	Mini Centrifuge N500C @10.000rpm (including 6x1.5/2.0ml angle rotor)	1pcs
NB-03-0160	Proteinase K (Powder)	100mg
NB-60-0001	NeoPrep mini	50columns
NB-12-8001-19	Combs for NeoPRO4 mini (1.5mm, 15 wells)	5pieces
NB-12-8001-20	Spacer glasses flat for NeoPRO4 mini (0.75mm, 100*83mm)	5pieces
NB-12-8001-04	Short glasses flat for NeoPRO4 mini (1.0mm, 100*73mm)	10pieces

***“Remember, if you ever get lost or have questions, don't hesitate to reach out to our customer support team. They're here to help you address any concerns you may have.”***



# CliniSciences

## Austria

Company: CliniSciences GmbH  
Address: Sternwartestrasse 76, A-1180  
Wien - Austria  
Telephone: +43 720 115 580  
Fax: +43 720 115 577  
Email: [oessterreich@clinisciences.com](mailto:oessterreich@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Belgium

Company: CliniSciences S.R.L  
Address: Avenue Stalingrad 52, 1000  
Brussels - Belgium  
Telephone: +32 2 31 50 800  
Fax: +32 2 31 50 801  
Email: [belgium@clinisciences.com](mailto:belgium@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Denmark

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [denmark@clinisciences.com](mailto:denmark@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Finland

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [suomi@clinisciences.com](mailto:suomi@clinisciences.com)  
Web: <https://www.clinisciences.com>



## France

Company: CliniSciences S.A.S  
Address: 74 Rue des Suisses, 92000  
Nanterre- France  
Telephone: +33 9 77 40 09 09  
Fax: +33 9 77 40 10 11  
Email: [info@clinisciences.com](mailto:info@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Germany

Company: Clinisciences GmbH  
Address: Wilhelm-Mauser-Str. 41-43,  
50827 Köln - Germany  
Telephone: +49 221 9498 320  
Fax: +49 221 9498 325  
Email: [deutschland@clinisciences.com](mailto:deutschland@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Iceland

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [island@clinisciences.com](mailto:island@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Ireland

Company: CliniSciences Limited  
Address: Ground Floor, 71 lower Baggot street  
Dublin D02 P593 - Ireland  
Telephone: +353 1 6971 146  
Fax: +353 1 6971 147  
Email: [ireland@clinisciences.com](mailto:ireland@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Italy

Company: CliniSciences S.r.l  
Address: Via Maremmana inferiore 378  
Roma 00012 Guidonia Montecelio - Italy  
Telephone: +39 06 94 80 56 71  
Fax: +39 06 94 80 00 21  
Email: [italia@clinisciences.com](mailto:italia@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Netherlands

Company: CliniSciences B.V.  
Address: Kraijenhoffstraat 137A,  
1018RG Amsterdam, - Netherlands  
Telephone: +31 85 2082 351  
Fax: +31 85 2082 353  
Email: [nederland@clinisciences.com](mailto:nederland@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Norway

Company: CliniSciences AS  
Address: c/o MerVerdi Munkeordunet 10  
1164 Oslo - Norway  
Telephone: +47 21 988 882  
Email: [norge@clinisciences.com](mailto:norge@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Poland

Company: CliniSciences sp.z.o.o.  
Address: ul. Rotmistrza Witolda Pileckiego 67  
lok. 200 - 02-781 Warszawa -Poland  
Telephone: +48 22 307 0535  
Fax: +48 22 307 0532  
Email: [polska@clinisciences.com](mailto:polska@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Portugal

Company: Quimigen Unipessoal LDA  
Address: Rua Almada Negreiros, Lote 5, Loja 14,  
2615-275 Alverca Do Ribatejo - Portugal  
Telephone: +351 30 8808 050  
Fax: +351 30 8808 052  
Email: [info@quimigen.com](mailto:info@quimigen.com)  
Web: <https://www.quimigen.pt>



## Spain

Company: CliniSciences Lab Solutions  
Address: C/ Hermanos del Moral 13  
(Bajo E), 28019, Madrid - Spain  
Telephone: +34 916 750 700  
Fax: +34 91 269 40 74  
Email: [espana@clinisciences.com](mailto:espana@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Sweden

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email:  [sverige@clinisciences.com](mailto: sverige@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Switzerland

Company: CliniSciences AG  
Address: Fracht Ost Flughafen Kloten  
CH-8058 Zürich - Switzerland  
Telephone: +41 (044) 805 76 81  
Fax: +41 (044) 805 76 75  
Email: [switzerland@clinisciences.com](mailto:switzerland@clinisciences.com)  
Web: <https://www.clinisciences.com>



## UK

Company: CliniSciences Limited  
Address: 11 Progress Business center, Whittle  
Parkway, SL1 6DQ Slough- United Kingdom  
Telephone: +44 (0)1753 866 511  
or +44 (0) 330 684 0982  
Fax: +44 (0)1753 208 899  
Email: [uk@clinisciences.com](mailto:uk@clinisciences.com)  
Web: <https://www.clinisciences.com>



## USA

Company: CliniSciences LLC  
Address: c/o Carr Riggs Ingram,  
500 Grand Boulevard, Suite 210 Miramar  
Beach, FL 32550- USA  
Telephone: +1 850 650 7790  
Fax: +1 850 650 4383  
Email: [usa@clinisciences.com](mailto:usa@clinisciences.com)  
Web: <https://www.clinisciences.com>

