

TRANSGEN

AQ7

AQ7 Series Probe-Based qPCR Products

Distributed by
cliniSciences

AQ7 Series Probe-Based qPCR Solutions – Precision for Every Testing Need

Quantitative Real-time PCR (qPCR) is a vital method for studying gene expression levels. It utilizes changes in fluorescence signals to monitor the amount of amplified product in real-time during each PCR cycle, enabling quantitative analysis of the initial template through the relationship between Ct values and standard curves. Currently, qPCR is widely applied across various fields of molecular biology, including: Diagnosis of human infectious diseases, Detection of pathogenic genes in animals, Inspection and quarantine of livestock and poultry products, Quality control and identification of biological products, etc.


The two most commonly used qPCR detection methods are the SYBR Green dye-based method and the fluorescent probe-based method. The non-specific SYBR Green dye method is more cost-effective but may yield false-positive results due to interference from non-specific PCR amplification products. In contrast, the probe-based method uses fluorescent probes to ensure synchronous accumulation of fluorescence signals with PCR product formation. This approach offers higher sensitivity, superior specificity, and simpler operation, making it the most widely used technique in clinical diagnostics today.

With 19 years of expertise in enzyme engineering and product development, TransGen Biotech has been consistently delivering high-quality, cost-effective solutions to IVD manufacturers, scientific institutions, and third-party testing organizations. Our products have been successfully applied in: Gene expression analysis, Pathogen detection, Genotyping, Copy number variation analysis, empowering researchers to achieve groundbreaking scientific advancements.

CONTENTS

02	10
Basic	Digital PCR
05	10
Extraction-free qPCR	Low Host Nucleic Acid Residue
06	12
Disease Diagnosis	Lyophilizable Series
08	
Genotyping	

AQ7 Series Probe-based qPCR Product Selection Guide



Basic

The *PerfectStart*[®] series probe-based qPCR reagents are high-performance premixed solutions specifically designed for precise and sensitive qPCR detection. The core enzyme is the *PerfectStart*[®] series hot-start DNA polymerase, a genetically engineered Taq DNA polymerase that utilizes antibody-mediated inactivation to effectively block polymerase activity at low temperatures, preventing non-specific amplification. When combined with an optimized reaction buffer, this system significantly enhances template affinity and qPCR amplification performance. It delivers improved sensitivity for low-copy targets and demonstrates excellent tolerance to impurities in nucleic acids (including PCR inhibitors), making it suitable for a wide range of detection scenarios.

PerfectStart[®] II Probe qPCR SuperMix (AQ711)

PerfectStart[®] II Probe qPCR SuperMix UDG (AQ712)


Features

- Blocking by 3 antibodies, with high specificity, high sensitivity, and high amplification efficiency.
- A specially optimized qPCR reaction buffer for enhanced sensitivity and specificity.
- Comes with a Passive Reference Dye (to minimize well-to-well differences) compatible with different qPCR instruments, guaranteeing reliable results.
- Wide range of applications: Successfully used for detecting African swine fever virus, Pseudorabies virus, *Vibrio harveyi*, *Enterocytozoon hepatopenaei*, and *Vibrio parahaemolyticus*, etc.
- Excellent stability: Maintains robust amplification performance even after repeated freeze-thaw cycles or storage in premix, at room temperature, or at 37°C.
- dUTP/UDG system: Effectively prevents carry-over contamination caused by PCR products, ensuring accurate results.

Application


Suitable for probe-based qPCR detection of DNA from animal, plant, and microbial, demonstrating superior species compatibility.

High sensitivity and amplification efficiency



Using different concentrations of PRRSV cDNA (10 pg, 0.1 pg, 0.01 pg) and BVDV cDNA (100 pg, 1 pg, 0.1 pg) as templates, amplification was performed with products from TransGen, Company V, and Company T. The results demonstrated that TransGen's product achieved a sensitivity of 0.1 pg, with superior amplification efficiency compared to Company V and Company T.

Wide range of applications



TransGen's product successfully detected DNA from *Vibrio harveyi*, *Enterocytozoon hepatopenaei* (EHP), and *Vibrio parahaemolyticus* across serial 10-fold dilutions (10^7 to 10^4 copies/mL)

Distributed by
CliniSciences

PerfectStart[®] III Fast Probe qPCR SuperMix UDG (AQ722)


Features

- Antibody-blocking technology with high specificity, sensitivity, and amplification efficiency, demonstrating strong inhibitor resistance and broad species compatibility.
- A specially optimized qPCR reaction buffer for enhanced sensitivity and specificity.
- Ultra-fast detection: Detection can be completed in as little as 35 minutes.
- High amplification capability: Suitable for multiplex PCR applications.
- Excellent stability—no performance degradation after 25 freeze-thaw cycles, 1-week storage at 37°C, or 15 days at room temperature.
- dUTP/UDG system: Effectively prevents carry-over contamination caused by PCR products, ensuring accurate results.
- Supports full premix formulation of primers and probes.
- With universal passive reference dye for signal calibration across wells.

Application


Suitable for complex sample detection, with validated applications in HPV, HBV, and animal pathogen testing.

Rapid amplification



Using serially diluted African swine fever virus (ASFV) genomic DNA (gDNA) as template, rapid qPCR was performed with TransGen's kit. Results demonstrated TransGen's ultra-fast amplification capability, requiring only 1-second extension time and completing the entire process in 35 minutes.

High GC tolerance



Using pseudorabies virus (PRV) gDNA (69% GC content) as template, qPCR was performed with TransGen and Company V kits. Results demonstrated TransGen's superior amplification performance for complex templates.

Distributed by
CliniSciences

PerfectStart® IV Fast Probe qPCR SuperMix UDG (AQ732)

Features


- High specificity: Utilizes next-generation hot-start Taq DNA polymerase with advanced antibody-blocking technology.
- High inhibitor tolerance: Enables direct amplification from blood, tissue homogenates, swabs, and other complex samples.
- Fully premix-compatible: Supports complete primer/probe premixing.
- Excellent stability: Maintains consistent performance after 50 freeze-thaw cycles, 30-day storage at 37°C, or 30-day stress testing at 37°C.
- Ultra-fast PCR: Completes 45 amplification cycles in 35 minutes.
- dUTP/UDG system: Effectively prevents carry-over contamination caused by PCR products, ensuring accurate results.

Application


Demonstrates superior performance with swab, blood, and impurity-containing samples. Successfully applied to challenging detection scenarios including African swine fever (ASFV) in low-quality extracts, HPV, lung cancer, and respiratory bacterial pathogens.

High impurity tolerance

Positive tissue homogenate



2% swine blood




— TransGen — Company V


Using ASFV-positive tissue homogenate and healthy swine blood added with ASFV gDNA as templates separately, qPCR targeting the ASFV national standard gene was performed with TransGen and Company V kits. Results demonstrated TransGen's superior inhibitor resistance.

Excellent direct PCR performance


M gene




Q gene



N gene



Z gene



— Swab — Nucleic acid release reagent — Water

Using nucleic acid release reagent-treated oral/nasal swabs samples from healthy pigs, nucleic acid release reagent and water as diluents, ASFV gDNA was diluted to equal nucleic acid concentrations. Quadruplex qPCR with TransGen's kit demonstrated reliable direct amplification from swab samples, proving compatibility with swab-based direct detection.

Extraction-Free qPCR

Extraction-Free qPCR is a reaction that enables direct amplification from crude samples without prior nucleic acid extraction. TransGen Biotech's extraction-free fluorescent quantitative PCR rapid detection products utilize highly robust enzymes and specially formulated buffers to perform nucleic acid amplification directly from lysed biological samples while ensuring result accuracy. This technology significantly reduces detection time, saves labor costs and experimental resources, and improves testing efficiency.


PerfectStart® V Fast Direct Probe qPCR SuperMix UDG (AQ742)

Features

- Simple workflow: Enables direct target gene detection without nucleic acid extraction or purification steps.
- Enhanced inhibitor resistance: Incorporates next-generation Taq DNA polymerase with improved template affinity and qPCR performance, demonstrating high tolerance for impurity-containing samples.
- Rapid amplification: Compatible with fast-cycling protocols, completing PCR in 35 minutes.
- dUTP/UDG system: Effectively prevents carry-over contamination caused by PCR products, ensuring accurate results.

Application

Extraction-free amplification for whole blood (anticoagulated), plasma, serum, and throat swab samples.



Using healthy swine whole blood of different concentration and water as templates (added with ASFV-positive nucleic acid) for qPCR, the results showed that AQ742 performed almost identically in samples containing 0.5% whole blood compared to blood-free samples. Even under conditions with 10% whole blood, it remained capable of effective detection, demonstrating strong resistance to inhibition and suitability for direct amplification of blood samples.


Disease Diagnosis

Probe-based qPCR technology offers rapid development, high sensitivity, and excellent specificity, making it well-suited for: Rapid screening and diagnosis of emerging infectious diseases, detection of viral infections with asymptomatic or early-stage presentations, and differential diagnosis of pathogens causing similar clinical symptoms. Multiplex nucleic acid detection technology enables simultaneous detection of multiple targets with high throughput and cost efficiency, meeting the demands of large-scale screening and quantitative analysis. TransGen Biotech has launched a series of detection kits for HPV, HBV, and animal pathogens, delivering fast, multiplex, and precise diagnostic solutions.

PerfectStart® III Fast Probe qPCR SuperMix UDG (P) (AQ722-S1)

Application

Suitable for ultra-multiplex HPV detection



Using AQ722-S1 for multiplex HPV testing, the results demonstrated that AQ722-S1 successfully achieved detection at 20 copies/reaction, exhibiting high sensitivity and excellent specificity.


PerfectStart® III Fast Probe qPCR SuperMix UDG(B) (AQ722-S2)

Application

HBV detection

Detection sensitivity: 1 IU/reaction

Using healthy human gDNA as template: 0/32



HBV virus was detected by AQ722-S2. The AQ722-S2 demonstrated reliable detection at 1 IU/reaction, with no false positives for healthy human gDNA, offering high sensitivity and superior specificity.

PerfectStart® IV Fast Probe qPCR SuperMix UDG(P) (AQ732-S1)

This product is a pre-mix reagent specially designed for ultra-multiplex detection of HPV using probe-based qPCR. It is supplied at a 2× concentration, simply add primers, probe and template to perform multiplex PCR reaction in one reaction well. This product contains PerfectStart® IV Taq hot-Start DNA polymerase (Genetically modified fourth generation Taq DNA polymerase), which significantly improves template affinity, qPCR amplification performance and sensitivity. In addition to amplifying extracted nucleic acids, this product facilitates direct amplification from crude lysates of cervical scrapings, eliminating the need for DNA purification. The dUTP/UDG in the SuperMix can function at room temperature to eliminate carry-over contamination caused by PCR products and aerosol to ensure the accuracy of results.

Application

Ultra-multiplex HPV detection directly from swab samples

PerfectStart® IV Fast Probe qPCR SuperMix UDG (AQ732-S2)

This product is a pre-mix reagent specially designed for animal disease detection using probe-based qPCR. It is supplied at a 2× concentration, simply add primers, probe and template to perform singleplex or multiplex PCR reaction in one reaction well. This product contains PerfectStart® IV Taq hot-start DNA polymerase (Genetically modified fourth generation Taq DNA polymerase), which significantly improves template affinity, qPCR amplification performance and detection sensitivity for inhibitor-containing or low-concentration samples. The dUTP/UDG in the SuperMix can function at room temperature to eliminate carry-over contamination caused by PCR products and aerosol to ensure the accuracy of results.

Application

Animal disease detection

Genotyping

Genotyping is the process of identifying genetic variations among individuals within a population. Single nucleotide polymorphisms (SNPs), or point mutations, are the most common type of genetic variation and account for a significant portion of phenotypic diversity among individuals. As third-generation genetic markers, SNPs can reveal differences in genetic material between individuals. Some SNP loci may also affect gene function, leading to phenotypic changes or even disease. Therefore, SNP detection is an essential component of genetic research and related fields. There are various technologies available for SNP detection, including PCR-based methods, sequencing, and microarrays. TransGen offers three methods for SNP detection: ARMS-qPCR, TaqMan SNP probe assay and high-resolution melt curve analysis.


PerfectStart® IV Multiplex Probe ARMS qPCR SuperMix UDG (AQ734)

Features

- High sensitivity and specificity: This product contains newly upgraded PerfectStart® IV hot-start Taq DNA polymerase, blocking by 3 antibodies, combined with the reaction buffer optimized for amplification-refractory mutation system (ARMS), enhancing both sensitivity and specificity.
- Easy to use: Supplied at a 2x concentration, simply add primers, probe and template to perform multiplex reaction.
- dUTP/UDG system: Effectively prevents carry-over contamination caused by PCR products, ensuring accurate results.

Application

Precision molecular diagnostics in oncology, like EGFR mutation detection in tumor targeted therapy



PerfectStart® V Probe SNP Genotyping qPCR SuperMix UDG (AQ744)


Features

- High specificity: The newly upgraded PerfectStart® IV hot-start Taq DNA polymerase, a next-generation genetically engineered Taq DNA polymerase, offers high specificity in genotyping detection across various sample types.
- Easy to use: Supplied at a 2x concentration, simply add primers, probe and template to perform multiplex reaction.
- dUTP/UDG system: Effectively prevents carry-over contamination caused by PCR products, ensuring accurate results.

Application

Widely used in large-scale SNP detection and gene expression analysis

Distributed by
CliniSciences



The human MTHFR C677T polymorphism was detected using AQ744. AQ744 can effectively distinguish all three genotypes (CC, CT, TT) at the MTHFR C677T locus. This enables more efficient and accurate support for folate metabolism screening and anomaly management.

PerfectStart® V Multiplex Probe Melt qPCR Mix UDG (AQ745)


Features

- High uniformity and specificity: Utilizing specially processed PerfectStart® V hot-start Taq DNA polymerase, this system delivers excellent performance in melting curve peak clarity, uniformity, and specificity in ultra-multiplex amplification.
- Easy to use: Supplied at a 2x concentration, simply add primers, probe and template to perform multiplex reaction.
- dUTP/UDG system: Effectively prevents carry-over contamination caused by PCR products, ensuring accurate results.

Application


Pathogen detection and identification (such as HPV genotyping), gene mutation analysis, drug resistance gene detection, etc.

ALDH2 gene polymorphism detection



AQ745 enables rapid, efficient, and accurate genotyping of ALDH2, providing critical guidance for the timely treatment of angina with nitroglycerin, alcohol consumption recommendations, and assessing cancer risks.

Ultra-multiplex HPV genotyping



AQ745 achieves highly multiplexed genotyping of various HPV subtypes, offering precise diagnostic support for HPV infection detection.

Distributed by
CliniSciences

Digital PCR

Digital PCR (dPCR) is a nucleic acid absolute quantification technology. Unlike qPCR, dPCR enables direct calculation of the number of DNA molecules, providing absolute quantification of the original sample. Unlike methods relying on the cycle threshold (Ct) of amplification curves, dPCR is independent of amplification efficiency, offering superior accuracy and reproducibility. It is widely applied in fields where Ct-based methods fall short, including copy number variation (CNV) analysis, mutation detection, gene expression studies, validation of next-generation sequencing (NGS) results, miRNA expression profiling, single-cell gene expression analysis.

PerfectStart[®] CHIPS dPCR SuperMix (AQ715)

This product is a 2x-premixed solution designed for chip-based digital PCR, enabling amplification and absolute quantification of DNA samples. Simply add primers, probe and template to perform detection. It exhibits high amplification efficiency and anti-interference capability, demonstrating broad applicability in tumor and various pathogen detection. The dUTP/UDG in the SuperMix can function at room temperature to eliminate carry-over contamination caused by PCR products and aerosol to ensure the accuracy of results.

Application: Chip-based dPCR

Ultra-Low Host Nucleic Acid Residue

Commercial molecular enzymes are typically produced through expression in engineered bacterial strains (e.g., *E. coli*), which often results in residual host genomic DNA contamination. Additionally, external nucleic acids may be introduced during reagent preparation from environmental or human sources, leading to potential DNA contamination in enzyme products. During pathogen detection, these contaminating background nucleic acids may be co-amplified with target sequences, compromising the reliability of qPCR results.

PerfectStart[®] Taq Low DNA qPCR SuperMix UDG (AQ714)

Features

- Minimal host DNA residue: No detectable genomic DNA from 16 common pathogenic bacteria (including *E. coli*, *P. aeruginosa*, and *S. aureus*) in NTC controls.
- Excellent amplification efficiency: Achieves 100% efficiency across a broad dynamic range (10-10⁵ copies/μl).
- Superior stability: Maintains stable performance under various conditions including repeated freeze-thaw cycles, pre-mixed reagent formulations, room temperature, and 37°C.
- dUTP/UDG system: Effectively prevents carry-over contamination caused by PCR products, ensuring accurate results.


Application

Specifically designed for applications requiring ultra-low host nucleic acid residue, including pathogen detection assays and quality control testing of biological products (particularly *E. coli*-derived preparations)


Distributed by
CliniSciences

Contamination control

E. coli residue



Human DNA contamination




— AQ712
— AQ714


NTC amplification of *E. coli* and human gDNA systems using AQ714 and AQ712 demonstrated undetectable host DNA contamination in AQ714, confirming its ultra-low residual host DNA content.

Wide range of applications


African swine fever virus (ASFV)




HPV-16



Group A streptococcus (GAS)



Brucella



AQ714 successfully amplifies ASFV, HPV-16, GAS, and Brucella, confirming wide template adaptability.

Distributed by
CliniSciences

Lyophilizable Series

TransGen Biotech can provide lyophilizable formulations for full range of IVD raw materials, supporting multiple lyophilization forms including powder/microspheres in 8-tube strips or microspheres in vials. Traditional PCR reaction systems require low-temperature storage, making them susceptible to temperature fluctuations and contamination. This results in high transportation costs and operational inconveniences, limiting their application in resource-limited settings and point-of-care testing. In contrast, lyophilized PCR reagents offer enhanced stability, reduced cold-chain dependence, improved operational efficiency and accuracy and lower costs. These advantages have made lyophilized reagents increasingly popular.

- Comprehensive offerings: Full range of IVD raw materials available.
- Compatible with multiple lyophilization forms: Powder/microspheres in 8-tube strips or microspheres in vials.
- Excellent stability: Maintains performance after 45 days at 45°C, stable for 60 days at 37°C.
- Customization: Custom lyophilization solutions.

PerfectStart® II Probe qPCR Lyo-SuperMix UDG (AQ713)

Application

Lyophilizable version of AQ712. Suitable for probe-based qPCR detection of DNA from animals, plants, and microorganisms.

PerfectStart® III Fast Probe qPCR Lyo-SuperMix UDG (AQ723)

Application

Lyophilizable version of AQ722. Suitable for the detection of complex samples, validated applications include HPV, HBV, and animal pathogen detection.

PerfectStart® IV Fast Probe qPCR Lyo-SuperMix UDG (AQ733)


Application

Lyophilizable version of AQ732. Excellent performance for swabs, blood and inhibitor-containing samples, validated applications include African Swine Vever (samples with poor extraction), HPV, lung cancer, and respiratory bacteria detection.


Compatible with multiple lyophilization forms



Microspheres in 8-tube strips




Powder in 8-tube strips




Microspheres in vials

Distributed by
CliniSciences

Excellent stability



The lyophilized microspheres of AQ723, after being subjected to thermal acceleration at 37°C for 60 days and 0 days, were reconstituted and used for quadruplex amplification with medium and low concentrations of ASFV-gDNA as templates. The results showed no performance change after 60 days of thermal acceleration at 37°C, indicating excellent stability.



The three forms of AQ713 (lyophilized microspheres, lyophilized powder and pre-lyophilization solution) stored at 45°C for 45 days were used for gene amplification. The results demonstrated no performance change after 45 days of storage at 45°C, confirming stable product performance.

Distributed by
CliniSciences

How to navigate through our website?

Homepage: When you first visit our website, you'll land on the homepage. This is where you can find an overview of our services, products, or any important announcements. From here, you can navigate to other sections of the website.

Consult our catalog and services

Choose your Country

Log in to your account

Choose your language

View Cart

Search for a product (description, reference, CAS number, target...)

upcyte® cells
Expanded primary cells

Categories

- Cell-based assay kits
- Culture flasks
- Culture supports for microscopy
- Detection kits for IHC
- FPLC Columns
- Histological special stains
- Immunoassays
- Lectins
- Next Generation Sequencing (NGS)
- Plant tissue culture
- Radiochemicals
- More categories...

All our products >

Search Bar: Look for specific Products and equipment

Search by selecting either RUO or IVD products

Refine your search by choosing from the following options

Apply Filter

Search for a product (description, reference, CAS number, target...)

Search result : 39126007 product found

Cat#	Description	Cond
NB-45-00042-100	Super Ni-NTA Agarose Resin	100ml
NB-45-00042-25	Super Ni-NTA Agarose Resin	25ml
NB-45-00058-4	Proteus 1 -step Batch Mini Spin Column Pack	40pc
NB-12-6001-3	NeoLine pipette 2-20 µl	1unit
NB-12-0023C	Mini Centrifuge N500C @10,000rpm (including 6x1.5/2.0ml angle rotor)	1pcs
NB-03-0160	Proteinase K (Powder)	100mg
NB-60-0001	NeoPrep mini	50columns
NB-12-8001-19	Combs for NeoPRO4 mini (1.5mm, 15 wells)	5pieces
NB-12-8001-20	Spacer glasses flat for NeoPRO4 mini (0.75mm, 100*83mm)	5pieces
NB-12-8001-04	Short glasses flat for NeoPRO4 mini (1.0mm, 100*73mm)	10pieces

Refine your search :

RUO CE / IVD

- Conjugate (100) +
- Reactivity species (100) +
- Class (66) +
- Host species (100) +
- Cell name (2) +
- Application (100) +
- Clone (100) +

APPLY FILTERS

REINITIALIZE

SHARE THIS SEARCH

“Remember, if you ever get lost or have questions, don't hesitate to reach out to our customer support team. They're here to help you address any concerns you may have.”

Distributed by
clinisciences

clinisciences

Austria

Company: CliniSciences GmbH
Address: Sternwartestrasse 76, A-1180
Wien - Austria
Telephone: +43 720 115 580
Fax: +43 720 115 577
Email: oesterreich@clinisciences.com
Web: <https://www.clinisciences.com>



Belgium

Company: CliniSciences S.R.L
Address: Avenue Stalingrad 52, 1000
Brussels - Belgium
Telephone: +32 2 31 50 800
Fax: +32 2 31 50 801
Email: belgium@clinisciences.com
Web: <https://www.clinisciences.com>



Denmark

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: denmark@clinisciences.com
Web: <https://www.clinisciences.com>



Finland

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: suomi@clinisciences.com
Web: <https://www.clinisciences.com>



France

Company: CliniSciences S.A.S
Address: 74 Rue des Suisses, 92000
Nanterre- France
Telephone: +33 9 77 40 09 09
Fax: +33 9 77 40 10 11
Email: info@clinisciences.com
Web: <https://www.clinisciences.com>



Germany

Company: Clinisciences GmbH
Address: Wilhelm-Mausier-Str. 41-43,
50827 Köln - Germany
Telephone: +49 221 9498 320
Fax: +49 221 9498 325
Email: deutschland@clinisciences.com
Web: <https://www.clinisciences.com>



Iceland

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: island@clinisciences.com
Web: <https://www.clinisciences.com>



Ireland

Company: CliniSciences Limited
Address: Ground Floor, 71 lower Baggot street
Dublin D02 P593 - Ireland
Telephone: +353 1 6971 146
Fax: +353 1 6971 147
Email: ireland@clinisciences.com
Web: <https://www.clinisciences.com>



Italy

Company: CliniSciences S.r.l
Address: Via Maremmana inferiore 378
Roma 00012 Guidonia Montecelio - Italy
Telephone: +39 06 94 80 56 71
Fax: +39 06 94 80 00 21
Email: italia@clinisciences.com
Web: <https://www.clinisciences.com>



Netherlands

Company: CliniSciences B.V.
Address: Kraaijenhoffstraat 137A,
1018RG Amsterdam, - Netherlands
Telephone: +31 85 2082 351
Fax: +31 85 2082 353
Email: nederland@clinisciences.com
Web: <https://www.clinisciences.com>



Norway

Company: CliniSciences AS
Address: c/o MerVerdi Munkerudtunet 10
1164 Oslo - Norway
Telephone: +47 21 988 882
Email: norge@clinisciences.com
Web: <https://www.clinisciences.com>



Poland

Company: CliniSciences sp.Z.o.o.
Address: ul. Rotmistrza Witolda Pileckiego 67
lok. 200 - 02-781 Warszawa -Poland
Telephone: +48 22 307 0535
Fax: +48 22 307 0532
Email: polska@clinisciences.com
Web: <https://www.clinisciences.com>



Portugal

Company: Quimigen Unipessoal LDA
Address: Rua Almada Negreiros, Lote 5, Loja 14,
2615-275 Alverca Do Ribatejo - Portugal
Telephone: +351 30 8808 050
Fax: +351 30 8808 052
Email: info@quimigen.com
Web: <https://www.quimigen.pt>



Spain

Company: CliniSciences Lab Solutions
Address: C/ Hermanos del Moral 13
(Bajo E), 28019, Madrid - Spain
Telephone: +34 916 750 700
Fax: +34 91 269 40 74
Email: espana@clinisciences.com
Web: <https://www.clinisciences.com>



Sweden

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: sverige@clinisciences.com
Web: <https://www.clinisciences.com>



Switzerland

Company: CliniSciences AG
Address: Fracht Ost Flughafen Kloten
CH-8058 Zürich - Switzerland
Telephone: +41 (044) 805 76 81
Fax: +41 (044) 805 76 75
Email: switzerland@clinisciences.com
Web: <https://www.clinisciences.com>



UK

Company: CliniSciences Limited
Address: 11 Progress Business center, Whittle
Parkway, SL1 6DQ Slough- United Kingdom
Telephone: +44 (0)1753 866 511
or +44 (0) 330 684 0982
Fax: +44 (0)1753 208 899
Email: uk@clinisciences.com
IWeb: <https://www.clinisciences.com>



USA

Company: CliniSciences LLC
Address: c/o Carr Riggs Ingram,
500 Grand Boulevard, Suite 210 Miramar
Beach, FL 32550- USA
Telephone: +1 850 650 7790
Fax: +1 850 650 4383
Email: usa@clinisciences.com
Web: <https://www.clinisciences.com>

