



RNAConnect

UltraMarathonRT Reverse Transcription Kit

A. INTRODUCTION

UltraMarathonRT (uMRT) is a highly optimized variant of MarathonRT with improved overall performance and sensitivity to low RNA input. uMRT is a group II intron encoded reverse transcriptase (RT) with ultra-high processivity and RNA template RNA unwinding capability. Unlike other RTs that require high temperature to disrupt stable structures, uMRT performance is maintained at ambient temperatures, from 20°C to 42°C, which prevents RNA template degradation.

B. KITS COMPONENTS

Part number	Description	Storage Recommendation
200001 or 200002	UltraMarathonRT (20 U/μL)	-20°C (up to 3mo) -80°C (long-term) Stable for 20x freeze/thaw
200011	RT Reaction buffer (2x)	-20°C
200012	uMRT Boost (20x)	-20°C
200013	Nuclease-free water	Room temperature, 4°C, or -20°C

C. ENZYME INFORMATION

- **Unit definition:** One (1) unit is equal to the amount of UltraMarathonRT that incorporates 1 nmole dTTP at 42°C in 30 minutes, when using poly(rA) as the template and Oligo(dT)₁₈ as the primer.
- **Enzyme concentration and storage buffer:** The enzyme is supplied at 20 U/μL in a buffer that contains 25 mM K-HEPES pH 7.5, 300 mM KCl, 10% glycerol and 1 mM DTT.

Distributed by:

CliniSciences Group

D. REACTION BUFFER (2X) COMPOSITION

Components and concentration in 2x:

- 100 mM Tris-HCl (pH 8.3 at 25°C)
- 400 mM KCl
- 8 mM MgCl₂
- 10 mM DTT
- 40% glycerol

E. ADDITIONAL REAGENTS REQUIRED AND NOT PROVIDED:

Reagents	Recommendation
dNTP mix	Use a stock that contains 10mM each of dATP, dCTP, dGTP, and dTTP of high purity (>99%) in a pH 7.0 solution (e.g. New England Biolab, Cat# N0447)
Primers	Oligo-(dT) ₁₈ primer (custom synthesis), random hexamer (e.g. Thermo Fisher, Cat# SO142) or gene-specific primer (custom synthesis)
RNA template	1 pg–2 µg total RNA or 1 pg–500 ng of poly(A)-RNA.

F. STEP-BY-STEP PROCEDURE

a. Anneal RT primers to RNA templates

- Combine the components as indicated in the table below to a nuclease-free microcentrifuge tube.
- Mix gently by tapping the tube. Collect the contents by brief centrifugation. Incubate at 95°C for 30 sec and snap cool on ice to anneal the primer to the template.

Components		Final amount	Volume
Primer	Oligo(dT) ₁₈ (5 µM)	5 pmol	1 µL
	Randomer (10 µM)	10 pmol	
	Gene-specific primer (2 µM)	2 pmol	
Template RNA	Total RNA	1 pg–2 µg	variable
	poly(A)-RNA	1 pg–500 ng	variable
dNTP mix stock, 10 mM each		0.5 mM final concentration	1 µL added to total
Nuclease-free water		Distributed by:	6 µL

CliniSciences Group

b. Prepare the RT reaction mix

- Combine the components as indicated in the table below in a nuclease-free microcentrifuge tube.
- Mix the contents gently by tapping the tube and collect the contents by brief centrifugation.

Components	Standard RNA Input	Low RNA Input
	(> 10 ng total RNA or PolyA RNA)	(< 10 ng total RNA or PolyA RNA)
Nuclease-free water	2 μ L	1 μ L
Reaction Buffer (2x)	10 μ L	
UltraMarathonRT (20 U/ μ L)	1 μ L	
RNaseOUT™ (40 U/ μ L) (optional) (e.g. Thermo Fisher, Cat# 10777019)	1 μ L	
uMRT Boost (20x)	Not needed	1 μ L
Total volume	14 μ L	

c. Carry out the reverse transcription reaction

- Add the RT reaction mix (14 μ L, prepared in Section **b**) to the annealed primer and RNA template (6 μ L, prepared in Section **a**) to make a 20 μ L reaction.
- Mix gently by tapping the tube.
- Incubate the mixture as indicated below.

	For RNA \leq 12 kb	For RNA \geq 12 kb
cDNA Synthesis	42°C for 15 min	42°C for 60 min
Inactivation	95°C for 1 min	

The cDNA should be stored at -20°C until use. (**Note:** for PCR amplification, it is recommended to keep the volume of unpurified reverse transcription product (cDNA) input at 10% or less of the PCR reaction volume.)

Distributed by:

CliniSciences Group

G. OPTIMIZATION AND TROUBLESHOOTING

Potential Issue	Likely Cause	Solution
cDNA yields are low	Enzyme concentration is too low	Enzyme concentration may be increased up to 60 units per reaction for a 20 μ L reaction volume
	Reaction time is too short	Consider extending the duration of the reaction – particularly if conducting it below 42°C
I did not get good results with my low abundant RNA samples	uMRT Boost was not utilized	For applications where the total amount of cellular input RNA is very low (10 ng or less), it is important to use uMRT Boost, an alternative uMRT buffer component that contains a specialized carrier molecule (RNAConnect catalog #R1003).
My PCR reaction did not work well	Too much uMRT Boost was used	<p>uMRT Boost is only needed for very low RNA template input (< 10 ng).</p> <p>High concentrations of uMRT Boost may inhibit some downstream PCR applications.</p> <p>1 μL is the appropriate amount of uMRT Boost (20X) for a 20 μL reaction volume.</p> <p>Clean up the cDNA samples before PCR reaction.</p>

REFERENCES

1. Chen Zhao and Anna Marie Pyle (2016) Crystal structures of a group II intron maturase reveal a missing link in spliceosome evolution. *Nat. Struct. Mol. Biol.* **23**, 558-565.
2. Chen Zhao, Fei Liu and Anna Marie Pyle (2018) An ultra-processive, accurate reverse transcriptase encoded by a metazoan group II intron. *RNA* **24**, 183-195.
3. Li-Tao Guo, Rebecca L Adams, Han Wan, Nicholas C Huston, Olga Potapova, Sara Olson, Christian M Gallardo, Brenton R Graveley, Bruce E Torbett, Anna Marie Pyle (2020) Sequencing and Structure Probing of Long RNAs Using MarathonRT: A Next-Generation Reverse Transcriptase. *J. Mol. Biol.* **432**, 3338-3352.
4. Li-Tao Guo, Sara Olson, Shivali Patel, Brenton R Graveley and Anna Marie Pyle (2022) Direct tracking of reverse-transcriptase speed and template sensitivity: implications for sequencing and analysis of long RNA molecules. *Nucleic Acids Res.* **50**, 6980-6989.

Distributed by:

CliniSciences Group

CliniSciences Group

Austria

Company: CliniSciences GmbH
Address: Sternwartestrasse 76, A-1180
Wien - Austria
Telephone: +43 720 115 580
Fax: +43 720 115 577
Email: osterreich@clinisciences.com
Web: <https://www.clinisciences.com>



Belgium

Company: CliniSciences S.R.L.
Address: Avenue Stalingrad 52, 1000
Brussels - Belgium
Telephone: +32 2 31 50 800
Fax: +32 2 31 50 801
Email: belgium@clinisciences.com
Web: <https://www.clinisciences.com>



Denmark

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: denmark@clinisciences.com
Web: <https://www.clinisciences.com>



Finland

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: suomi@clinisciences.com
Web: <https://www.clinisciences.com>



France

Company: CliniSciences S.A.S
Address: 74 Rue des Suisses, 92000
Nanterre- France
Telephone: +33 9 77 40 09 09
Fax: +33 9 77 40 10 11
Email: info@clinisciences.com
Web: <https://www.clinisciences.com>



Germany

Company: Biotrend Chemikalien GmbH
Address: Wilhelm-Mauser-Str. 41-43,
50827 Köln - Germany
Telephone: +49 221 9498 320
Fax: +49 221 9498 325
Email: info@biotrend.com
Web: <https://www.biotrend.com>



Iceland

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: island@clinisciences.com
Web: <https://www.clinisciences.com>



Ireland

Company: CliniSciences Limited
Address: Ground Floor, 71 lower Baggot street
Dublin D02 P593 - Ireland
Telephone: +353 1 6971 146
Fax: +353 1 6971 147
Email: ireland@clinisciences.com
Web: <https://www.clinisciences.com>



Italy

Company: CliniSciences S.r.l
Address: Via Maremmana inferiore 378
Roma 00012 Guidonia Montecelio - Italy
Telephone: +39 06 94 80 56 71
Fax: +39 06 94 80 00 21
Email: italia@clinisciences.com
Web: <https://www.clinisciences.com>



Netherlands

Company: CliniSciences B.V.
Address: Kraijenhoffstraat 137A,
1018RG Amsterdam, Netherlands
Telephone: +31 85 2082 351
Fax: +31 85 2082 353
Email: nederland@clinisciences.com
Web: <https://www.clinisciences.com>



Norway

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: norge@clinisciences.com
Web: <https://www.clinisciences.com>



Poland

Company: CliniSciences sp.Z.o.o.
Address: ul. Rotmistrza Witolda Pileckiego 67
lok. 200 - 02-781 Warszawa -Poland
Telephone: +48 22 307 0535
Fax: +48 22 307 0532
Email: polska@clinisciences.com
Web: <https://www.clinisciences.com>



Portugal

Company: Quimigen Unipessoal LDA
Address: Rua Almada Negreiros, Lote 5, Loja 14,
2615-275 Alverca Do Ribatejo - Portugal
Telephone: +351 30 8808 050
Fax: +351 30 8808 052
Email: info@quimigen.com
Web: <https://www.quimigen.pt>



Spain

Company: CliniSciences Lab Solutions
Address: C/ Hermanos del Moral 13
(Bajo E), 28019, Madrid - Spain
Telephone: +34 91 269 40 65
Fax: +34 91 269 40 74
Email: espana@clinisciences.com
Web: <https://www.clinisciences.com>



Sweden

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: sverige@clinisciences.com
Web: <https://www.clinisciences.com>



Switzerland

Company: CliniSciences Limited
Address: Marktgasse 18 8302 Kloten -
Switzerland
Telephone: +41 (044) 805 76 81
Fax: +41 (044) 805 76 75
Email: switzerland@clinisciences.com
Web: <https://www.clinisciences.com>



UK

Company: CliniSciences Limited
Address: 11 Progress Business center, Whittle
Parkway, SL1 6DQ Slough- United Kingdom
Telephone: +44 (0)1753 866 511
or +44 (0) 330 684 0982
Fax: +44 (0)1753 208 899
Email: uk@clinisciences.com
Web: <https://www.clinisciences.com>



USA

Company: Biotrend Chemicals LLC
Address: c/o Carr Riggs Ingram,
500 Grand Boulevard, Suite 210 Miramar
Beach, FL 32550- USA
Telephone: +1 850 650 7790
Fax: +1 850 650 4383
Email: info@biotrend-usa.com
Web: <https://www.biotrend-usa.com>

