

The Serum-Free System With Recombinant RocketFuel

Cat# NB-11-0087

The Serum-Free Recombinant Growth Factor System. Complete cell growth and passage system including CSC Media, Attachment Factor, CultureBoost-R, Passage Reagent Group and Cell Freezing Medium. Contains Recombinant Growth Factor only.

Media Component:

Serum-Free Medium - does not contain bovine derived components. Human serum albumin and other components of human origin are specified by the formulation. This medium becomes complete once activated with the included RocketFuel supplement containing Recombinant Growth Factors only. Certified and intended for experimental application.

Growth Supplement Component:

CultureBoost-R is the broad-spectrum supplement used to activate CSC Complete Medium. CultureBoost-R contains Recombinant Growth Factor and Porcine heparin.

PRG Component:

Passage Reagent Group™ (PRG) is a matched set of CSC Certified™ reagents for releasing cells from culture for subculture or freezing. The PRG contains three parts: PRG-1 (EDTA -dPBS Solution), PRG-2 (Trypsin/EDTA -dPBS Solution) and PRG-3 (Trypsin Inhibitor-dPBS Solution). The chelating agent EDTA in PRG-1 prepares for PRG-2, which contains highly purified trypsin. PRG-3 inactivates the protease in PRG-2 and stabilizes the cell membranes.

Cell Freezing Medium Component:

CSC Cell Freezing Medium is a specialized media, when used in conjunction with CSC Passage Reagent Group, provides a beneficial environment for the freeze/thaw cycle of cell cultures, assisting in the minimization of cellular damage during the process.

Attachment Factor Component:

Attachment Factor™ is an extracellular matrix (ECM) product that promotes cell attachment to tissue culture surfaces and encourages correct polarity and cytoskeletal organization. Attachment Factor™ also contributes to correct *in vitro* growth factor and biological response modifier presentation. Use of Attachment Factor™ is critically important when cultures are to be initiated, grown, passaged, or used within the CSC Medium family.

CSC media and reagents are Sterile, made with WFI and all components are cGMP and ISO Compliant.