



1871 N GAFFEY STREET
 SUITE A & B
 SAN PEDRO CA 90731
 www.celprogen.com

Phone: 310 547 3975
 Fax: 310 547 2975
 Email: info@celprogen.com
 stemcells@celprogen.com

DATA SHEET

Human Bone Marrow Derived Stem Cell Culture Extracellular Matrix (ECM)
--

Catalog number: E36093-20

Description: Human Bone Marrow Derived Stem Cell Culture Extracellular Matrix for growth, expansion, and proliferation of Human Bone Marrow Derived Stem Cell in tissue culture with their appropriate media. The Human Bone Marrow Derived Stem Cell ECM is available in the following formats: You may choose the tissue culture ware from the table below; each tissue culture ware is packaged with the number of individual items per package. One can only purchase each tissue culture ware as complete packages and will not be able to mix and match individual tissue culture ware formats. For individual catalog items can be picked from the table below or you may please visit our website at www.celprogen.com and choose the formats from the drop down menu. Each format is tissue culture tested with their appropriate media.

Catalog	Description
E36093-20-T25	Human Bone Marrow Derived Stem Cell Culture Extracellular Matrix - T25 Flask (10/Pk)
E36093-20-T75	Human Bone Marrow Derived Stem Cell Culture Extracellular Matrix - T75 Flask (5/Pk)
E36093-20-T150	Human Bone Marrow Derived Stem Cell Culture Extracellular Matrix - T150 Flasks (5/pk)
E36093-20-T225	Human Bone Marrow Derived Stem Cell Culture Extracellular Matrix - T225 Flasks (5/pk)
E36093-20-4Well	Human Bone Marrow Derived Stem Cell Culture ECM - 4 Well Plates (5/Pk)
E36093-20-4MCS	Human Bone Marrow Derived Stem Cell Culture Extracellular Matrix - 4 Well Microscope Chamber Slides
E36093-20-MS4	Human Bone Marrow Derived Stem Cell Culture Extracellular Matrix 4 Well Microscope Slide
E36093-20-6Well	Human Bone Marrow Derived Stem Cell Culture Extracellular Matrix - 6 Well Plates (5/Pk)
E36093-20-8MCS	Human Bone Marrow Derived Stem Cell Culture Extracellular Matrix - 8 Well Microscope Chamber Slides
E36093-20-MS8	Human Bone Marrow Derived Stem Cell Culture Extracellular Matrix 8 well Microscope Slide
E36093-20-MS10	Human Bone Marrow Derived Stem Cell Culture- 10 Well Microscope Slides
E36093-20-12Well	Human Bone Marrow Derived Stem Cell Culture Extracellular Matrix - 12 Well Plates (5/Pk)
E36093-20-MS18	Human Bone Marrow Derived Stem Cell Culture- 18 Well Microscope Slides
E36093-20-24Well	Human Bone Marrow Derived Stem Cell Culture Extracellular Matrix - 24 Well Plates (5/Pk)
E36093-20-48Well	Human Bone Marrow Derived Stem Cell Culture Extracellular Matrix - 48 Well Plates (5/Pk)



1871 N GAFFEY STREET
 SUITE A & B
 SAN PEDRO CA 90731
 www.celprogen.com

Phone: 310 547 3975
 Fax: 310 547 2975
 Email: info@celprogen.com
 stemcells@celprogen.com

E36093-20-96Well	Human Bone Marrow Derived Stem Cell Culture Extracellular Matrix - 96 Well Plates (5/Pk)
E36093-20-96WW	Human Bone Marrow Derived Stem Cell Culture Extracellular Matrix - 96 White Well Plates (5/pk)
E36093-20-96BW	Human Bone Marrow Derived Stem Cell Culture Extracellular Matrix - 96 Black Well Plates (5/pk)
E36093-20-384WW	Human Bone Marrow Derived Stem Cell Culture Extracellular Matrix - 384 White Well Plates (5/pk)
E36093-20-384BW	Human Bone Marrow Derived Stem Cell Culture Extracellular Matrix - 384 Black Well Plates (5/pk)
E36093-20-CS12	Human Bone Marrow Derived Stem Cell Culture Extracellular Matrix- 12mm coverslips (10/pk)
E36093-20-CS15	Human Bone Marrow Derived Stem Cell Culture Extracellular Matrix- 15mm coverslips (10/pk)
E36093-20-CS18	Human Bone Marrow Derived Stem Cell Culture Extracellular Matrix- 18mm coverslips (10/pk)
E36093-20-CS22	Human Bone Marrow Derived Stem Cell Culture Extracellular Matrix- 22mm coverslips (10/pk)
E36093-20-CS25	Human Bone Marrow Derived Stem Cell Culture Extracellular Matrix- 25mm coverslips (10/pk)
E36093-20-PD6	Human Bone Marrow Derived Stem Cell Culture Extracellular Expansion Matrix - 6cm Petri Dish (5/pk)
E36093-20-PD10	Human Bone Marrow Derived Stem Cell Culture Extracellular Expansion Matrix - 10cm Petri Dish (5/pk)
E36093-20-6W-GB	Human Bone Marrow Derived Stem Cell Culture Extracellular Matrix 6 Well Glass Bottom Plates (5/pk)
E36093-20-96W-GB	Human Bone Marrow Derived Stem Cell Culture Extracellular Matrix 96 Well Glass Bottom Plates (5/pk)
E36093-20-3D-6Well	Human Bone Marrow Derived Stem Cell Culture 3D Cell Culture System 6 Well Format (1/pk)
E36093-20-3D-12Well	Human Bone Marrow Derived Stem Cell Culture 3D Cell Culture System 12 Well Format (1/pk)
E36093-20-3D-24Well	Human Bone Marrow Derived Stem Cell Culture 3D Cell Culture System 24 Well Format (1/pk)

Storage Conditions: 2- 8 °C long term storage 18 months

Room temperature for short term storage 12 months

Product Orders: Before submitting an order, you will be asked to read and accept the terms and conditions of Celprogen's Material Transfer Agreement (MTA).

Permits & Forms: In addition to the MTA mentioned above, other CELPROGEN and/or regulatory permits may be required for the transfer of this CELPROGEN material. Anyone purchasing CELPROGEN material is ultimately responsible for obtaining the permits.

Notices &



1871 N GAFFEY STREET
SUITE A & B
SAN PEDRO CA 90731
www.celprogen.com

Phone: 310 547 3975
Fax: 310 547 2975
Email: info@celprogen.com
stemcells@celprogen.com

Disclaimers:

CELPROGEN products are intended for laboratory research purposes only. They are not intended for use in Humans. This product, Human Bone Marrow Derived Stem Cell Culture Expansion Extracellular Matrix, is established and manufactured by CELPROGEN Inc., and is for Research Use Only. This product is not for re-sale or may not be transferred to third party prior to written request and approval by CELPROGEN Inc.